









Catalog # 15-0136 (R)-(+)-5,5'-Bis(diphenylphosphino)-4,4'-bi-1,3-benzodioxole, min. 98% (R)-SEGPHOS®




Note: Manufactured under license of Takasago patent. Takasago SEGPHOS® Ligand Kit component.


Technical Notes:

1. Biaryl bisphosphine ligand with narrow dihedral angle. The SEGPHOS® ligand has been applied to a variety of metal catalyzed reactions. In many cases, yields and enantioselectivities, exceed results obtained earlier using BINAP.^{1,2}
2. As ruthenium complex, SEGPHOS® generally gives higher levels of chiral induction in asymmetric hydrogenations of α,β , and γ -functionalized ketones. See ruthenium complexes [44-0096](#), [44-0518](#), [44-0168](#).
3. Used in Rh-catalyzed transformations such as: (a) 1,4-addition of boronic acids to coumarins,⁴ (b) addition of titanium reagents to imines,⁶ (c) cotrimerization of alkenes and acetylenes,¹⁰ (d) double [2+2] cycloaddition,¹¹ (e) indanone formation.^{12a,b}
4. Used in Pd-catalyzed transformations such as: (a) cycloaddition of 1,6-enyne,⁵ (b) arylyative cyclization of allenyl aldehydes with boronic acids,¹³ (c) synthesis of chromans.¹⁴
5. Used in Cu-catalyzed transformations such as: (a) nitroso Diels-Alder,⁷ (b) reductive aldol condensation,⁹ (c) conjugate reduction of unsaturated sulfones,¹⁵ and phosphonates.¹⁶
6. Iridium-catalyzed asymmetric hydrogenation of quinolines activated by chloroformates.
7. Iridium-catalyzed asymmetric transfer hydrogenation used in polyketide construction.¹⁷
8. Rhodium-catalyzed asymmetric hydroarylation of 3-pyrrolines.¹⁸
9. Palladium-catalyzed regio- and enantioselective dearomatization of pyrroles to 2H-pyrroles.¹⁹
10. Rhodium-catalyzed asymmetric synthesis of cyclopentanols.²⁰
11. Silver-catalyzed asymmetric Mannich-type reaction.²¹


Tech. Note (3)
Ref. (12a)Tech. Note (4)
Ref. (5)Tech. Note (4)
Ref. (14)Tech. Note (5)
Ref. (7)Tech. Note (5)
Ref. (9)Tech. Note (6)
Ref. (8)Tech. Note (8)
Ref. (18)



Tech. Note (9)
Ref. (19)



Tech. Note (10)
Ref. (20)



Tech. Note (11)
Ref. (21)

References:

1. *Tetrahedron*, **2005**, 61, 5405.
2. *Topics Organometal. Chem.*, **2004**, 6, 63.
3. *Org. Lett.*, **2003**, 5, 5043.
4. *Org. Lett.*, **2005**, 7, 2285.
5. *Angew. Chem. Int. Ed.*, **2001**, 40, 249.
6. *Angew. Chem. Int. Ed.*, **2004**, 43, 6125.
7. *J. Am. Chem. Soc.*, **2004**, 126, 4128.
8. *Angew. Chem. Int. Ed.*, **2006**, 45, 2260.
9. *Org. Lett.*, **2005**, 7, 4225.
10. *Org. Lett.*, **2008**, 10, 2849.
11. *Org. Lett.*, **2008**, 10, 2825.
12. (a) *Angew. Chem. Int. Ed.*, **2007**, 46, 3735. (b) *Org. Lett.*, **2006**, 8, 3379.
13. *Org. Lett.*, **2008**, 10, 1047.
14. *Org. Lett.*, **2009**, 11, 4752.
15. *Angew. Chem.*, **2007**, 119, 3393.
16. *Tetrahedron. Lett.*, **2009**, 50, 6720.
17. *J. Am. Chem. Soc.*, **2011**, 133, 12795.
18. *J. Am. Chem. Soc.*, **2013**, 135, 10990.
19. *J. Am. Chem. Soc.*, **2014**, 136, 6590.
20. *J. Am. Chem. Soc.*, **2014**, 136, 7217.
21. *Org. Lett.*, **2014**, 16, 86.