

Sold in collaboration with HepatoChem

96-7560 EvoluChem™ Photocatalytic Alkylation Kit

1 kit

Product overview:

The trifluoroborate alkylation reaction (Minisci reaction)¹ is a powerful late stage functionalization tool. Our kit allows convenient, one-step production of eight different analogues of a lead compound in mg quantities. Each reaction vial contains 75 μmol of trifluoroborate alkylation reagent (pre-weighed) and a stirring bar to react with 50 μmol of substrate. C-H functionalization will primarily occur on electron-deficient heteroarenes at one or several positions.

Benefits


- Facilitates screen of photochemical reaction conditions
- Enables up to 32 reaction conditions simultaneously
- Save substrate using low scale reaction conditions
- Save time on optimization

Material required, but not supplied

- Customer supplied substrate
- Customer supplied reaction solvent(s)
- Trifluoroacetic acid
- EvoluChem™ PhotoRedOx Box
- EvoluChem™ Light Source 18W-450 nm
- EvoluChem™ Vial Holder (98-7650)
- Nitrogen or Argon line for sparging solvents with two needles
- DMSO
- Stirring plate
- 500 μL - 1.0 ml syringe

Recommendations

- Safety personal protection such as gloves, safety glasses and lab coat should be worn at all times.
- Always use a clean and dry syringe to add and transfer solution.



Storage and Stability

- Store at 2-8°C in dark.
- Stable for 12 months.

Kit Contents

| Description | Label | Quantity | Amount |
|--|--|-----------|-----------------------|
| Potassium cyclopropyltrifluoroborate / Potassium persulfate | cyclopropyl-BF ₃ K/ K ₂ S ₂ O ₈ | 2 x vials | 75 μmol / 100 μmol |
| Potassium cyclobutyltrifluoroborate / Potassium persulfate | cyclobutyl-BF ₃ K/ K ₂ S ₂ O ₈ | 2 x vials | 75 μmol / 100 μmol |
| Potassium cyclopentyltrifluoroborate / Potassium persulfate | cyclopentyl-BF ₃ K/ K ₂ S ₂ O ₈ | 2 x vials | 75 μmol / 100 μmol |
| Potassium cyclohexyltrifluoroborate / Potassium persulfate | cyclohexyl-BF ₃ K/ K ₂ S ₂ O ₈ | 2 x vials | 75 μmol / 100 μmol |
| Potassium ethyltrifluoroborate / Potassium persulfate | ethyl-BF ₃ K/ K ₂ S ₂ O ₈ | 2 x vials | 75 μmol / 100 μmol |
| Potassium isopropyltrifluoroborate / Potassium persulfate | isopropyl-BF ₃ K/ K ₂ S ₂ O ₈ | 2 x vials | 75 μmol / 100 μmol |
| Potassium methoxymethyltrifluoroborate / Potassium persulfate | MOM-BF ₃ K/ K ₂ S ₂ O ₈ | 2 x vials | 75 μmol / 100 μmol |
| <i>tert</i> -butyl peracetate | TBA | 2 x vials | 250 μmol |
| Ir[dF(CF ₃)ppy] ₂ (dtbbpy)[PF ₆] (Strem# 77-0425) / | Ir-cat-1 | 2 x vials | 8 μmol |

Visit www.strem.com for new product information and a searchable catalog.

Strem Chemicals, Inc.
7 Mulliken Way
Newburyport, MA 01950
U.S.A.
Tel: 978.499.1600
Fax: 978.465.3104
Email: info@strem.com

Strem Chemicals, Inc.
15, rue de l'Atome
Zone Industrielle
67800 BISCHHEIM France
Tel: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.
Postfach 1215
77672 KEHL
Germany
Tel: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.
An Independent Distributor of Strem Chemicals Products
Newton Hall, Town Street
Newton, Cambridge
England CB22 7ZE
Tel: +44 (0)1223 873 028
Fax: +44 (0)1223 870207
Email: enquiries@strem.co.uk

| Reagent Information | | | |
|---------------------|--|--------------|---------|
| Strem Item# | Vial | CAS | MW |
| 77-0425 | Ir[dF(CF ₃)ppy] ₂ (dtbbpy)][PF ₆] | 870987-63-6 | 1121.91 |
| --- | Potassium cyclopropyltrifluoroborate | 1065010-87-8 | 147.98 |
| --- | Potassium cyclobutyltrifluoroborate | 1065010-88-9 | 162.00 |
| --- | Potassium cyclopentyltrifluoroborate | 1040745-70-7 | 176.03 |
| --- | Potassium ethyltrifluoroborate | 44248-07-9 | 135.97 |
| --- | Potassium isopropyltrifluoroborate | 1041642-13-0 | 149.99 |
| --- | Potassium methoxymethyltrifluoroborate | 910251-11-5 | 151.97 |
| --- | <i>tert</i> -butyl peracetate solution TBPA (50% mineral spirits) | 107-71-1 | 132.16 |
| --- | Potassium persulfate | 7727-21-1 | 270.32 |

Typical Protocol

- The typical protocol is performed at 0.1 mol/l of substrate prepared as solution in DMSO or H₂O/ACN with 5 equiv. of trifluoroacetic acid. Each sealed reaction vial contains either 75 μmol of BF₃K reagent with 100 μmol of K₂S₂O₈ or 250 μmol of *tert*-butyl peracetate.
- Typically, the reaction is performed with 6200K white LED or blue LED (450 nm).
- The alkylation production kit contains 2 sets of vials allowing the screening of two different substrates or 1 substrate in two solvents.
- Sparging reaction solvents with nitrogen or argon while transferring reagents is important to achieve highest conversions of product. See protocol diagram for instructions.

Protocol at 500 μl volume reaction condition

- In the vial provided containing the Ir catalyst, prepare 4.0 ml of substrate solution at 0.1 mol/L of substrate (400 μmol substrate needed) in either DMSO or H₂O/ACN with 153 μL of trifluoroacetic acid (5 equiv.). Mix thoroughly.
- Degas the substrate/catalyst solution with subsurface sparging via N₂ or Ar line with exit needle for 5 minutes.
- Using a clean and dry syringe, add 500 μL of the substrate solution to each reaction vial (8 reactions).
- Place samples in vial holder **98-7650**.
- Turn on lamp (6200K white LED or blue LED 450 nm) and stir vials for 2 to 4 hours (or longer if necessary and there is substrate remaining). Be sure to plug in fan to maintain RT.
- Upon completion of reaction, prepare analytical sample for each reaction condition with 5 μl sample diluted into 200 μl in either DMSO or water/acetonitrile 50/50.
- Analyze resulting analytical samples by LC/MS.

Visit www.strem.com for new product information and a searchable catalog.


Strem Chemicals, Inc.
7 Mulliken Way
Newburyport, MA 01950
U.S.A
Tel: 978.499.1600
Fax: 978.465.3104
Email: info@strem.com

Strem Chemicals, Inc.
15, rue de l'Atome
Zone Industrielle
67800 BIRSCHHEIM France
Tel: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.
Postfach 1215
77672 KEHL
Germany
Tel: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.
An Independent Distributor of Strem Chemicals Products
Newton Hall, Town Street
Newton, Cambridge
England CB22 7ZE
Tel: +44 (0)1223 873 028
Fax: +44 (0)1223 870207
Email: enquiries@strem.co.uk

Protocol Diagram



Visit www.strem.com for new product information and a searchable catalog.

Strem Chemicals, Inc.
7 Mulliken Way
Newburyport, MA 01950
U.S.A
Tel: 978.499.1600
Fax: 978.465.3104
Email: info@strem.com

Strem Chemicals, Inc.
15, rue de l'Atome
Zone Industrielle
67800 BISCHEIM France
Tel: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.
Postfach 1215
77672 KEHL
Germany
Tel: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.
An Independent Distributor of Strem Chemicals Products
Newton Hall, Town Street
Newton, Cambridge
England CB22 7ZE
Tel: +44 (0)1223 873 028
Fax: +44 (0)1223 870207
Email: enquiries@strem.co.uk