



A Publication of Ascensus Specialties

THE STREM CHEMIKER

Vol. XXXII No.1

October, 2021



Air-Stable tris-(stilbene)Ni(0) Complexes

by Dr. Josep Cornella and Rakan Saeb

Table of Contents

Biographical Sketches:

Dr. Josep Cornella	1
Rakan Saeb	1

Article: Air-Stable tris-(stilbene)Ni(0) Complexes	2-14
---	-------------

Nickel Catalysts Related to the Article	15-20
--	--------------

New Products Since Chemiker XXXI (November 2019)	21-95
---	--------------

Ascensus Legacy Products	96-97
---------------------------------------	--------------

Ascensus Legacy Products - Coming Soon	98-101
---	---------------

cGMP & Custom Synthesis	102
--	------------

Available Booklets	103
---------------------------------	------------

American Chemical Society Award for Distinguished Service in the Advancement of Inorganic Chemistry	2021	Prof. Kenneth D. Karlin, Johns Hopkins University
	2020	Prof. Bruce E. Bursten, Worcester Polytechnic Institute
	2019	Prof. Debbie C. Crans, Colorado State University
	2018	Prof. Thomas B. Rauchfuss, University of Illinois
	2017	Prof. William B. Tolman, University of Minnesota

Canadian Society for Chemistry Award for Pure or Applied Inorganic Chemistry	2021	Laura Turculet, MCIC, Dalhousie University
	2020	Joe Gilroy, MCIC, University of Western Ontario
	2019	Prof. Tim Storr, Simon Fraser University
	2018	Prof. Eric Rivard, University of Alberta
	2017	Prof. Dwight Seferos, University of Toronto

USA Office

Strem
7 Mulliken Way
Newburyport, MA 01950
USA

Tel: (978) 499-1600
Fax: (978) 465-3104
Toll free numbers below US & Canada only
Tel: (800) 647-8736
Fax: (800) 517-8736
Email: info@ascensuspecialties.com

European Offices

Strem
15, rue de l'Atome
Zone Industrielle
F-67800 BISCHHEIM, France

Tel: +33 (0) 3 88 62 52 60
Fax: +33 (0) 3 88 62 26 81
Email: info-europe@ascensuspecialties.com

Strem Chemicals UK, Ltd.
An Independent Distributor of Strem Products
Newton Hall
Town Street, Newton
Cambridge, CB22 7ZE UK

Tel: +44 (0)1223 873 028
Fax: +44 (0)1223 870 207
Email: enquiries@strem.co.uk

Ascensus Specialties Global Headquarters

Ascensus Specialties
Global Headquarters
2821 Northup Way, Suite 275
Bellevue, WA 98004

Tel: (425) 448 1679 / TF
Tel: (888) 797 9362
Email: ask@ascensuspecialties.com
www.ascensuspecialties.com

©Copyright 2021 by



The Strem Chemiker
Vol. XXXII No.1
October, 2021

www.strem.com





IMPORTANT ANNOUNCEMENT!

We were acquired by Ascensus Specialties

Quote from Michael Strem, founder & former President of Strem Chemicals, "We thoughtfully chose Ascensus as the partner for the next phase of Strem's growth. Our employees' interests and well-being have been my top priority, and the strong cultural fit will provide opportunities for our team. The additional resources and global reach of Ascensus will deliver Strem's strong capabilities to an even broader audience of customers."

About Ascensus Specialties

For over 70 years, Ascensus Specialties has been one of the world's leading providers of sodium borohydride, specialty alcoholates, and boranes for a variety of industries.

Our specialty chemistries improve the way our clients create their products by making what they do easier, more efficient, safer, and more selective. From our world-class production plants in Elma, Washington, Evans City, Pennsylvania, and Newburyport, Massachusetts, we have a global reach that ensures our clients can consistently manufacture their products to the highest standards. What does all this mean? Our clients can make the best possible products to change people's lives for the better, and we are able to make the world a better place with chemistry made for what matters.

Our Expanded Product Lines

We offer an unparalleled range of high-quality specialty reagents that enable complex transformations so you can continue to create the products that matter. Our products are available in sample, intermediate, and bulk packages.

[Borohydrides](#) | [Boranes](#) | [Specialty Alcoholates](#) | [Boronic Acids & Esters](#) |
[Alkali Metals](#) | [APIs](#) | [Custom Synthesis](#) | [Catalysis](#) | [Thin Film Deposition](#)



CHEMISTRY MADE FOR WHAT MATTERS

Glossary of Terms

[α]_D	Specific rotation
AAS	Atomic Absorption Standard
ACS	Conforms to American Chemical Society specifications
air sensitive	Product may chemically react with atmospheric oxygen or carbon dioxide at ambient conditions. Handle and store under an inert atmosphere of nitrogen or argon.
amp	Ampouled
b.p.	Boiling point in °C at 760mm of Hg (Torr), unless otherwise noted
d.	Density
dec.	Decomposes
elec. gr.	Electronic Grade, suitable for electronic applications
f.p.	Flash point in °F
gran.	Granular
heat sensitive	Product may chemically degrade if stored for prolonged periods of time at ambient temperatures or higher. Store at 5°C or lower.
hydrate	Unspecified water content which may vary slightly from lot to lot
hygroscopic	Product may absorb water if exposed to the atmosphere for prolonged periods of time (dependent on humidity and temperature). Handle and store under an inert atmosphere of nitrogen or argon.
light sensitive	Product may chemically degrade if exposed to light
liq.	Liquid
m.p.	Melting point in °C
moisture sensitive	Product may chemically react with water. Handle and store under an inert atmosphere of nitrogen or argon.
NMR grade	Suitable as a Nuclear Magnetic Resonance reference standard
optical grade	For optical applications
pwdr.	Powder
primary standard	Used to prepare reference standards and standardize volumetric solutions
PURATREM	Product has a minimum purity of 99.99% (metals basis)
purified	A grade higher than technical, often used where there are no official standards
P. Vol.	Pore volume
pyrophoric	Product may spontaneously ignite if exposed to air at ambient conditions
reagent	High purity material, generally used in the laboratory for detecting, measuring, examining or analyzing other substances
REO	Rare Earth Oxides. Purity of a specific rare-earth metal expressed as a percentage of total rare-earths oxides.
SA	Surface area
store cold	Product should be stored at -18°C or 4°C, unless otherwise noted (see product details)
subl.	Sublimes
superconductor grade	A high purity, analyzed grade, suitable for preparing superconductors
tech. gr.	Technical grade for general industrial use
TLC	Suitable for Thin Layer Chromatography
v.p.	Vapor pressure mm of Hg (Torr)
xtl.	Crystalline

About Purity

Chemical purity	is reported after the chemical name, e.g. Ruthenium carbonyl, 99%
Metals purity	is reported in parentheses with the respective element, e.g. Gallium (III) bromide, anhydrous, granular (99.999%-Ga) PURATREM where 100% minus the metal purity is equal to the maximum allowable percentage of trace metal impurity

Biographical Sketches



Dr. Josep Cornella

Josep Cornella (Pep) was born in La Bisbal del Penedès, a small town in south Catalunya. He graduated in chemistry in 2008 from the University of Barcelona and carried MSc studies in the Department of Organic Chemistry studying the chemistry of allylboron reagents. After completing his masters thesis, he moved to the United Kingdom to pursue doctoral studies in the group of Prof. Igor Larrosa (QMUL). In early 2012, he earned his PhD working on the use of aromatic carboxylic acids as aryl donors in metal-catalyzed decarboxylative reactions. He then moved back to Catalunya, where he joined the group of Prof. Ruben Martin (ICIQ) with a COFUND Fellowship and later as Marie Curie Postdoctoral Fellow. There, he developed novel transformations involving Ni-catalyzed C–O bond activation and carbon dioxide insertion into organic molecules. In 2015, Pep obtained a Beatriu de Pinós Fellowship to carry out further postdoctoral studies in the group of Prof. Phil S. Baran at The Scripps Research Institute, California, USA. During this time, he worked on the discovery and implementation of new transformations based on the concept of “redox-active esters” as practical and readily available partners for Ni- and Fe-catalyzed C–C bond forming reactions. In summer 2017, he was selected as Max Planck Research Group Leader (MPRGL) to create and lead the Laboratory for Sustainable Catalysis at the Max-Planck-Institut für Kohlenforschung (Germany). His work focuses on the development of sustainable catalytic strategies based on the design of ligands, reagents and catalysts to streamline organic synthesis.



Rakan Saeb


Rakan grew up in Munich and studied Chemistry at the Ludwig Maximilian University of Munich (LMU Munich). He carried out undergraduate research in the group of Dr. Oliver Thorn-Seshold, working on the synthesis of novel topoisomerase inhibitors, as well as in the group of Prof. Paul Knochel, studying the chemistry of organolithium species. During his undergraduate studies Rakan worked in the biotech company 4SC AG, synthesizing potential drug candidates for cancer treatment. Driven by his interest in organometallic chemistry and catalysis, he joined Dr. Josep Cornella at the Max-Planck-Institut für Kohlenforschung, focusing on studies on Ni(0)-olefin complexes. In 2021 he graduated from LMU Munich, obtaining a M.Sc. in chemistry and is currently pursuing PhD studies with Dr. Cornella, further investigating the chemistry of Ni(0)-olefin complexes.

Air-Stable *tris*-(stilbene)Ni(0) Complexes¹

Dr. Josep Cornella and Rakan Saeb
Max-Planck-Institut für Kohlenforschung
Kaiser-Wilhelm-Platz, 1
45470 Mülheim an der Ruhr
Germany


Ni-olefin chemistry: a brief historical perspective

In 1897, Sabatier and Senderens reported on the gas-phase hydrogenolysis² and hydrogenation³ of substances aided by nickel (Ni) salts, a feature considered as the beginning of Ni catalysis. Despite the power of this seminal milestone, the application of Ni in the context of organic molecules remained underexplored for several years. In a remarkable example, Reppe reported on the tri- and tetramerization of acetylene (**1**) to generate benzene⁴ and COT (cyclooctatetraene, **2**),⁵ respectively. Up until then, the synthesis of tetraene **2** had been limited to the modification of *pseudopelletierine* (**3**), a naturally occurring compound obtained from the root-bark of the pomegranate tree. At that time, the synthesis of COT required a 13-step sequence, affording an overall yield of 1-2%.⁶ Yet, the application of Ni(CN)₂ on simple acetylene (**1**), led to the synthesis of COT in an outstanding 70% yield in one catalytic step (Scheme 1).^{5,7} This example illustrates the power of catalysis and highlights the paradigm shift in the construction of organic molecules.



Scheme 1: The power of catalysis: catalytic preparation of COT vs. traditional synthetic route.


Although the work of Reppe truly epitomizes the remarkable synergy between alkynes and Ni, the studies on the interactions between Ni and unsaturated compounds by Prof. Karl Ziegler and Prof. Günther Wilke laid the foundations of modern organonickel chemistry and industrial applications thereof. Yet, these contributions initiated with an unexpected and serendipitous observation during the study of the catalytic synthesis of higher molecular weight olefins from simple ethylene (**4**). The process under investigation at that time at the Max-Planck-Institut für Kohlenforschung was the “*Aufbaureaktion*” (addition reaction), where triethylaluminum (**5**) adds into terminal olefins to generate higher order alkylaluminum compounds (Scheme 2). At that time, this process received a lot of attention as it permitted the rapid synthesis of surfactants. Raising the reaction temperature leads to further addition of the alkylaluminium reagent into ethylene (“*Verdrängungsreaktion*”, displacement reaction) which, after elimination resulted in the formation of high order α -olefins (Scheme 2, path a).⁸ In a routine experiment in 1952, Erhard Holzcamp, a researcher working with Ziegler, did not obtain the expected polymerized olefins when heating 100 atm of ethylene with triethylaluminum to 100 °C. Instead, the reaction yielded mainly 1-butene (**6**). This unexpected outcome was systematically repeated by Holzcamp and led to a focused investigation into what were the critical factors that changed the outcome of the well-established “*Verdrängungsreaktion*”.⁹ Scientists concluded that traces of metallic particles, formed after exhaustive acidic washings of the reactor walls, were responsible for catalyzing the formation of 1-butene. This hypothesis led scientists to examine all elements of the periodic table to identify which element was responsible for the observed reactivity. Finally, it was found that Ni salts were the catalysts responsible for the formation of 1-butene, and this reactivity became known as the “*nickel effect*” (Scheme 2, path b).⁹ It is important to mention that the exhaustive scrutiny of a plethora of salts of different elements led to the discovery that titanium/aluminum catalysts efficiently catalyze the polymerization of olefins under mild conditions,⁹ an achievement that led Karl Ziegler together with Giulio Natta to be awarded with the Nobel Prize in Chemistry in 1963.¹⁰




Scheme 2: The "Aufbaureaktion", the "Verdrängungsreaktion", and the "nickel effect".

The discovery of the "nickel effect" opened up a new field of research and spurred a myriad of programs studying the interaction and reactivity of unsaturated hydrocarbons with nickel compounds.^{7, 11-13} Reed reported that Reppe catalysts based on $\text{Ni}(\text{CO})_4$ were able to catalyze the dimerization of butadiene (**7**) to generate COD (1,5-cyclooctadiene, **8**) in 30% – 40% yield.¹⁴ However, treatment of other Ni(II) salts with organometallic reagents would lead to the precipitation of metallic nickel with poor catalytic reactivity.¹⁵ Nevertheless, Wilke and co-workers studied further the "nickel effect" responsible for the abovementioned dimerization. Hence, when $\text{Ni}(\text{acac})_2$ was mixed with phenylacetylene (**9**) (as stabilizer) and reduced with organoaluminum compounds, 24% COD and 63% *t,t*-CDT (*trans,trans*-1,5,9-cyclododecatriene, **10**) were obtained (Scheme 3A).¹⁵⁻¹⁶ By altering the compound that stabilizes the reduced nickel, Wilke was capable of shifting the course of the reaction, thus favoring either the trimerization or the dimerization products (Scheme 3B).^{15, 17-18} This procedure was of high industrial relevance at that time, as CDT could be used as a starting material for the synthesis of dodecandioic acid and azacyclotridecan-2-one; precursors in the production of polyesters and super polyamides.¹⁹ To the best of our knowledge, this is one of the first examples in which, by changing the "stabilizer" that affects the Ni center, a completely different reaction outcome is obtained.^{15, 17-18} We currently understand this as the effect of the ligand on the Ni, which influences the energy of transition states, leading to different reaction pathways.

A. Wilke's reaction of butadiene with Ni and aluminium reagents




B. Wilke's reaction of butadiene with Ni and aluminium reagents with different ligands



Scheme 3: (A) Dimerization of butadiene by employing the "nickel effect". (B) By altering the stabilizer, *t,t*-CDT and 1,5-COD could be generated with high selectivity and excellent yields.^{15, 17-18}


Studies on the mechanism of the cyclooligomerization of butadiene led to the synthesis, isolation and characterization of the 16-electron complex (*trans,trans,trans*-cyclododeca-1,5,9-triene)nickel(0) (Ni(*t*-CDT), **13**), which represents an early example of *homoleptic metal-olefin complexes*.^{15-16, 20-21} Addition of 1,5-cyclooctadiene (COD) led to displacement of CDT from the nickel center and concomitant formation of the 18-electron bis(1,5-cyclooctadiene)nickel(0) (Ni(COD)₂, **14**). Both complexes turned out to be applicable in the catalytic trimerization of butadiene,¹⁵ and can be easily prepared from Ni(acac)₂ and AlR₃ (Scheme 4).^{16, 21}



Scheme 4: Synthesis of the first homoleptic Ni(0)-olefin complexes.^{15-16, 20-21}

The facile access to Ni(COD)₂ ([Strem 28-0010](#)) and Ni(*t*-CDT) laid the foundations for the development of other Ni(0)-olefin complexes and accelerated the research on organonickel chemistry. Since then, the 18-electron Ni(COD)₂ has reigned sovereign as a source of Ni(0) due to its relative stability compared to other Ni(0) precursors, and facile exchange of its COD ligands. Despite the vast utility, both, Ni(*t*-CDT) and Ni(COD)₂ suffer from instability towards air, thus decomposing rapidly to Ni black when exposed to traces of oxygen. The lability of the COD and *t,t,t*-CDT ligands provides great reactivity for ligand exchange; nevertheless, this property is tightly connected to thermal instability, requiring storage of these compounds under cryogenic conditions to avoid decomposition.²² Despite the great success of Ni(COD)₂ as a Ni(0) source, efforts have been placed in developing alternative Ni(0) sources to circumvent such practical caveats. Over the years, other 18-, 16-, and even 14-electron Ni(0)-olefin complexes were developed. The latter, represent extremely sensitive compounds which have been shown to be highly pyrophoric.¹¹⁻¹² A solution to their stability is the addition of external σ -donating phosphines¹² or NHC ligands, thus turning an unstable 14-electron Ni(0)-olefin complexes into stable heteroleptic 16-electron complexes.²³ This strategy was recently applied by Montgomery and Hazari, showing that such complexes can be applied as synthons for other bis-olefin-Ni(0) complexes as well as some air-stable pre-catalysts.²⁴⁻²⁷ In addition, Ni(II) pre-catalysts bearing various σ -donating ligands,²⁸⁻³² NHC-Ni(0)-olefin complexes,²⁴⁻²⁷ as well as encapsulation of Ni(COD)₂ with paraffin³³ have been developed and represent only a few examples of providing practical access to Ni(0) for catalysis. Some of these pre-catalysts are indeed extremely useful and are now commercially available in the Strem catalog. While being extremely valuable solutions for specific problems, different issues associated with such complexes still exist. For example, these complexes usually lack modularity, and thus restrict extensive ligand survey. Additionally, the necessity of reducing Ni(II) complexes *via* ligand exchange/reduction/disproportionation pathways in order to release a Ni(0) species can result in undefined ligand to metal ratios and formation of undesired salt byproducts.³⁴ Therefore, a modular Ni(0) source that enables rapid and broad ligand screening and would not suffer from the issues connected to Ni(COD)₂ was still highly coveted.

A fundamental discovery in the area of Ni(0)-olefin complexes was reported by Steigerwald and Nuckolls. While searching for a stable tetraalkyl Ni(IV) complex, the authors synthesized a unique 16-electron Ni(0)-olefin complex bearing three (5*Z*,11*E*)-dibenzo[*a,e*][8]annulene molecules as ligands (**15**). The authors obtained crystallographic evidence for its connectivity and pointed out a rather exceptional stability towards oxidation under air; however, no catalytic application was reported (Figure 1).³⁵



15


Nuckolls' and Steigerwalds complex

Figure 1: Air-stable 16-electron Ni(0)-olefin complex.³⁵

Despite this isolated precedent, Ni(COD)₂ still dominates the field of organonickel chemistry for 60 years as Ni(0) precursor. Recently, Schrauzer's air-stable Ni(0)-olefin complex Ni(COD)(DQ) ([Strem 28-0040](#)) (DQ = duroquinone, **16**)³⁶ was reported by Engle to be catalytically active in certain catalytic domains where a strongly nucleophilic ligand is present.³⁷

Ni(^{4-CF₃}stb)₃ – 1st generation of air-stable Ni(0)-olefin pre-catalysts

During our investigations on the mechanism of a low-temperature Ni-olefin-Li-catalyzed Kumada-Corriu cross-coupling, the preparation of various 16-electron Ni(0)-olefin complexes was of great interest.³⁸⁻³⁹ We turned our attention to a rather forgotten complex reported by Wilke in a patent: Ni(*t*-stb)₃ (tris-(*trans*-stilbene)Ni(0), **17**). This complex was originally prepared by an olefin exchange reaction, using Ni(*t*-CDT) and *trans*-stilbene.²¹ Similar to Steigerwald and Nuckolls' complex (**15**),³⁵ **17** exhibited certain stability when exposed to air (6 hours at 25 °C or 3 days at 0 °C without decomposition). Yet, Ni(*t*-stb)₃ was temperature sensitive and slow ligand dissociation eventually led to decomposition (<1 week).²² When trifluoromethyl groups were introduced at the *para*-position of the stilbene (*trans*-1,2-bis(4-trifluoromethyl)phenyl)ethene, ^{4-CF₃}stb, **18**), Ni(^{4-CF₃}stb)₃ (**19**) was obtained in 70% yield (Scheme 5, left).²² Surprisingly, Ni(^{4-CF₃}stb)₃ exhibited a remarkable air-stability for several months when kept in a freezer at -18 °C, and could be manipulated at ambient temperature on a benchtop setting for several days before decomposition was observed.²²



Scheme 5: Syntheses of **19**.

X-ray analysis revealed that, similarly to **15** and **13**, Ni(^{4-CF₃}stb)₃ consists of three olefin ligands wrapped around the Ni center in a propeller-like arrangement, thus rendering a slightly distorted trigonal-planar coordination geometry (Figure 2). Donation of electron density from the Ni center towards the olefin π* orbitals results in an elongation of the stilbene olefin bond lengths (average of 1.389 Å), which coincides with Nuckolls' and Steigerwald's complex (**15**) (average of 1.39 Å)³⁵ and with Ni(COD)₂ (average of 1.39 Å).^{22,40} It is noteworthy to point out that Ni(^{4-CF₃}stb)₃ is chiral, resulting in a racemate in the unit cell. With the goal of facilitating the synthesis of **19**, a practical method was developed utilizing the traditional Wilke recipe: reduction of Ni(acac)₂ with AlEt₃ in the presence of **18** generated 23 g of the complex in high yields in a single run (95%, Scheme 5, right). Besides the X-ray analysis, additional information on the 1:3 metal:ligand ratio of **19** was obtained by elemental analysis, coordination with external ligands as well as by degradation studies. However, no conclusive NMR analysis was obtained for this complex and thus, it could not be discarded that the powder form of the complex consists of alternative coordination modes around the Ni center. Studies toward fully elucidating the structural features of **19** are currently ongoing in our laboratory.²²





Figure 2: ORTEP diagram of complex **19**. Ar = $p\text{-CF}_3\text{-C}_6\text{H}_4$. Crystallized from toluene. Two molecules of toluene are present in the unit cell (omitted for clarity). Colors: Green: Ni; Black: C; Yellow: F. Hydrogen atoms were omitted for clarity. Left: side view. Right: top view along the principal axis.

Ligand exchange studies

In order to evaluate the potential of $\text{Ni}(\text{}^{4\text{-CF}_3\text{stb}})_3$ as a viable Ni(0) source for catalysis, its ligand exchange abilities were explored. Hence, $\text{Ni}(\text{}^{4\text{-CF}_3\text{stb}})_3$ was mixed with a set of different ligands commonly applied in nickel catalysis, namely phosphines and diamines. Interestingly, complex **19** underwent smooth ligand exchange with bidentate phosphines such as dppe (1,1'-bis-(diphenylphosphino)ferrocene, **20**),⁴¹⁻⁴² monodentate triphenyl phosphine (PPh_3 , **21**) and amines such as 2,2'-bipyridine (bipy, **22**) (Scheme 6).⁴³⁻⁴⁴ In all cases, 16-electron complexes were formed bearing one $^{4\text{-CF}_3\text{stb}}$ unit. The products formed were characterized by crystallographic and spectroscopic techniques, and identified as (dppe) $\text{Ni}(\text{}^{4\text{-CF}_3\text{stb}})$ (**23**), $(\text{PPh}_3)_2\text{Ni}(\text{}^{4\text{-CF}_3\text{stb}})$ (**24**) and $(\text{bipy})\text{Ni}(\text{}^{4\text{-CF}_3\text{stb}})$ (**25**) respectively.²² The latter example exhibits similar features as Nuckolls' and Steigerwald's complex, when two ligands in **15** are displaced by bipy.³⁵




Scheme 6: Ligand exchange reactions investigated for $\text{Ni}(\text{}^{4\text{-CF}_3\text{stb}})_3$ with dppe, PPh_3 , bipy and THF.

As shown in Scheme 6, ligand exchange with complex **19** is highly facile for a variety of mono- and bidentate ligands. Yet, such facile ligand exchange turned out to be its Achilles heel: when dissolved in THF at room temperature in the absence of any supporting ligand such as diamine, phosphine or NHC, **19** would rapidly decompose into Ni black. Attempts to intercept intermediates of this decomposition resulted in the isolation of the extremely sensitive compound $\text{Ni}(\text{}^{4\text{-CF}_3\text{stb}})_2(\text{THF})$ (**26**) at a very low temperature, which corresponds to the substitution of one $^{4\text{-CF}_3\text{stb}}$ for a THF molecule. Although still speculative, formation of **26** suggests that decomposition of **19** in THF solutions proceeds via ligand dissociation facilitated by the nucleophilic ethereal solvent.²² Decomposition of tris-Ni(0)-olefin complexes in ethereal solvents has been noted in the case of the 16-electron complex $(\mu\eta^2, \eta^2\text{C}_7\text{H}_{12})_2[\text{Ni}(\eta^2, \eta^2\text{C}_7\text{H}_{12})_2]_2$. While no evidence was provided, ligand displacement by THF at 20 °C was believed to be the reason for decomposition.⁴⁵ Due to that rapid exchange of ligands in solution, complete characterization of the $\text{Ni}(\text{}^{4\text{-CF}_3\text{stb}})_3$ in solution still poses severe challenges.²²

Catalytic transformations


The facile ligand exchange in **19** observed in solution immediately revealed that this behavior should be of high value in catalytic contexts. Moreover, complex **19** has been shown to deliver well-defined Ni(0):L ratios, thus permitting a control over the species at the onset of the reaction. In order to explore the catalytic abilities of **19**, a survey of important catalytic transformations, which utilize Ni as catalyst was performed. The results were benchmarked with the yields reported using either Ni(COD)₂ or an alternative Ni(II) precursor.

Hartwig and co-workers reported that Ni(COD)₂ was inefficient in the formation of heterobiaryls, making the application of a Ni(II) species essential for a successful coupling.²⁸ Ni(4-CF₃stb)₃ turned out to be a good pre-catalyst in this catalytic transformation, resulting in the generation of high yields (>99%) of the coupling product (**27**). When Ni(4-CF₃stb)₃ was stored for more than 180 days in a freezer under air, an 82% yield was obtained (Scheme 7A).²² Good yields were also observed in the aryl coupling with acetals developed by Doyle (Scheme 7B).^{22, 46} Furthermore, **19** demonstrated to be also applicable in C(sp²)-C(sp²) Negishi cross-coupling reactions between an aryl bromide and a vinyl zinc reagent (Scheme 7C).^{22, 47} In the context of alkyl-alkyl Negishi coupling, **19** also performed well although higher temperatures were required (Scheme 7D).^{22, 48} The requirement of higher temperature might indicate that differences between Ni(COD)₂ and Ni(4-CF₃stb)₃ might exist in generating the catalytically active species (*vide infra*). Heck coupling also were amenable and **19** afforded good yields of the C-C product **31** (Scheme 7E).^{22, 49}



Scheme 7: C–C bond forming reactions catalyzed by **19**. (A) Heteroaryl Suzuki cross-coupling. (B) C–O arylation of acetals. (C) Aryl-vinyl Negishi cross-coupling. (D) Alkyl-Alkyl Negishi cross-coupling. (E) A Heck-type reaction. ^aThis reaction was performed with a batch of **19** that was stored in a freezer at –18 °C under air for >180 d. ^bThe yield was determined by gas chromatography with a flame ionization detector using an internal standard. ^c25 °C. ^dUsing a catalyst stored under air at 25 °C for 14 d.


The cycloaddition of benzonitrile with 2,3-dimethylbuta-1,3-diene, using catalytic amounts of PCy_3 and $\text{Ni}(\text{COD})_2$, was reported to generate the pyridine derivative **32** in 88% yield.⁵⁰ When $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$ was employed, a yield of 84% was achieved. Additionally, when displacing the air-sensitive ligand by its protonated air-stable counterpart $\text{PCy}_3\bullet\text{HBF}_4$ together with amounts of K_3PO_4 , equimolar to the ligand, an only slightly reduced yield (77%) was obtained, permitting the performance of this catalytic transformation without the need of a glovebox (Scheme 8A).²²



Scheme 8: Additional examples of $\text{Ni}(0)$ -catalysis. (A) Oxidative cyclization. (B) C–H activation. (C) Buchwald-Hartwig C–N bond formation with alkylamines. (D) C–SMe reduction with silanes. ^aThis reaction was performed using PCy_3 . ^bThe yield was determined by gas chromatography with a flame ionization detector using an internal standard.

While first examples of $\text{C}(\text{sp}^2)\text{--H}$ bond activation were already reported in 1963,⁵¹ it was not until recently that Ni-catalyzed C–H activation became a routinely explored strategy. Following a recent protocol by Chatani,⁵² isoquinolone derivative **33** (Scheme 8B) could be generated in an excellent yield of 94% (86% with $\text{Ni}(\text{COD})_2$).^{22, 52} Another cross-coupling reaction of paramount importance is the Buchwald-Hartwig amination reaction to forge $\text{C}(\text{sp}^2)\text{--N}$ bonds.⁵³ Amination of 4-chlorobenzotrifluoride with morpholine proceeded smoothly with the aid of SiPr as the ligand for Ni (91%, Scheme 8C, **34**).^{22, 32} $\text{Ni}(0)$ -olefin complexes were demonstrated to be a valuable source for the formation of heterogeneous $\text{Ni}(0)$ nanoparticles by Martin and co-workers. Reduction of 2-(methylthio)naphthalene with silanes afforded 91% yield of **35** (Scheme 8D).^{22, 54}

Finally, we demonstrated that $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$ was able to have orthogonal reactivity to that of $\text{Ni}(\text{COD})_2$. For example, in the C–H activation of electron-poor arenes, the use of $\text{Ni}(\text{COD})_2$ results in the formation of catalytically inactive species as a result of hydrometalation of the COD. Yet, the electron-deficient nature of the double bond in $\text{}^4\text{-CF}_3\text{stb}$ prevents such hydrometalation in the ligand, and permits the hydrometalation to occur in the triple bond, thus generating **36** in excellent yields (Scheme 9).^{22, 26}



Scheme 9: $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$ prevents detrimental deactivation pathways associated to the COD in $\text{Ni}(\text{COD})_2$. ^aYield determined by ^{19}F NMR.

Shortcomings of the $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$ pre-catalyst

Although $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$ is applicable to a variety of Ni-catalyzed organic transformations, several limitations still prevailed and required further detailed analysis. For example, $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$ is unstable at high temperatures when heated in the solid state. Moreover, the complex decomposes within a week when exposed to air at room temperature. More importantly, $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$ is also unstable in solution without supporting ligands, which prevents the preparation of stock solutions or its use in certain homogenous "ligand-less" reactions. Finally, the electron-deficient nature of the $\text{}^4\text{-CF}_3\text{stb}$ ligand affects elementary steps in certain catalytic domains, thus requiring harsher conditions (Figure 3).²²





Figure 3: Drawbacks of the air-stable $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$.

These limitations have consequences when attempting certain catalytic reactions (Scheme 10). For example, the Buchwald-Hartwig amination using anilines required slightly elevated reaction temperatures (130 °C) to generate **37** in 90% yield (Scheme 10A).⁵⁵ When performing the reaction at 100 °C, as reported for $\text{Ni}(\text{COD})_2$, traces of product were observed.²² A similar behavior was observed in the Ni-catalyzed ester formation from amides reported by Garg. While ester **38** was generated in a yield of 85% when using $\text{Ni}(\text{COD})_2$ at 100 °C,⁵⁶ only 39% was obtained when using $\text{Ni}(\text{}^4\text{-CF}_3\text{stb})_3$. Heating up the reaction to 130 °C led to an increase in yield up to 65% (Scheme 10B).²² We believe that in these cases, the high energetic cost to dissociate the stilbene ligand from L-Ni(0)-stilbene complexes might affect the energetic profile of the reaction (*vide infra*).²²


A.



Ni(COD)₂: 96%
 Ni(4-CF₃stb)₃: <5% (90% at 130 °C)

37


B.



Ni(COD)₂: 85%
 Ni(4-CF₃stb)₃: 39% (65% at 130 °C)

38


C.



Ni(COD)₂: 92%
 Ni(4-CF₃stb)₃: <5%

39

D.



Ni(COD)₂: 70%
 Ni(4-CF₃stb)₃: <5%


40

Scheme 10: Limitations in the catalytic reactivity of Ni(4-CF₃stb)₃. (A) Buchwald-Hartwig C–N bond formation with arylamines. (B) Ester formation through C–N bond activation of amides. (C) Alkyl-alkyl Kumada cross-coupling. (D) Conversion of vinyl triflates into vinyl iodides.

Examples of the instability of Ni(4-CF₃stb)₃ in solution can be found when attempting the Kambe coupling⁵⁷ or the Reisman triflate-iodide exchange.⁵⁸ In both cases Ni(4-CF₃stb)₃ does not survive the time required to generate the active species, thus decomposing without catalyzing the reaction (Scheme 10C and D).²²

Ni(4-tBu^{stb})₃ – 2nd generation of air-stable pre-catalysts

In order to tackle the problematics associated with Ni(4-CF₃stb)₃, various Ni(0) stilbene complexes bearing different substituents in the stilbene unit were synthesized, isolated and characterized.⁵⁹ Eventually, solid-state analysis led us to hypothesize that in addition to its electronic influence, the CF₃-moiety also exerted a severe steric constraint at the edges of the complexes. We then envisaged that in order to discriminate between electronics and steric effects in the stability toward oxidation, the CF₃ group could be replaced by the isostere ^tBu group. Indeed, a complex bearing three *p*-Bu-stilbene ligands (**41**) was successfully synthesized following the optimized protocol for **19** (Scheme 11A). This orange complex (**42**, Ni(4-^tBustb)₃) was found to be remarkably air stable when compared to all the previously synthesized complexes, including **19**. A direct comparison of the physical properties between **19** and **42** is highlighted in Scheme 11B. It is clear that **42** exhibits both an increased stability in solid state as well as in solution.⁵⁹ While **19** suffers from decomposition in solution at ambient temperature when no external ligand is added, **42** does not suffer from these drawbacks, and the Ni species remain stable in solution in various solvents. Although both complexes are stable for at least 3 months in the freezer (–18 °C), **42** did not show signs of decomposition for one month at room temperature.⁵⁹ Once again, the composition of this complex as well as its structure was determined by elemental analysis, X-ray crystallography, XPS, titration and degradation studies. Yet, complete structural analysis by NMR was prevented due to the exchange of the ligands in solution.⁵⁹



Scheme 11: (A) Preparation of $\text{Ni}(\text{4-tBu}\text{stb})_3$. (B) Comparison between **19** and **42** and ORTEP structures of these complexes compared to each other. Colors: Green: Ni, Black: C, Yellow: F. Hydrogen atoms are omitted for clarity.

Catalytic transformations

The different stabilities in solution between **19** and **42** suggested that different ligand exchange rates occur, which could dramatically affect the catalytic properties of each complex. An example of such differences was noted in the Buchwald-Hartwig amination reaction with anilines (*vide supra*).⁵⁵ When the reaction profile was monitored by NMR, severe differences were encountered. While complex **19** was inefficient at the reported temperature (100°C), the newly prepared complex **42** catalyzed the reaction smoothly (Figure 4).⁵⁹





Figure 4: Kinetic profiles of the consumption of the aryl chloride in the Buchwald-Hartwig amination.

The differences in reactivity observed, led us to hypothesize that $\text{Ni}(\text{4-}^{\text{tBu}}\text{stb})_3$ could be amenable in transformations where $\text{Ni}(\text{4-}^{\text{CF}_3}\text{stb})_3$ failed to react. Therefore, reactions in Scheme 10 were repeated using $\text{Ni}(\text{4-}^{\text{tBu}}\text{stb})_3$ instead (Scheme 12). As shown, the Buchwald-Hartwig amination (Scheme 12A),^{55, 59} the activation of amides (Scheme 12B).^{56, 59} Kambe's alkyl-alkyl cross-coupling using isoprene as ligand (Scheme 12C),⁵⁷ as well as the catalytic conversion of vinyl triflates into vinyl iodides (Scheme 12D),⁵⁸ recovered the reactivity when using $\text{Ni}(\text{4-}^{\text{tBu}}\text{stb})_3$, affording similar yields as seen with $\text{Ni}(\text{COD})_2$.⁵⁹



Scheme 12: Catalytic transformations performed. (A) Buchwald-Hartwig C–N bond formation with arylamines. (B) Ester formation through C–N bond activation of amides. (C) Alkyl-alkyl Kumada cross-coupling. (D) Conversion of vinyl triflates into vinyl iodides.

The differences in reactivity between different stilbene ligands are certainly of great interest, as tuning the electronics on the olefin should now be considered as an additional parameter when optimizing Ni-catalyzed transformations. Investigations on the physical properties, catalytic activity, stereoelectronic effects and structural differences rendered by the simple change of substituents in the stilbene unit in $\text{Ni}(\text{O})$ complexes are currently ongoing in our laboratory.

Summary

In summary, complexes $\text{Ni}(\text{4-}^{\text{CF}_3}\text{stb})_3$ (Strem 28-0060) and $\text{Ni}(\text{4-}^{\text{tBu}}\text{stb})_3$ (Strem 28-0070) represent more air-stable alternatives to the current $\text{Ni}(\text{O})$ sources, which allows practical incorporation of $\text{Ni}(\text{O})$ into routine experimentation. These two complexes permit the practitioner to perform a plethora of various catalytic transformations commonly applied in organic synthesis without the need of using gloveboxes or complex Schlenk techniques. It has been shown that $\text{Ni}(\text{4-}^{\text{CF}_3}\text{stb})_3$ can be as catalytically competent as $\text{Ni}(\text{COD})_2$, and in some instances, surpasses its catalytic activity. The instability of $\text{Ni}(\text{4-}^{\text{CF}_3}\text{stb})_3$ in solution without the addition of external ligands together with the decomposition after prolonged exposure to air, led to the design and synthesis of $\text{Ni}(\text{4-}^{\text{tBu}}\text{stb})_3$. This last complex represents a robust source of $\text{Ni}(\text{O})$, with broader scope of reactivity which approaches the levels of $\text{Ni}(\text{COD})_2$.

References

1. This article has been elaborated utilizing parts of the masters thesis of Mr. Rakan Saeb (Ludwig Maximilian University of Munich, Germany, **2021**), Doctoral Dissertation from Dr. Lukas Nattmann (Ruhr Universität Bochum, Germany, **2020**) as well as the two articles in references 22 and 59.
2. Sabatier, P.; Senderens, J. B. Action du nickel sur l'ethylene. *C. r. hebd. séances Acad. sci.* **1897**, *124*, 616-618.
3. Sabatier, P.; Senderens, J. B. Action du nickel sur l'ethylene. Synthèse de l'ethane. *C. r. hebd. séances Acad. sci.* **1897**, *124*, 1358-1361.
4. Reppe, W.; Schweckendiek, W. Cyclisierende Polymerisation von Acetylen. III Benzol, Benzolderivate und hydroaromatische Verbindungen. *Justus Liebig's Ann. Chem.* **1948**, *560*, 104-116.
5. Reppe, W.; Schlichting, O.; Klager, K.; Toepel, T. Cyclisierende Polymerisation von Acetylen I Über Cyclooctatetraen. *Justus Liebig's Ann. Chem.* **1948**, *560*, 1-92.
6. Willstätter, R.; Waser, E. Über Cyclo-octatetraen. *Ber. Dtsch. Chem. Ges.* **1911**, *44*, 3423-3445.
7. Jolly, P. W.; Wilke, G. The Organic Chemistry of Nickel - Organic Synthesis. In *Organometallic Chemistry*; Maitlis, P. M.; Stone, F. G. A.; West, R., Eds.; Academic Press: New York **1975**; Vol. 2.
8. Ziegler, K. Aluminium-organische Synthese im Bereich olefinischer Kohlenwasserstoffe. *Angew. Chem.* **1952**, 323-350.
9. Ziegler, K.; Holzkamp, E.; Breil, H.; Martin, H., Das Mülheimer Normaldruck-Polyäthylen-Verfahren. *Angew. Chem.* **1955**, 541-547.
10. Wilke, G. Fifty Years of Ziegler Catalysts: Consequences and Development of an Invention. *Angew. Chem., Int. Ed.* **2003**, *42*, 5000-5008.
11. Schrauzer, G. N. Some Advances in the Organometallic Chemistry of Nickel. In *Advances in Organometallic Chemistry*; Stone, F. G. A.; West, R., Eds.; Academic Press: New York **1965**; Vol. 2, pp 1-48.
12. Wilke, G.; Jolly, P. W. The Organic Chemistry of Nickel - Organonickel Complexes. In *Organometallic Chemistry*; Maitlis, P. M.; Stone, F. G. A.; West, R., Eds.; Academic Press: New York; **1974**; Vol. 1.
13. Tamaru, Y., *Modern Organonickel Chemistry*; WILEY-VCH Verlag GmbH & Co. KGaA: Weinheim, **2005**.
14. Reed, H. W. B. The catalytic cyclic polymerisation of butadiene. *J. Chem. Soc.* **1954**, 1931-1941.
15. Wilke, G.; Bogdanović, B.; Heimbach, P.; Kröner, M.; Müller, E. W. Reactions of New Pi-Complex Catalysts. In *POLYMERIZATION AND POLYCONDENSATION PROCESSES*; Platzner, N. A. J.; Eds.; American Chemical Society: Washington, DC **1962**; Vol. 34, pp 137-144.
16. Bogdanović, B.; Kröner, M.; Wilke, G., Übergangsmetallkomplexe, I. Olefin-Komplexe des Nickels(0). *Justus Liebig's Ann. Chem.* **1966**, *699*, 1-23.
17. Wilke, G.; Müller, E. W. Verfahren zur katalytischen Di- bzw. Trimerisation von 1,3-Diolefinen. DE1140569B, **1962**.
18. Brenner, W.; Heimbach, P.; Hey, H.; Müller, E. W.; Wilke, G., Über die katalytische Umwandlung von Olefinen, III. Synthese von *cis,cis*-Cyclooctadien-(1.5) und *cis*-1.2-Divinyl-cyclobutan. *Justus Liebig's Ann. Chem.* **1969**, *727*, 161-182.
19. Wilke, G. Tagung der Nordwestdeutschen Chemiedozenten in Münster/Westf. *Angew. Chem.* **1957**, 392-400.
20. Wilke, G. Hauptversammlung der Gesellschaft Deutscher Chemiker. *Angew. Chem.* **1960**, *72*, 565-593.
21. Wilke, G.; Müller, E. W.; Kröner, M.; Heimbach, P.; Breil, H. Verfahren zur Herstellung von CO- und NO-freien Komplexverbindungen der Übergangsmetalle. DE000001191375, **1965**.
22. Nattmann, L.; Saeb, R.; Nöthling, N.; Cornella, J. An air-stable binary Ni(0)-olefin catalyst. *Nat. Catal.* **2020**, *3*, 6-13.
23. Haynes II, M. T.; Jackson, E. P.; Montgomery, J. Nickel Complexes of *N*-Heterocyclic Carbenes. In *N-Heterocyclic Carbenes*; Nolan, S. P.; Eds.; Wiley-VCH: Weinheim; **2014**; pp 371-396.
24. Todd, D. P.; Thompson, B. B.; Nett, A. J.; Montgomery, J. Deoxygenative C-C Bond-Forming Processes via a Net Four-Electron Reductive Coupling. *J. Am. Chem. Soc.* **2015**, *137*, 12788-12791.
25. Nett, A. J.; Cañellas, S.; Higuchi, Y.; Robo, M. T.; Kochkodan, J. M.; Haynes, M. T.; Kampf, J. W.; Montgomery, J. Stable, Well-Defined Nickel(0) Catalysts for Catalytic C-C and C-N Bond Formation. *ACS Catal.* **2018**, *8*, 6606-6611.
26. Nett, A. J.; Zhao, W.; Zimmerman, P. M.; Montgomery, J. Highly Active Nickel Catalysts for C-H Functionalization Identified through Analysis of Off-Cycle Intermediates. *J. Am. Chem. Soc.* **2015**, *137*, 7636-7639.
27. Wu, J.; Faller, J. W.; Hazari, N.; Schmeier, T. J. Stoichiometric and Catalytic Reactions of Thermally Stable Nickel(0) NHC Complexes. *Organometallics* **2012**, *31*, 806-809.
28. Ge, S.; Hartwig, J. F. Highly Reactive, Single-Component Nickel Catalyst Precursor for Suzuki-Miyaura Cross-Coupling of Heteroaryl Boronic Acids with Heteroaryl Halides. *Angew. Chem., Int. Ed.* **2012**, *51*, 12837-12841.
29. Park, N. H.; Teverovskiy, G.; Buchwald, S. L. Development of an Air-Stable Nickel Precatalyst for the Amination of Aryl Chlorides, Sulfamates, Mesylates, and Triflates. *Org. Lett.* **2014**, *16*, 220-223.
30. Standley, E. A.; Jamison, T. F. Simplifying Nickel(0) Catalysis: An Air-Stable Nickel Precatalyst for the Internally Selective Benzoylation of Terminal Alkenes. *J. Am. Chem. Soc.* **2013**, *135*, 1585-1592.
31. Shields, J. D.; Gray, E. E.; Doyle, A. G. A Modular, Air-Stable Nickel Precatalyst. *Org. Lett.* **2015**, *17*, 2166-2169.
32. Magano, J.; Monfette, S. Development of an Air-Stable, Broadly Applicable Nickel Source for Nickel-Catalyzed Cross-Coupling. *ACS Catal.* **2015**, *5*, 3120-3123.
33. Dander, J. E.; Weires, N. A.; Garg, N. K. Benchtop Delivery of Ni(COD)₂ using Paraffin Capsules. *Org. Lett.* **2016**, *18*, 3934-3936.
34. Beromi, M. M.; Banerjee, G.; Brudvig, G. W.; Charboneau, D. J.; Hazari, N.; Lant, H. M. C.; Mercado, B. Q. Modifications to the Aryl Group of dppe-Ligated Ni σ -Aryl Precatalysts: Impact on Speciation and Catalytic Activity in Suzuki-Miyaura Coupling Reactions. *Organometallics* **2018**, *37*, 3943-3955.


References (continued)

35. Carnes, M.; Buccella, D.; Chen, J. Y.-C.; Ramirez, A. P.; Turro, N. J.; Nuckolls, C.; Steigerwald, M. A Stable Tetraalkyl Complex of Nickel(IV). *Angew. Chem., Int. Ed.* **2009**, *48*, 290-294.
36. Schrauzer, G. N.; Thyret, H. Neuartige „Sandwich“-Verbindungen des Nickel(0). Zur Kenntnis von Durochinon-Nickel(0)-Komplexen mit cyclischen Dienen. *Z. Naturforsch. B* **1962**, *17*, 73-76.
37. Tran, V. T.; Li, Z.-Q.; Apolinar, O.; Derosa, J.; Joannou, M. V.; Wisniewski, S. R.; Eastgate, M. D.; Engle, K. M. Ni(COD)(DQ): An Air-Stable 18-Electron Nickel(0)–Olefin Precatalyst. *Angew. Chem., Int. Ed.* **2020**, *59*, 7409-7413.
38. Nattmann, L.; Lutz, S.; Ortsack, P.; Goddard, R.; Cornella, J. A Highly Reduced Ni–Li–Olefin Complex for Catalytic Kumada–Corriu Cross-Couplings. *J. Am. Chem. Soc.* **2018**, *140*, 13628-13633.
39. Lutz, S.; Nattmann, L.; Nöthling, N.; Cornella, J. The Role of 16-electron Nickel(0)–Olefin Complexes in Low-Temperature C(sp²)–C(sp³) Kumada Cross-Couplings. *Organometallics*, Article ASAP.
40. Diebks, H.; Dietrich, H. Die Kristallstruktur von bis-Cyclooctadien-(1,5)-Nickel(0) Z. *Kristallogr. – Cryst. Mater.* **1965**, *122*, 1-23.
41. Yin, G.; Kalvet, I.; Englert, U.; Schoenebeck, F. Fundamental Studies and Development of Nickel-Catalyzed Trifluoromethylthiolation of Aryl Chlorides: Active Catalytic Species and Key Roles of Ligand and Traceless MeCN Additive Revealed. *J. Am. Chem. Soc.* **2015**, *137*, 4164-4172.
42. Birkholz, M.-N.; Freixa, Z.; van Leeuwen, P. W. N. M. Bite angle effects of diphosphines in C–C and C–X bond forming cross coupling reactions. *Chem. Soc. Rev.* **2009**, *38*, 1099-1118.
43. Tasker, S. Z.; Standley, E. A.; Jamison, T. F. Recent advances in homogeneous nickel catalysis. *Nature* **2014**, *509*, 299-309.
44. Choi, J.; Fu, G. C. Transition metal–catalyzed alkyl–alkyl bond formation: Another dimension in cross-coupling chemistry. *Science* **2017**, *356*, eaaf7230.
45. Proft, B.; Pörschke, K.-R.; Lutz, F.; Krüger, C. 1,6-Heptadien-Nickel(0)-Komplexe: rac/meso-(μ-η²,η²-C₇{Ni(η²-C₇H₁₂)₂ und L-Ni(η²,η²-C₇H₁₂). *Chem. Ber.* **1991**, *124*, 2667-2675.
46. Graham, T. J. A.; Doyle, A. G. Nickel-Catalyzed Cross-Coupling of Chromene Acetals and Boronic Acids. *Org. Lett.* **2012**, *14*, 1616-1619.
47. Staudaher, N. D.; Stolley, R. M.; Louie, J. Synthesis, mechanism of formation, and catalytic activity of Xantphos nickel π-complexes. *Chem. Commun.* **2014**, *50*, 15577-15580.
48. Zhou, J.; Fu, G. C. Cross-Couplings of Unactivated Secondary Alkyl Halides: Room-Temperature Nickel-Catalyzed Negishi Reactions of Alkyl Bromides and Iodides. *J. Am. Chem. Soc.* **2003**, *125*, 14726-14727.
49. Matsubara, R.; Gutierrez, A. C.; Jamison, T. F. Nickel-Catalyzed Heck-Type Reactions of Benzyl Chlorides and Simple Olefins. *J. Am. Chem. Soc.* **2011**, *133*, 19020-19023.
50. Ohashi, M.; Takeda, I.; Ikawa, M.; Ogoshi, S. Nickel-Catalyzed Dehydrogenative [4 + 2] Cycloaddition of 1,3-Dienes with Nitriles. *J. Am. Chem. Soc.* **2011**, *133*, 18018-18021.
51. Kleiman, J. P.; Dubeck, M. The Preparation of Cyclopentadienyl [o-(Phenylazo)Phenyl]Nickel. *J. Am. Chem. Soc.* **1963**, *85*, 1544-1545.
52. Shiota, H.; Ano, Y.; Aihara, Y.; Fukumoto, Y.; Chatani, N. Nickel-Catalyzed Chelation-Assisted Transformations Involving Ortho C–H Bond Activation: Regioselective Oxidative Cycloaddition of Aromatic Amides to Alkynes. *J. Am. Chem. Soc.* **2011**, *133*, 14952-14955.
53. Roughley, S. D.; Jordan, A. M. The Medicinal Chemist's Toolbox: An Analysis of Reactions Used in the Pursuit of Drug Candidates. *J. Med. Chem.* **2011**, *54*, 3451-3479.
54. Barbero, N.; Martin, R. Ligand-Free Ni-Catalyzed Reductive Cleavage of Inert Carbon–Sulfur Bonds. *Org. Lett.* **2012**, *14*, 796-799.
55. Wolfe, J. P.; Buchwald, S. L. Nickel-Catalyzed Amination of Aryl Chlorides. *J. Am. Chem. Soc.* **1997**, *119*, 6054-6058.
56. Hie, L.; Baker, E. L.; Anthony, S. M.; Desrosiers, J.-N.; Senanayake, C.; Garg, N. K. Nickel-Catalyzed Esterification of Aliphatic Amides. *Angew. Chem., Int. Ed.* **2016**, *55*, 15129-15132.
57. Terao, J.; Watanabe, H.; Ikumi, A.; Kuniyasu, H.; Kambe, N. Nickel-Catalyzed Cross-Coupling Reaction of Grignard Reagents with Alkyl Halides and Tosylates: Remarkable Effect of 1,3-Butadienes. *J. Am. Chem. Soc.* **2002**, *124*, 4222-4223.
58. Hofstra, J. L.; Poremba, K. E.; Shimozone, A. M.; Reisman, S. E. Nickel-Catalyzed Conversion of Enol Triflates into Alkenyl Halides. *Angew. Chem., Int. Ed.* **2019**, *58*, 14901-14905.
59. Nattmann, L.; Cornella, J. Ni(⁴⁻¹⁰stb)₃: A Robust 16-Electron Ni(0) Olefin Complex for Catalysis. *Organometallics* **2020**, *39*, 3295-3300.

NICKEL (Compounds)

28-0010


Bis(1,5-cyclooctadiene)nickel (0), 98+% (1295-35-8)
 $(C_8H_{12})_2Ni$; FW: 275.08; yellow xtl.; m.p. 60° dec. (under N_2)
air sensitive, (store cold)




2g
10g

Technical Notes:


1. Pre-catalyst for the coupling of benzoxazole with N-tosylhydrazones.
2. Catalyst precursor for heterogeneously nickel-catalysed hydrogenolysis of aryl ethers without arene hydrogenation.
3. Pre-catalyst for reductive cleavage of C-OMe bonds with silanes as reducing agents.
4. Pre-catalyst for the cross-coupling reactions of benzylic pivalates with arylboroxines.
5. Pre-catalyst for the cross-coupling of benzylic carbamates with arylboronic esters.
6. Pre-catalyst for the direct arylation of $C(sp^3)-H$ Bonds in aliphatic amides via bidentate-chelation.
7. Pre-catalyst for the cross-coupling reactions of potassium alkoxyalkyl- and benzyltrifluoroborates with an array of aryl bromides and co-catalyzed by iridium photoredox catalyst under visible light at ambient temperature.
8. Pre-catalyst for highly regioselective indoline synthesis, co-catalyzed by Ruthenium photoredox catalyst.
9. Catalyst for conversion of amides to esters via the activation of amide C-N bonds.
10. Catalyst for borylation of aryl fluorides via C-F cleavage.
11. Catalyst for Suzuki-Miyaura coupling of amides.
12. Catalyst for the cross-coupling reaction of the aryl methyl ether alkylation.




Tech. Note (1)
Ref. (1)




Tech. Note (2)
Ref. (2)




Tech. Note (3)
Ref. (3)



Tech. Note (4)
Ref. (4)








Tech. Note (5)
Ref. (5)



Tech. Note (6)
Ref. (6)

NICKEL (Compounds)


28-0010 Bis(1,5-cyclooctadiene)nickel (0), 98+% (1295-35-8)
(continued)



References:

1. *Angew. Chem. Int. Ed.*, **2012**, 51,775.
2. *J. Am. Chem. Soc.*, **2012**, 134, 20226.
3. *J. Am. Chem. Soc.*, **2013**, 135, 1997.
4. *J. Am. Chem. Soc.*, **2013**, 135, 3307.
5. *J. Am. Chem. Soc.*, **2013**, 135, 3303.
6. *J. Am. Chem. Soc.*, **2014**, 136, 898.
7. *Science*, **2014**, 345, 433.
8. *J. Am. Chem. Soc.*, **2015**, 137, 9531.
9. *Nature*, **2015**, 524, 79.
10. *J. Am. Chem. Soc.*, **2015**, 137, 12470.
11. *Nature Chemistry*, **2016**, 8, 75.
12. *Angew. Chem. Int. Ed.*, **2016**, 55, 6093.

28-0040 (1,5-Cyclooctadiene)(duroquinone) nickel(0),
NEW min. 98% (40759-64-6)
C₁₈H₂₄NiO₂; FW: 331.08; orange to red solid;
m.p. 227(dec)




250mg
1g

Technical Notes:

1. Catalyst for C-CN bond activation in late-stage carbon isotope exchange of aryl nitriles.¹
2. Air stable Nickel(0)-olefin precatalyst for a variety of synthetic methods including Suzuki-Miyaura cross coupling reactions and amination of aryl chlorides.²


NICKEL (Compounds)

28-0040 (1,5-Cyclooctadiene)(duroquinone) nickel(0), min. 98% (40759-64-6)
(continued)




(Y = CH or N)


Late-stage incorporation of ^{13}C and ^{14}C labels
First carbon isotope exchange method for nitriles
Over 30 examples, with 10 complex pharmaceuticals
Aryl, heteroaryl, and alkenyl nitriles



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)




Tech. Note (2)
Ref. (2)

References:


1. *J. Am. Chem. Soc.* **2021**, *143*, 4817.
2. *Angew. Chem. Int. Ed.* **2020**, *59*, 7409.

28-0070 **Tris(trans-1,2-bis(4-tert-butylphenyl)ethene)nickel(0)**, min. 97% (2468315-70-8)
NEW $\text{C}_{66}\text{H}_{64}\text{Ni}$; FW: 936.07; red-orange powdr.
(store cold)
Note: Prolonged storage under inert atmosphere recommended.
Patents: EP19189236.3, DE 102019214138.2. Product sold under, use subject to, terms and conditions of Limited Use License found at www.strem.com/sgk1



Technical Note:


1. Temperature and air-stable Ni(0) source. Applications include a fast kinetic profile in Buchwald-Hartwig amination reactions. Other uses include various catalytic transformations with and without external ligands, and many other industrially relevant transformations.




Tech. Note (1)
Ref. (1)

NICKEL (Compounds)


28-0070 Tris(trans-1,2-bis(4-tert-butylphenyl)ethene)nickel(0), min. 97% (2468315-70-8)
(continued)




Tech. Note (1)
Ref. (1)




Tech. Note (1)
Ref. (1)



Tech. Note (1)
Ref. (1)



Tech. Note (1)
Ref. (1)



Tech. Note (1)
Ref. (1)

References:

1. *Organometallics*, **2020**, *39*, 3295-3300.

28-0060

NEW


Tris(trans-1,2-bis(4-(trifluoromethyl)phenyl)ethene)nickel(0), min. 97% (2413906-36-0)

C₄₈H₃₀F₁₈Ni; FW: 1007.41; red solid

(store cold)

Note: Prolonged storage under inert atmosphere recommended.


Patents: EP19189236.3, DE 102019214138.2. Product sold under, use subject to, terms and conditions of Limited Use License found at www.strem.com/sgk1



1g
5g

Technical Notes:


1. Air-stable Ni(0) source. Catalyst can be weighed out in air. Decomposition of the catalyst occurs at 25 °C eventually under air after ca. 1 week.
2. Applications include various catalytic transformations.
3. Electron-deficient stilbene prevents ligand hydrometallation.




Tech. Note (1)
Ref. (1)

NICKEL (Compounds)


28-0060 Tris(trans-1,2-bis(4-(trifluoromethyl)phenyl)ethene)nickel(0), min. 97% (2413906-36-0)
(continued)




Tech. Note (1)
Ref. (1)




Tech. Note (1)
Ref. (1)




Tech. Note (1,2)
Ref. (1)




Tech. Note (1,2)
Ref. (1)



Tech. Note (1,2)
Ref. (1)



Tech. Note (1,2)
Ref. (1)




Tech. Note (1,2)
Ref. (1)

NICKEL (Compounds)

28-0060 Tris(trans-1,2-bis(4-(trifluoromethyl)phenyl)ethene)nickel(0), min. 97% (2413906-36-0)
(continued)



Tech. Note (1,2)
Ref. (1)



Tech. Note (1,2,3)
Ref. (1)

References:

1. *Nat. Catal.* **2020**, 3, 6-13.

BERYLLIUM (Compounds)

04-4000

NEW

HAZ

Beryllium(II) acetylacetonate min. 97%


(10210-64-7)

 $\text{Be}(\text{CH}_3\text{COCHCOCH}_3)_2$; FW: 207.22; white

powdr.; m.p. 108; d. 1.116

moisture sensitive

Note: sublimates at 90 °C at 0.2 Torr



2g

10g


CERIUM (Compounds)

58-5500

NEW

Tetrakis[1-(methoxy)-2-methyl-2-propanolato]**cerium, [Ce(mmp)4], 98% (1262520-82-0)** $\text{C}_{20}\text{H}_{44}\text{CeO}_8$; FW: 552.68; light yellow to orange

powdr.

air sensitive, moisture sensitive

500mg

2g

Technical Note:

1. CVD and ALD alkoxide precursor for preparation of cerium thin films:

Film	Reactants/Conditions	Substrate/Temperature	Ref
CeO ₂	[Ce(mmp)4] in toluene; 1 mbar; Oxidants: O ₂ (MOCVD) or H ₂ O (ALD)	Si(100); 250-600°C (MOCVD); 150-350°C (ALD)	[1-3]
CeAlO ₃	[Ce(mmp)4] in toluene; 1 mbar; Oxidant: O ₂ (MOCVD); Et ₃ Al(OEt)	Si(100), Si(100)/TiN; 400-450°C	[4]
Pt@CeO ₂	[Ce(mmp)4] in cyclohexane; 10 hPa Oxidant: O ₂ (MOCVD); [MeCpPtMe ₃]	Si(100); 400°C	[5]
CeO ₂	[Ce(mmp)4] in cyclohexane; 10 mbar Oxidant: O ₂ (MOCVD)	Si(100), Carbon foil; 400 °C	[6]
CeO ₂	[Ce(mmp)4] in toluene; 1 mbar Oxidant: O ₂ (ALD)	Si(100), Si(100)/TiN; 250°C	[7]

Thermal behavior

- Under atmospheric pressure [Ce(mmp)4] shows a continual, gradual loss of mass over the temperature range 25-180°C, followed by rapid loss of mass at 180-250°C. Mass loss is complete at 300°C with remaining residue of 16.5% (probably CeO₂).
- [Ce(mmp)4] sublimates intact at 120°C without deposition of residues at low pressure (0.8 Torr).
- Thermal decomposition at ~275 °C [1].
- TGA data is available in [1] and [2].

References:

1. *Chem. Vap. Deposition* **2009**, 15, 259.
2. *Inorg. Chem.* **2011**, 50, 11644.
3. *Nanoscale Res. Lett.* **2013**, 8, 456.
4. *Thin Solid Films* **2013**, 536, 68.
5. *Thin Solid Films* **2015**, 589, 246.
6. *Surf. Coat. Tech.* **2015**, 280, 148.
7. *Beilstein J. Nanotechnol.* **2018**, 9, 890.

CHROMIUM (Compounds)

24-0600

Methylbenzoate chromium tricarbonyl, min.

NEW


HAZ

98% (12125-87-0)

C₁₁H₈CrO₅; FW: 272.18; orange solid;

m.p. 95-96°C

air sensitive, moisture sensitive





500mg

2g

Technical Note:

1. Chromium carbonyl based catalyst for the hydrogenation of conjugated dienes^{2,3} and alkynes.¹

Tech. Note (1)
Ref. (3)Tech. Note (1)
Ref. (1)

References:

1. *Synthesis* **1993**, 7, 643.
2. *Russ. Chem. Bull.* **2002**, 51, 1341.
3. *Russ. Chem. Bull.* **2018**, 67, 923.

ELECTROPOLISHED STAINLESS STEEL BUBBLERS (Vertical)

95-1050

Stainless steel bubbler, 600ml, vertical,
electropolished with fill-port, pneumatic
actuator valves, DOT 4B

NEW



1cyl

GADOLINIUM (Compounds)

64-4250

Gadolinium(III) bromide, anhydrous (99.99%-Gd) (REO) PURATREM
(13818-75-2)

NEW

GdBr₃; FW: 396.96; white xtl.; m.p. 770

air sensitive, moisture sensitive, hygroscopic

5g

25g

64-3575

Tris(N,N'-di-i-propylformamidinato)gadolinium(III), (99.999+%-Gd) PURATREM Gd-FMD
C₂₁H₄₅GdN₆; FW: 538.87; white to off-white powdr.


NEW

air sensitive, moisture sensitive

Note: Product sold under, use subject to, terms
and conditions of label license at
www.strem.com/harvard2.

1g

5g



KITS (Compounds)

96-0280

Chromalite® SpectraChrom Kit

NEW

Sold in collaboration with Purolite for research purposes only.

1 kit

Contains 25ml of the following:

Chromalite MIDA/M/Ni
Chromalite MIDA/M/Co
Chromalite MIDA/M/Zn
Chromalite MIDA/M/Cu
Chromalite MIDA/M/Fe
Chromalite MIDA/M

METALS SCAVENGING AGENTS (Compounds)

08-4225

Glucamine Alkyl Silica, borate scavenger

10g


50g

NEW

yellow-white powdr.

*moisture sensitive*Note: Particle size 37-200 μm

Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.



16-0800

Imidazole Sulfide Mercapto Alkyl Silica, metals and organic impurity scavenger

10g


50g

NEW

yellow-white powdr.

*moisture sensitive*Note: Particle size 37-200 μm

Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.



16-0805

Octadecyl Sulfide Propyl Silica, non-polar impurity scavenger

10g


50g

NEW

yellow-white powdr.

*moisture sensitive*Note: Particle size 37-200 μm

Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.



08-4230

Poly Carboxylic Acid Alkyl Silica, metals and organic impurity scavenger

10g


50g

NEW

yellow-white powdr.

*moisture sensitive*Note: Particle size 37-200 μm

Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.



16-0815

Poly Mercaptoalkyl Silica, metals and organic impurity scavenger

10g


50g

NEW

yellow-white powdr.

*moisture sensitive*Note: Particle size 20-200 μm

Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.



07-9255

Poly Quaternary Ammonium Alkyl Silica, metals and organic impurity scavenger

10g


50g

NEW

yellow-white powdr.


*moisture sensitive*Note: Particle size 37-200 μm

Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.




METALS SCAVENGING AGENTS (Compounds)

16-0810 Polyamine Sulfide Alkyl Silica, metals and organic impurity scavenger
NEW yellow-white powdr.
moisture sensitive
 Note: Particle size 37-200 μm
 Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.




10g
50g

16-0820 Tetraamine Sulfide Alkyl Silica, metals and organic impurity scavenger
NEW yellow-white powdr.
moisture sensitive
 Note: Particle size 37-200 μm
 Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.



10g
50g

16-0825 Thiourea Sulfide Mercapto Silica, metals and organic impurity scavenger
NEW yellow-white powdr.
moisture sensitive
 Note: Particle size 37-200 μm
 Sold in collaboration with Si-Novation for research purposes only. Patent: PCT/CN2016/092663.



10g
50g

MOFS AND LIGANDS FOR MOF SYNTHESIS (Compounds)

29-3150 Copper benzene-1,3,5-tricarboxylate MOF (HKUST-1(Cu))/PVDF membrane
NEW (35/65 wt.%) (51937-85-0)
 blue membrane; SA: BET 500 m^2/g
moisture sensitive


1x1cm

40-1125 Zirconium 1,4-dicarboxybenzene MOF (UiO-66)/PVDF membrane (60/40 wt.%) (1072413-89-8)
NEW white membrane; SA: BET 650-700 m^2/g
moisture sensitive

1x1cm

NICKEL (Compounds)

28-0040 (1,5-Cyclooctadiene)(duroquinone) nickel(0),
NEW min. 98% (40759-64-6)
 $\text{C}_{18}\text{H}_{24}\text{NiO}_2$; FW: 331.08; orange to red solid;
 m.p. 227(dec)




250mg
1g

Technical Notes:

1. Catalyst for C-CN bond activation in late-stage carbon isotope exchange of aryl nitriles.¹
2. Air stable Nickel(0)-olefin precatalyst for a variety of synthetic methods including Suzuki-Miyaura cross coupling reactions and amination of aryl chlorides.²


NICKEL (Compounds)

28-0040 (1,5-Cyclooctadiene)(duroquinone) nickel(0), min. 98% (40759-64-6)
(continued)




(Y = CH or N)


Late-stage incorporation of ^{13}C and ^{14}C labels
First carbon isotope exchange method for nitriles
Over 30 examples, with 10 complex pharmaceuticals
Aryl, heteroaryl, and alkenyl nitriles



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)



Tech. Note (2)
Ref. (2)

References:

1. *J. Am. Chem. Soc.* **2021**, *143*, 12, 4817-4823.
2. *Angew. Chem. Int. Ed.* **2020**, *59*, 7409.

28-0070

NEW

Tris(trans-1,2-bis(4-tert-butylphenyl)ethene)

nickel(0), min. 97% (2468315-70-8)


$\text{C}_{66}\text{H}_{84}\text{Ni}$; FW: 936.07; red-orange powd.

(store cold)

Note: Prolonged storage under inert atmosphere recommended.

Patents: EP19189236.3, DE 102019214138.2.

Product sold under, use subject to, terms and conditions of Limited Use License found at www.strem.com/sgk1




1g

5g

Technical Note:

1. Temperature and air-stable Ni(0) source. Applications include a fast kinetic profile in Buchwald-Hartwig amination reactions. Other uses include various catalytic transformations with and without external ligands, and many other industrially relevant transformations.




Tech. Note (1)
Ref. (1)


NICKEL (Compounds)

28-0070
(continued)


Tris(trans-1,2-bis(4-tert-butylphenyl)ethene)nickel(0), min. 97% (2468315-70-8)




Tech. Note (1)
Ref. (1)




Tech. Note (1)
Ref. (1)



Tech. Note (1)
Ref. (1)



Tech. Note (1)
Ref. (1)



Tech. Note (1)
Ref. (1)

References:

1. *Organometallics*, **2020**, 39, 3295-3300.


NITROGEN (Compounds)

07-1070

NEW

(1Z,3Z)-1,3-Bis[[(4S)-4,5-dihydro-4-phenyl-2-oxazolyl]methylene]-2,3-dihydro-5,6-dimethyl-1H-isoindole, 95% (1358991-79-3)
C₃₀H₂₇N₃O₂; FW: 461.6; light yellow to yellow powd.
moisture sensitive, (store cold)
Note: Sold in collaboration with Daicel for research purposes only

25mg
100mg




NITROGEN (Compounds)


07-1070 (1Z,3Z)-1,3-Bis[[(4S)-4,5-dihydro-4-phenyl-2-oxazolyl]methylene]-2,3-dihydro-5,6-dimethyl-1H-isoindole, 95% (1358991-79-3)

Technical Notes:


1. Ligand for the Ni-catalyzed enantioselective fluorination of oxindoles.
2. Ligand used in Nozaki-Hiyama-Kishi reaction.
3. Ligand for the highly enantioselective Cu-catalyzed alkylation of β -ketoesters with benzyl and allylic alcohols.
4. Ligand for the highly enantioselective Cu-catalyzed electrophilic trifluoromethylation and trifluoromethylthiolations of β -ketoesters.




Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (1)



Tech. Note (3)
Ref. (2)



Tech. Note (4)
Ref. (3-4)

References:

1. *Chem. Eur. J.* **2011**, 17, 14922.
2. *J. Am. Chem. Soc.* **2012**, 134, 2946.
3. *J. Am. Chem. Soc.* **2012**, 134, 10769.
4. *Chem. Eur. J.* **2014**, 20, 93.


07-1065**NEW**

(1Z,3Z)-1,3-Bis[[(4R)-4,5-dihydro-4-phenyl-2-oxazolyl]methylene]-2,3-dihydro-1H-isoindole, 95% (1429056-54-1)

$C_{28}H_{23}N_3O_2$; FW: 433.5; light yellow to yellow powd.

moisture sensitive, (store cold)

Note: Sold in collaboration with Daicel for research purposes only.




25mg
100mg

Technical Notes:



1. Ligand for the enantioselective Fe-catalyzed azidation of β -ketoesters and oxindoles.
2. Ligand for the enantioselective Fe-catalyzed hydrosilylation of ketones.
3. Ligand for the enantioselective Zn-catalyzed allylation of oxindole derivatives with allyl bromide.
4. Ligand for the Zn-catalyzed asymmetric electrophilic cyanation of 3-substituted oxindoles.
5. Ligand for the α -selective Co-catalyzed hydroborations of alkynes.

NITROGEN (Compounds)


07-1065 (1Z,3Z)-1,3-Bis[[(4R)-4,5-dihydro-4-phenyl-2-oxazolyl]methylene]-2,3-dihydro-1H-isoidole, 95% (1429056-54-1)




Tech. Note (1)
Ref. (1)




Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)



Tech. Note (4)
Ref. (4)



Tech. Note (5)
Ref. (5)

References:

1. *J. Am. Chem. Soc.* **2013**, 135, 5356.
2. *J. Am. Chem. Soc.* **2015**, 137, 2456.
3. *Angew. Chem. Int. Ed.* **2016**, 55, 7852.
4. *Org. Lett.* **2017**, 19, 4018.
5. *Angew. Chem. Int. Ed.* **2020**, 59, 23010.

07-1068 (1Z,3Z)-1,3-Bis[[(4S)-4,5-dihydro-4-(tert-butyl)-2-oxazolyl]methylene]-2,3-dihydro-1H-isoidole, 98% (1361563-41-8)


NEW

$C_{24}H_{31}N_3O_2$; FW: 393.5; light yellow to yellow powdr.

moisture sensitive, (store cold)

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg




Technical Notes:

1. Ligand for the Ni-catalyzed enantioselective fluorination of oxindoles.
2. Ligand used in Nozaki-Hiyama-Kishi reaction.


NITROGEN (Compounds)

07-1068
(continued)

(1Z,3Z)-1,3-Bis[[[(4S)-4,5-dihydro-4-(tert-butyl)-2-oxazolyl]methylene]-2,3-dihydro-1H-isoin-
dole, 98% (1361563-41-8)



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (1)


References:

1. *Chem. Eur. J.* **2011**, *17*, 14922.

07-8230

NEW


3,6-Bis(dimethylamino)-9-(2,6-dimethyl-phenyl)-10-methylacridin-10-ium bromide, min. 95% (2380163-04-0)
C₂₈H₃₀BrN₃; FW: 464.45; orange to red to dark red solid
air sensitive
Note: Sold in collaboration with Solvias for research purposes only.




50mg
250mg

Technical Notes:

1. Organophotoredox catalyst used in Ni-assisted C-C cross-coupling without C-O bond formation.
2. Acridinium catalysts used for photoredox deuteration of clomipramine.



Tech. Note (1)
Ref. (1, 3)



Tech. Note (2)
Ref. (2)


References:

1. *Angew. Chem. Int. Ed.* **2018**, *57*, 2436.
2. *Synthesis* **2019**, *51*, 4359.
3. *ACS Catal.* **2020**, *10*, 210.

07-8225

NEW

3,6-Bis(dimethylamino)-9-mesityl-10-methylacridin-10-ium bromide, min. 95% (2180894-90-8)
C₂₇H₃₂BrN₃; FW: 478.48; orange to red to dark red solid
air sensitive
Note: Sold in collaboration with Solvias for research purposes only.



50mg
250mg

Technical Note:

1. See 07-8230 (page 29)


NITROGEN (Compounds)

07-1058

NEW


2,2-Bis[(4R)-4-tert-butyl-2-oxazolin-2-yl]propane, 98%, (99% ee) (131833-97-1)
 $C_{17}H_{30}N_2O_2$; FW: 294.4; off white to light yellow powdr.
moisture sensitive, (store cold)
 Note: Sold in collaboration with Daicel for research purposes only.

100mg




Technical Notes:


1. Ligand for the enantioselective H-atom transfer Reactions to generate β 2-amino acids.
2. Ligand for the Cu-catalyzed asymmetric alkylation of β -keto esters with xanthydrols.
3. Ligand used in Zn-catalyzed asymmetric hetero-Diels-Alder reaction of trans-1-methoxy-3-trimethylsilyloxy-buta-1,3-diene.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)


References:

1. *Angew. Chem. Int. Ed.* **2004**, 43, 1235.
2. *Adv. Synth. Catal.* **2013**, 355, 2815.
3. *Eur. J. Org. Chem.* **2020**, 5388.

07-1076

NEW


(4R,4'R,5S,5'S)-2,2'-Cyclobutylidenebis[4,5-dihydro-4,5-diphenyloxazole], 98%
 $C_{34}H_{30}N_2O_2$; FW: 498.6; off white to light yellow powdr.
moisture sensitive, (store cold)
 Note: Sold in collaboration with Daicel for research purposes only

25mg
100mg

07-1079


NEW

(4S,4'S)-2,2'-Cyclobutylidenebis[4,5-dihydro-4-tert-butyloxazole], min. 95% (298693-02-4)
 $C_{18}H_{30}N_2O_2$; FW: 306.4; off white to light yellow powdr.
moisture sensitive, (store cold)
 Note: Sold in collaboration with Daicel for research purposes only

25mg
100mg

Technical Note:

1. Ligand used in the enantioselective organolithium additions to imines.



Tech. Note (1)
Ref. (1)

References:

1. *Adv. Synth. Catal.* **2008**, 350, 1023.

NITROGEN (Compounds)


07-1062

NEW

(4S,4'S)-2,2'-Cyclopropylidenebis[4-tert-butyl-4,5-dihydrooxazole], 95% (195379-09-0)C₁₇H₂₈N₂O₂; FW: 292.4; off white to light yellow powdr.






moisture sensitive, (store cold)

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

Technical Notes:

1. Ligand for the highly enantioselective Cu-catalyzed Friedel-Crafts alkylations of pyrroles and indoles with α -hydroxy enones.
2. Ligand used in Lewis acid catalyzed Aldol reactions to generate optically active pyridine derivatives from prochiral pyridine-N-oxides.
3. Ligand used in Cu-catalyzed tandem Nazarov cyclization/Wagner-Meerwein rearrangement to synthesize Enokipodin B.
4. Ligand used in enantioselective Cu-catalyzed oxy-alkynylation of diazo compounds.
5. Ligand used in Cu-catalyzed oxyvinylation of diazo compounds.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)Tech. Note (4)
Ref. (4)Tech. Note (5)
Ref. (5)

References:

1. *J. Am. Chem. Soc.* **2005**, 127, 4154.
2. *Chem. Eur. J.* **2006**, 12, 3472.
3. *Chem. Eur. J.* **2013**, 19, 4835.
4. *J. Am. Chem. Soc.* **2017**, 139, 8420.
5. *Org. Lett.* **2020**, 22, 3884.

NITROGEN (Compounds)


07-1073

NEW

(4*S*,4'*S*)-2,2'-(4,6-Dibenzofurandiyl)bis[4,5-dihydro-4-phenyloxazole], 98% (246040-77-7)
 $C_{30}H_{22}N_2O_3$; FW: 458.5; off white to light yellow powdr.

moisture sensitive, (store cold)


Note: Sold in collaboration with Daicel for research purposes only




25mg
100mg

Technical Notes:


1. Ligand for the Ni-catalyzed enantioselective addition of nitrones to activated cyclopropanes.
2. Ligand for the Zn-catalyzed Friedel-Crafts alkylation with α -substituted acrylates.
3. Ligand for the Lewis acid mediated enantioselective intramolecular redox reactions to generate ring-fused tetrahydroquinolines.
4. Ligand for the Lewis acid-catalyzed highly enantioselective α -cyanation with 4-acetylphenyl cyanate.
5. Ligand for the Ni-catalyzed asymmetric 1,3-dipolar cycloaddition of β -fluoroalkylated α,β -Unsaturated 2-pyridylsulfones with nitrones to generate chiral fluoroalkylated isoxazolidines and γ -amino alcohols.




Tech. Note (1)
Ref. (1)




Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)



Tech. Note (4)
Ref. (4)



Tech. Note (5)
Ref. (5)

References:

1. *J. Am. Chem. Soc.* **2005**, 127, 5764.
2. *Angew. Chem. Int. Ed.* **2008**, 47, 9913.
3. *J. Am. Chem. Soc.* **2009**, 131, 13226.
4. *Chem. Eur. J.* **2017**, 23, 1775.
5. *Angew. Chem. Int. Ed.* **2017**, 56, 1510.

NITROGEN (Compounds)

07-1085


NEW

(4R,4'R,5S,5'S)-2,2'-(1,3-Dihydro-2H-inden-2-ylidene)bis[4,5-dihydro-4,5-diphenyloxazole], 98% (1656253-81-4)

C₃₉H₃₂N₂O₂; FW: 560.7; off white to light yellow powdr.

moisture sensitive, (store cold)


Note: Sold in collaboration with Daicel for research purposes only



25mg
100mg

Technical Note:

1. Ligand used in Cu(II)-catalyzed enantioselective intramolecular cyclization of N-alkenylureas.



Tech. Note (1)
Ref. (1)

References:

1. *Org. Lett.* **2015**, 17, 1018.

07-1434


NEW

(4R,4'R)-2,2'-(1,3-Dihydro-2H-inden-2-ylidene) bis[4,5-dihydro-4-isopropoxyloxazole], 98%, (99% ee)

C₂₁H₂₈N₂O₂; FW: 340.5; off white to light yellow powdr.

moisture sensitive, (store cold)

Note: Sold in collaboration with Daicel for research purposes only



100mg

07-1435


NEW

(4S,4'S)-2,2'-(1,3-Dihydro-2H-inden-2-ylidene) bis[4,5-dihydro-4-isopropoxyloxazole], 98%, (99% ee)

C₂₁H₂₈N₂O₂; FW: 340.5; off white to light yellow powdr.

moisture sensitive, (store cold)

Note: Sold in collaboration with Daicel for research purposes only



100mg

07-1082


NEW

(4R,4'R)-2,2'-(1,3-Dihydro-2H-inden-2-ylidene) bis[4,5-dihydro-4-tert-butyloxazole], 98%, (99% ee)

C₂₃H₃₂N₂O₂; FW: 368.5; off white to light yellow powdr.

moisture sensitive, (store cold)

Note: Sold in collaboration with Daicel for research purposes only



100mg

07-1049


NEW

2-[(4R)-4,5-Dihydro-4-isopropyl-2-oxazolyl]pyridine, 98% (99% ee) (132187-16-7)

C₁₁H₁₄N₂O; FW: 190.2; light yellow to yellow powdr.

moisture sensitive, (store cold)


Note: Sold in collaboration with Daicel for research purposes only.



100mg

Technical Note:

1. Ligand for the Pd-catalyzed enantioselective diene cyclization/hydrosilylation.



Tech. Note (1)
Ref. (1)

References:

1. *J. Org. Chem.* **2000**, 65, 3836.

NITROGEN (Compounds)

07-1048

NEW

2-[(4S)-4,5-Dihydro-4-isopropyl-2-oxazolyl]


pyridine, 98%, (99% ee) (108915-04-4)

C₁₁H₁₄N₂O; FW: 190.2; off white to light yellow powdr.

moisture sensitive, (store cold)





Note: Sold in collaboration with Daicel for research purposes only.

100mg




Technical Notes:

1. Ligand used in the Pd-catalyzed asymmetric aerobic Aza-Wacker-Type cyclization.
2. Ligand used in Pd-catalyzed conjugate addition of arylboronic acids to β -substituted cyclic enones.
3. Ligand for enantioselective Ir-catalyzed silylation of aromatic C–H bonds.
4. Ligand for the Co-catalyzed enantioselective Markovnikov 1,2-hydrosilylation of conjugated dienes.
5. Ligand used for the Ni-catalyzed intramolecular desymmetrization addition of aryl halides to 1,3-diketones.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)Tech. Note (4)
Ref. (4)

R = Ar, Alk; R' = H, Ar, Alk

Tech. Note (5)
Ref. (5)

References:

1. *Angew. Chem. Int. Ed.* **2012**, 51, 9141.
2. *J. Am. Chem. Soc.* **2013**, 135, 14996.
3. *Angew. Chem. Int. Ed.* **2017**, 56, 1092.
4. *ACS Catal.* **2019**, 9, 1612.
5. *Chem. Commun.*, **2020**, 56, 8194.

NITROGEN (Compounds)

07-1340

NEW



Dodecabenzylbambus[6]uril, Bn-BU[6]

(1308315-95-8)

$C_{114}H_{108}N_{24}O_{12}$; FW: 2006.23; white powdr.

Note: Patents: CZ302710B6, EP2501699B1, US8779120, WO2011057590A1.

100mg



Technical Notes:

1. Exhibits exceptionally strong binding of anions with different sizes in chloroform.
2. Macrocyclic compound used in the real-time analysis of anion mixtures by NMR methods.
3. Used as a transmembrane Cl-/HCO3- antiporter.
4. Used as a carrier in tailor-made liquid membranes for highly selective electromembrane extractions of inorganic anions.
5. Used for anion complexation in the study of the influence of counterions on the [N-I-N]+ halogen bond in solution, in the solid state and in silico.

References:

1. *ChemPlusChem.*, **2015**, *80*, 1601-1606.
2. *Chem. Commun.*, **2015**, *51*, 4666.
3. *Chem.*, **2019**, *5*, 429-444.
4. *Anal. Chim. Acta*, **2017**, *950*, 49-56.
5. *Chem. Sci.*, **2015**, *6*, 3746-3756.

07-1342

NEW



Dodecamethylbambus[6]uril hydrate, BU[6]

(1227292-62-7)

$C_{42}H_{60}N_{24}O_{12} \cdot xH_2O$; FW: 1093.08 ; white powdr.

Note: Patents: CZ302710B6, EP2501699B1, US8779120, WO2011057590A1.

100mg



Technical Note:

1. A new macrocyclic compound which is able to bind anions with high affinity and selectivity.

References:

1. *Chem. Eur. J.*, **2011**, *17*, 5605-5612.
2. *Org. Lett.*, **2011**, *13*(15), 4000.

07-1041


NEW

**2-[(4S)-4-Isobutyl-4,5-dihydro-2-oxaz-
olyl]-5-(trifluoromethyl)pyridine, 98%, (99%
ee)**

$C_{13}H_{15}F_3N_2O$; FW: 272.3; off white to light yellow
powdr.

moisture sensitive, (store cold)

Note: Sold in collaboration with Daicel for
research purposes only



100mg

07-7059


NEW

**2-[(4S)-4-Isopropyl-4,5-dihydro-2-oxaz-
olyl]-5-(trifluoromethyl)pyridine, 98%, (99%
ee)**

$C_{12}H_{13}F_3N_2O$; FW: 258.2; off white to light yellow
powdr.

moisture sensitive, (store cold)

Note: Sold in collaboration with Daicel for
research purposes only.




100mg

NITROGEN (Compounds)

07-1088


NEW

(4*R*,4'*R*)-2,2'-(1-Phenylpropane-2,2-diyl)bis(4,5,5-triphenyl-4,5-dihydrooxazole), 98% (2409652-69-1)
 $C_{51}H_{42}N_2O_2$; FW: 714.9; off white to light yellow powdr.
moisture sensitive, (store cold)
 Note: Sold in collaboration with Daicel for research purposes only

25mg
100mg

Technical Note:

- Ligand for the Ni-catalyzed enantioconvergent reductive hydroalkylation of olefins with α -heteroatom phosphorus alkyl electrophiles.

Tech. Note (1)
Ref. (1)


References:

- J. Am. Chem. Soc.* **2020**, *142*, 214.

07-1055


NEW

1-[(4*R*)-4-*tert*-Butyl-4,5-dihydro-2-oxazolyl]isoquinoline, 98%, (99% ee) (1402851-53-9)
 $C_{16}H_{18}N_2O$; FW: 254.3; off white to light yellow powdr.
moisture sensitive, (store cold)
 Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

Technical Note:

- Ligand for the Ni-catalyzed asymmetric cross-coupling of secondary electrophiles with secondary nucleophiles, specifically, stereoconvergent Negishi reactions of racemic benzylic bromides with achiral cycloalkylzinc reagents.

Tech. Note (1)
Ref. (1)


References:

- J. Am. Chem. Soc.* **2012**, *134*, 17003.

07-1045

NEW

2-[(4*R*)-4-*tert*-Butyl-4,5-dihydro-2-oxazolyl]pyridine, 95%, (99% ee) (242482-28-6)
 $C_{12}H_{16}N_2O$; FW: 204.3; off white to light yellow powdr.
moisture sensitive, (store cold)
 Note: Sold in collaboration with Daicel for research purposes only

25mg
100mg

Technical Notes:




- Ligand for asymmetric Pd-catalyzed directed intermolecular fluoroarylation of styrenes.
- Ligand used in Ir visible light photoredox/Ni-catalyzed cross-coupling of acyl chlorides with potassium alkoxymethyltrifluoroborates.
- Ligand for Ni-catalyzed asymmetric intramolecular reductive Heck reaction of unactivated alkenes.
- Ligand used in Au(III)-catalyzed intramolecular alkoxycyclization of a 1,6-enyne.
- Ligand for the Ni-catalyzed asymmetric reductive arylbenzylation of unactivated alkene.

NITROGEN (Compounds)



07-1045

2-[(4R)-4-tert-Butyl-4,5-dihydro-2-oxazolyl]pyridine, 95%, (99% ee) (242482-28-6)

(continued)



R₃= H or D; X=CH₂, O, NTs



X=CH₂, O, NTs

References:

1. *J. Am. Chem. Soc.* **2014**, 136, 4101.
2. *Org. Lett.* **2016**, 18, 732.
3. *Org. Lett.* **2019**, 21, 6989.
4. *J. Am. Chem. Soc.* **2019**, 141, 18221.
5. *Org. Lett.* **2020**, 22, 2724.

NITROGEN (Compounds)

07-1044


NEW

2-[(4S)-4-tert-Butyl-4,5-dihydro-2-oxazoly]pyridine, 98%, (99% ee) (117408-98-7)C₁₂H₁₆N₂O; FW: 204.3; off white to light yellow powdr.

moisture sensitive, (store cold)





Note: Sold in collaboration with Daicel for research purposes only

100mg




Technical Notes:

1. Ligand used in Pd-catalyzed conjugate addition of arylboronic acids to β -substituted cyclic enones.
2. Ligand for the enantioselective Pd-catalyzed oxidative β,β -fluoroarylation of α,β -unsaturated carbonyl derivatives.
3. Ligand for the Pd-catalyzed asymmetric Aza-Wacker-Type cyclization of N-Ts hydrazine-tethered tetrasubstituted olefins.
4. Ligand for the Co-catalyzed enantioselective Markovnikov 1,2-hydrosilylation of conjugated dienes.
5. Ligand for the Pd-catalyzed enantioselective desymmetrization of polycyclic cyclohexenediones.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)EWG = COR, CO₂R, CONR₁R₂, CO₂H, CH₂CN, CHOTech. Note (3)
Ref. (3)Tech. Note (4)
Ref. (4)

R = Ar, Alk; R' = H, Ar, Alk

Tech. Note (5)
Ref. (5)

References:


1. *J. Am. Chem. Soc.* **2013**, 135, 14996.
2. *Angew. Chem. Int. Ed.* **2016**, 55, 9045
3. *J. Am. Chem. Soc.* **2018**, 140, 7587.
4. *ACS Catal.* **2019**, 9, 1612.
5. *Org. Lett.* **2019**, 21, 8689.

NITROGEN (Compounds)

07-1050

NEW


2-[(4S)-4-tert-Butyl-4,5-dihydro-2-oxazolyl]-4-(tri-fluoromethyl)pyridine, 98% (1257527-14-2)
 $C_{13}H_{15}F_3N_2O$; FW: 272.3; off white to light yellow powdr.
moisture sensitive, (store cold)
 Note: Sold in collaboration with Daicel for research purposes only.




100mg

Technical Notes:


1. Ligand used in enantioselective Pd-catalyzed alkene difunctionalization reaction.
2. Ligand used in enantioselective Heck arylations of acyclic alkenyl alcohols.
3. Ligand for enantioselective Pd-catalyzed conjugate addition of arylboronic acids to β -substituted cyclic enones.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)

References:


1. *J. Am. Chem. Soc.* **2010**, 132, 17471.
2. *Science* **2012**, 338, 1455.
3. *J. Am. Chem. Soc.* **2013**, 135, 14996.

PHOSPHORUS (Compounds)

15-8406

NEW


[S(R)]-N-((1R)-(Adamantan-1-yl)(5-(diphenyl-phosphanyl)-9,9-dimethyl-9H-xanthen-4-yl)methyl)-2-methyl-2-propanesulfonamide, 95% (2162939-91-3)
 $C_{42}H_{48}NO_2PS$; FW: 661.9; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



25mg

Technical Note:

1. Ligand for the Pd-catalyzed enantioselective arylation of general sulfonate anions.



Tech. Note (1)
Ref. (1)

R¹ = Alkyl, Aryl

X = I, Br

References:

1. *J. Am. Chem. Soc.* **2018**, 140, 3467.

PHOSPHORUS (Compounds)

15-8404

NEW

[S(R)]-N-[(1S)-(Adamantan-1-yl)(5-(diphenylphosphanyl)-9,9-dimethyl-9H-xanthen-4-yl)methyl]-2-methyl-2-propanesulfonamide, 95% (2183514-08-9)
 $C_{42}H_{48}NO_2PS$; FW: 661.9; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg

Technical Note:


1. See 15-8406 (page 39)

15-8370

NEW


[S(R)]-N-[(1S)-1-[2-(9-Anthracenyl)phenyl]-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95% (1936438-22-0)
 $C_{38}H_{36}NOPS$; FW: 585.7
air sensitive, moisture sensitive, (store cold)
 Note: Sold in collaboration with Daicel for research purposes only.

50mg




Technical Notes:


1. Organocatalyst used in enantioselective intermolecular cross Rauhut-Currier reactions of activated alkenes with acrolein.
2. Used in phosphine-catalyzed asymmetric umpolung addition of trifluoromethyl ketimines to Morita-Baylis-Hillman carbonates.
3. Used in the enantioselective phosphine-catalyzed allylic alkylations of mix-indene with MBH carbonates.
4. Used in phosphine-catalyzed asymmetric synthesis of α -quaternary amine via umpolung γ -addition of ketimines to allenates.
5. Ligand for the direct Pd-catalyzed cross-coupling reaction of easily accessible secondary phosphine oxides and aryl bromides, which provides rapid access to P-chiral phosphine oxides.




Tech. Note (1)
Ref. (1)




Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)



Tech. Note (4)
Ref. (4)

EWG = CF₃, ester

Tech. Note (5)
Ref. (5)

PHOSPHORUS (Compounds)


15-8370 [S(R)]-N-[(1S)-1-[2-(9-Anthracenyl)phenyl]-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95% (1936438-22-0)

References:

1. *Chem. Commun.*, **2016**, 52, 7612.
2. *Angew. Chem. Int. Ed.* **2016**, 55, 13316.
3. *Org. Lett.* **2017**, 19, 6080.
4. *Org. Lett.* **2017**, 19, 6550.
5. *J. Am. Chem. Soc.* **2019**, 141, 20556.

15-8380**NEW**


[S(R)]-N-[(1S)-1-[2-(9-Anthracenyl)phenyl]-2-(diphenylphosphino)ethyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{38}H_{38}NO_2PS$; FW: 599.8; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8478**NEW**


[S(R)]-N-[(S)-[(3-(Benzyloxy)-2-(dicyclohexylphosphino)phenyl)-(1,1'-biphenyl)methyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{43}H_{54}NO_2PS$; FW: 679.9; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8482**NEW**


[S(R)]-N-[(S)-[(3-(Benzyloxy)-2-(dicyclohexylphosphino)phenyl)-(2-naphthalenyl)methyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{41}H_{52}NO_2PS$; FW: 653.9; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8480**NEW**


[S(R)]-N-[(S)-[(3-(Benzyloxy)-2-(dicyclohexylphosphino)phenyl)-(2-naphthalenyl)methyl]-2-methyl-2-propanesulfonamide, 95%
 $C_{40}H_{50}NO_2PS$; FW: 639.9; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8476**NEW**

[S(R)]-N-[(R)-[(3-(Benzyloxy)-2-(dicyclohexylphosphino)phenyl)phenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{37}H_{50}NO_2PS$; FW: 603.8; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



25mg

PHOSPHORUS (Compounds)

15-8474


NEW

[S(R)]-N-[(S)-[(3-(Benzyloxy)-2-(dicyclohexylphosphino)phenyl)phenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%

C₃₇H₅₀NO₂PS; FW: 603.8; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8472


NEW

[S(R)]-N-[(R)-[(3-(Benzyloxy)-2-(dicyclohexylphosphino)phenyl)phenylmethyl]-2-methyl-2-propanesulfonamide, 95%

C₃₆H₄₈NO₂PS; FW: 589.8; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8470


NEW

[S(R)]-N-[(S)-[(3-(Benzyloxy)-2-(dicyclohexylphosphino)phenyl)phenylmethyl]-2-methyl-2-propanesulfonamide, 95%

C₃₆H₄₈NO₂PS; FW: 589.8; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8506


NEW

N-[(1S)-1-[1',1'-Biphenyl]-2-yl-2-(diphenylphosphino)ethyl]-3,5-bis(trifluoromethyl)-benzamide, 95% (2089424-10-0)

C₃₆H₂₆F₆NOP; FW: 621.6; white to off-white powdr.

air sensitive, moisture sensitive


Note: Sold in collaboration with Daicel for research purposes only.




50mg

Technical Notes:


1. Organocatalyst used in aza-[4+2] cycloaddition of N-sulfonyl-1-aza-1,3-dienes with vinyl ketones.
2. Catalyst for the enantioselective [3+2] cycloadditions of γ -substituted allenates with β -perfluoroalkyl enones.
3. Catalyst for asymmetric intermolecular cross Rauhut-Currier reaction of β -perfluoroalkyl-substituted enones and vinyl ketones.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)

References:

1. *Org. Lett.* **2017**, 19, 1710.
2. *Chem. Sci.*, **2017**, 8, 4660.
3. *Adv. Synth. Catal.* **2017**, 359, 3347.

PHOSPHORUS (Compounds)

15-8368

NEW


[S(R)]-N-[(1S)-1-[1,1'-Biphenyl]-2-yl-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95% (1936438-14-0)

C₃₀H₃₂NOPS; FW: 485.6; white to off-white powder.

air sensitive, moisture sensitive


Note: Sold in collaboration with Daicel for research purposes only.

100mg




Technical Notes:


1. Organocatalyst used in enantioselective intermolecular cross Rauht-Currier reactions of activated alkenes with acrolein.
2. Organocatalyst used in aza-[4+2] cycloaddition of N-sulfonyl-1-aza-1,3-dienes with vinyl ketones.
3. Catalyst for the enantioselective [3+2] cycloadditions of γ -substituted allenates with β -perfluoroalkyl enones.
4. Ligand for the direct Pd-catalyzed cross-coupling reaction of easily accessible secondary phosphine oxides and aryl bromides, which provides rapid access to P-chiral phosphine oxides.




Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)



Tech. Note (4)
Ref. (4)

References:


1. *Chem. Commun.*, **2016**, 52, 7612.
2. *Org. Lett.* **2017**, 19, 1710.
3. *Chem. Sci.*, **2017**, 8, 4660.
4. *J. Am. Chem. Soc.* **2019**, 141, 20556.

PHOSPHORUS (Compounds)

15-8520




NEW

(11bR)-2,6-Bis[4-(9-anthracenyl)-2,6-bis(isopropyl)phenyl]-4-hydroxy-4-oxide-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin, 95%, (99% ee)
(1236191-19-7)
 $C_{72}H_{61}O_4P$; FW: 1021.2; off white to light yellow powdr.
Note: Sold in collaboration with Daicel for research purposes only.

10mg
25mg

Technical Notes:

1. Organocatalyst for enantioselective Friedel-Crafts alkylation reaction of indoles with α,β -unsaturated acyl phosphonates.
2. Used in the Lewis base/chiral Brønsted acid catalyzed enantioselective bromocycloetherification.
3. Chiral phosphoric acid catalyst for the kinetic resolution of indolines based on a self-redox reaction.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)

References:

1. *Chem. Commun.* **2010**, 46, 4112.
2. *Org. Lett.*, **2012**, 14, 256.
3. *Angew. Chem. Int. Ed.* **2016**, 55, 3148.

15-8521


NEW

(11bS)-2,6-Bis[4-(9-anthracenyl)-2,6-bis(isopropyl)phenyl]-4-hydroxy-4-oxide-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin, 95%, (99% ee)
(1051435-82-5)
 $C_{72}H_{61}O_4P$; FW: 1021.2; off white to light yellow powdr.
Note: Sold in collaboration with Daicel for research purposes only.

10mg
25mg


Technical Notes:

1. Used in direct Brønsted acid catalyzed asymmetric three-component Kabachnik-Fields reaction.
2. Used in direct catalytic asymmetric synthesis of cyclic aminals from aldehydes.
3. Brønsted acid catalyst for asymmetric acetalization of aliphatic aldehydes.


Tech. Note (1)
Ref. (1)

PHOSPHORUS (Compounds)

15-8521 (11bS)-2,6-Bis[4-(9-anthracenyl)-2,6-bis(isopropyl)phenyl]-4-hydroxy-4-oxide-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphocin, 95%, (99% ee) (1051435-82-5)



Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)


References:

1. *Angew. Chem. Int. Ed.* **2013**, 52, 4474.
2. *J. Am. Chem. Soc.* **2008**, 130, 15786.
3. *Angew. Chem. Int. Ed.* **2008**, 47, 5079.

15-0564

NEW

(11aS)-3,7-Bis[3,5-bis(trifluoromethyl)phenyl]phenyl]-10,11,12,13-tetrahydro-5-hydroxy-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee)
C₆₁H₃₁F₂₄O₄P; FW: 1314.8; off white to light yellow powdr.
Note: Sold in collaboration with Daicel for research purposes only.




25mg
100mg

15-8588

NEW

(11aR)-3,7-Bis[3,5-bis(tert-butyl)-4-methoxyphenyl]-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee)
C₅₁H₆₇O₆P; FW: 807; off white to light yellow powdr.
Note: Sold in collaboration with Daicel for research purposes only.



25mg
100mg

15-8589

NEW

(11aS)-3,7-Bis[3,5-bis(tert-butyl)-4-methoxyphenyl]-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee)
C₅₁H₆₇O₆P; FW: 807; off white to light yellow powdr.
Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

PHOSPHORUS (Compounds)


15-0578

NEW

(11aS)-3,7-Bis[3,5-bis(tert-butyl)phenyl]-10,11,12,13-tetrahydro-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee)

$C_{46}H_{55}O_4P$; FW: 690.9; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.



25mg


15-8592

NEW

(11aR)-3,7-Bis[3,5-bis(tert-butyl)phenyl]-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee)

$C_{48}H_{63}O_4P$; FW: 747; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.



25mg

100mg


15-8593

NEW

(11aS)-3,7-Bis[3,5-bis(tert-butyl)phenyl]-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee)

$C_{48}H_{63}O_4P$; FW: 747; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.



25mg

100mg


15-8575

NEW

(11aR)-3,7-Bis(4-chlorophenyl)-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%

$C_{33}H_{29}Cl_2O_4P$; FW: 591.5; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.



25mg

100mg

15-3360


NEW

1-[2,6-Bis(cyclohexyloxy)phenyl]-2-(di-tert-butylphosphaneyl)-1H-imidazole, min. 95% (2179272-79-6)

$C_{29}H_{45}N_2O_2P$; FW: 484.66; white to off-white to pale yellow solid

air sensitive

Note: Sold in collaboration with Solvias for research purposes only. Patents: EP 1625133, US 2007/0123707.




100mg

500mg

Technical Notes:


- Ligand for the Pd-catalyzed Barbier-Negishi coupling of secondary alkyl bromides with aryl and alkenyl triflates and nonaflates.
- Used in Pd-catalyzed enantioselective synthesis of β 2- and β 3- amino acids via asymmetric and regiodivergent C-H Functionalization of Boc-1,3-oxazinanones.



Tech. Note (1)
Ref. (1)

PHOSPHORUS (Compounds)

15-3360 **1-[2,6-Bis(cyclohexyloxy)phenyl]-2-(di-tertbutylphosphaneyl)-1H-imidazole, min. 95%**
(continued) (2179272-79-6)



References:


1. *Angew. Chem. Int. Ed.* **2018**, 57, 1982.
2. *Nat. Catal.* **2019**, 2, 882.

15-8508**NEW**

N-[(1S)-1-[3',5'-Bis(1,1-dimethylethyl)[1,1'-biphenyl]-2-yl]-2-(diphenylphosphino)ethyl]-3,5-bis(trifluoromethyl)-benzamide, 95%
 (2089424-11-1)

C₄₃H₄₂F₆NOP; FW: 733.8; white to off-white powdr.
air sensitive, moisture sensitive


Note: Sold in collaboration with Daicel for research purposes only.




50mg

Technical Notes:


1. Organocatalyst used in aza-[4+2] cycloaddition of N-sulfonyl-1-aza-1,3-dienes with vinyl ketones.
2. Catalyst for the enantioselective [3+2] cycloadditions of γ-substituted allenates with β-perfluoroalkyl enones.
3. Phosphine catalyst for difunctionalization of β-fluoroalkyl α, β-Enones to generate β-amino α-diazo carbonyl compounds.



Tech. Note (1)
 Ref. (1)



Tech. Note (2)
 Ref. (2)



Tech. Note (3)
 Ref. (3)

References:

1. *Org. Lett.* **2017**, 19, 1710.
2. *Chem. Sci.*, **2017**, 8, 4660.
3. *Angew. Chem. Int. Ed.* **2018**, 57, 15787.


PHOSPHORUS (Compounds)

15-8372

NEW


[S(R)]-N-[(1S)-1-[3',5'-Bis(1,1-dimethylethyl)-4'-methoxy[1,1'-biphenyl]-2-yl]-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95% (1936438-26-4)
 $C_{39}H_{50}NO_2PS$; FW: 627.9; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg




Technical Notes:


1. Organocatalyst used in enantioselective intermolecular cross Rauhut-Currier reactions of activated alkenes with acrolein.
2. Used in phosphine-catalyzed asymmetric umpolung addition of trifluoromethyl ketimines to Morita-Baylis-Hillman carbonates.
3. Organocatalyst for [4+2] annulation of electron-deficient dienes and vinyl ketones.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)

References:


1. *Chem. Commun.*, **2016**, 52, 7612.
2. *Angew. Chem. Int. Ed.* **2016**, 55, 13316.
3. *Adv. Synth. Catal.* **2018**, 360, 682.

15-8438

NEW


[S(R)]-N-[(R)-[3,5-Bis(1,1-dimethylethyl)-4-methoxyphenyl][2-(dicyclohexylphosphino)phenyl]methyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{39}H_{62}NO_2PS$; FW: 640.0; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg



Technical Notes:

1. Ligand for the Pd-catalyzed enantioselective reductive Heck reactions to generate 3,3-disubstituted 2,3-dihydrobenzofuran.
2. Ligand for the Pd-catalyzed enantioselective carboidination of olefin-tethered aryl iodides.
3. Ligand for the Pd-catalyzed enantioselective dicarbonylization of unactivated alkenes via Heck-Suzuki coupling reaction.




Tech. Note (1)
Ref. (1)


PHOSPHORUS (Compounds)

15-8438
(continued)

[S(R)]-N-[(R)-[3,5-Bis(1,1-dimethylethyl)-4-methoxyphenyl][2-(dicyclohexylphosphino)phenyl]methyl]-N,2-dimethyl-2-propanesulfonamide, 95%



X = O; N-Boc

Tech. Note (2)
Ref. (2)

X = C, O, N; R = Aryl, Alkyl, Vinyl

Tech. Note (3)
Ref. (3-4)

References:

1. *Angew. Chem. Int. Ed.* **2018**, 57, 10373.
2. *J. Am. Chem. Soc.* **2019**, 141, 8110.
3. *Angew. Chem. Int. Ed.* **2019**, 58, 14653.
4. *Angew. Chem. Int. Ed.* **2020**, 59, 2769.

15-8436

NEW

[S(R)]-N-[(S)-[3,5-Bis(1,1-dimethylethyl)-4-methoxyphenyl][2-(dicyclohexylphosphino)phenyl]methyl]-N,2-dimethyl-2-propanesulfonamide, 95%

50mg

C₃₈H₆₂NO₂PS; FW: 640.0; white to off-white powd.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.

Technical Note:

1. See 15-8438 (page 48)

15-8434

NEW


[S(R)]-N-[(R)-[3,5-Bis(1,1-dimethylethyl)-4-methoxyphenyl][2-(dicyclohexylphosphino)phenyl]methyl]-2-methyl-2-propanesulfonamide, 95%

50mg

C₃₈H₆₀NO₂PS; FW: 625.9; white to off-white powd.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



15-8320

NEW


[S(R)]-N-[(R)-[3,5-Bis(1,1-dimethylethyl)-4-methoxyphenyl][2-(diphenylphosphino)phenyl]methyl]-2-methyl-2-propanesulfonamide, 95% (1616688-64-2)

50mg

C₃₈H₄₈NO₂PS; FW: 613.80; white to off-white powd.


air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



Technical Note*

1. Ligand for the enantioselective Au-catalyzed [3+3] cycloaddition reaction of 2-(1-alkynyl)-alk-2-en-1-ones with nitrones.

Tech. Note (1)
Ref. (1)

References:

1. *Angew. Chem. Int. Ed.* **2014**, 53, 4350.

PHOSPHORUS (Compounds)

15-8318	[S(R)]-N-[(S)-[3,5-Bis(1,1-dimethylethyl)-4-methoxyphenyl][2-(diphenylphosphino)phenyl]methyl]-2-methyl-2-propanesulfonamide, 95% (1616688-63-1) $C_{38}H_{48}NO_2PS$; FW: 613.80; white to off-white pwdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.	50mg
---------	---	------


NEW

Technical Note:

- See 15-8320 (page 49)


15-8464	[S(R)]-N-[(S)-3,5-Bis(dimethylethyl)phenyl][2-(dicyclohexylphosphanyl)phenyl]-N,2-dimethyl-2-propanesulfonamide, 95% $C_{38}H_{60}NOPS$; FW: 607.9; light-yellow to yellow pwdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.	50mg
---------	---	------

NEW




15-8462	[S(R)]-N-[(S)-3,5-Bis(dimethylethyl)phenyl][2-(dicyclohexylphosphanyl)phenyl]-2-methyl-2-propanesulfonamide, 95% $C_{37}H_{58}NOPS$; FW: 595.9; white to off-white pwdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.	50mg
---------	---	------

NEW




15-8578	(11aR)-3,7-Bis((4-(1,1-dimethylethyl)phenyl)-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee) $C_{47}H_{47}O_4P$; FW: 634.8; off white to light yellow pwdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
---------	--	---------------

NEW




15-8562	(11aR)-3,7-Bis(3,5-dimethylphenyl)-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee) $C_{37}H_{39}O_4P$; FW: 578.7; off white to light yellow pwdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
---------	--	---------------

NEW



15-8563	(11aS)-3,7-Bis(3,5-dimethylphenyl)-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee) $C_{37}H_{39}O_4P$; FW: 578.7; off white to light yellow pwdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
---------	--	---------------

NEW




PHOSPHORUS (Compounds)

15-0165

NEW

(2R,3R)-(+)-Bis(diphenylphosphino)butane,
min. 98%, (R,R)-CHIRAPHOS (74839-84-2)
(C₆H₅)₂PCH(CH₃)CH(CH₃)P(C₆H₅)₂; FW: 426.48;
white xtl.; m.p. 104-109 °C
air sensitive




250mg

1g

5g

Technical Notes:

1. A bidentate chiral phosphine ligand used in chiral catalysis.^{1,2}
2. Also used in metal catalysis chemistry.³


Fe(acac)₃ (cat.)

RMgBr

ZnX₂.TMEDA

THF, 23-85°C

DG = directing group,
chiral ligand 15-0165,
up to 46% ee



R = alkyl, aryl,
up to 94% yield.

Tech. Note (1)

Ref. (1)

References:


1. *Organometallics*, **2017**, 4979-4989.
2. *J. Am. Chem. Soc.* **1996**, 118, 4723-4724.
3. *J. Am. Chem. Soc.* **2008**, 130, 12874-12875.

15-8250

NEW

amp

Bis(3,5-di-tert-butylphenyl)cyclohexylphosphine,
min. 98% CyTyrannoPhos (2097604-67-4)
C₃₄H₅₃P; FW: 492.76; white to off-white solid
air sensitive, moisture sensitive








250mg

1g

Technical Note:

1. Ligand used for Ni-catalysed Suzuki coupling of reactions of benzylic acetals.



Tech. Note (1)

Ref. (1)


References:

1. *Nat. Chem.*, **2017**, 9, 779.

15-1269

NEW

(-)-4,6-Bis((S)-(2-methoxyphenyl)(phenyl)
phosphinyl)dibenzo[b,d]furan, min. 97%
(2119686-75-6)
C₃₈H₃₀O₃P₂; FW: 596.60; white to off-white solid
air sensitive, moisture sensitive, (store cold)
Note: Sold under license of WO2017/191310 for
research purposes only.







25mg

100mg

Technical Note:

1. P-Chirogenic Xantphos ligand used for the synthesis and application in Rhodium-catalyzed asymmetric hydrogenation.



Tech. Note (1)

Ref. (1)

References:

1. *ACS Catal.*, **2017**, 7, 6162-6169.


PHOSPHORUS (Compounds)

15-8524

NEW

(11bR)-2,6-Bis[2-methoxyphenyl-5-(tert-butyl)-4-hydroxy-4-oxide-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin, 98%, (99% ee)
 $C_{42}H_{41}O_6P$; FW: 672.7; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

15-8525

NEW

(11bS)-2,6-Bis[2-methoxyphenyl-5-(tert-butyl)-4-hydroxy-4-oxide-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin, 98%, (99% ee)

$C_{42}H_{41}O_6P$; FW: 672.7; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

15-3365


NEW

2-[Bis(1-methylethyl)phosphino]-1-phenyl-1H-pyrrole, min. 95% (1257847-61-2)

$C_{16}H_{22}NP$; FW: 259.33; white to off-white to pale yellow solid




air sensitive

Note: Sold in collaboration with Solvias for research purposes only. Patents: EP 1625133, US 2007/0123707.

100mg
500mg

Technical Notes:

1. Ligand for the Ru-catalyzed direct amination of secondary alcohols using ammonia.
2. Ligand for the Pd-catalyzed α - and β -arylation of acyclic N-Boc amines.
3. Used in Pd-catalyzed enantioselective synthesis of β^2 - and β^3 - amino acids via asymmetric and regiodivergent C-H Functionalization of Boc-1,3-oxazinanes.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)

References:


1. *Angew. Chem. Int. Ed.* **2010**, 49, 8130.
2. *Angew. Chem. Int. Ed.* **2014**, 53, 2678.
3. *Nat. Catal.* **2019**, 2, 882.

PHOSPHORUS (Compounds)

15-1285



NEW

1,3-Bis[(1R)-1-(1-naphthalenyl)ethyl]-2,3-dihydro-1H-1,3,2-diazaphosphol-2-yl trifluoromethanesulfonate, min. 98% (2377935-66-3)
 $C_{27}H_{24}N_2O_3PS$; FW: 544.53; light beige solid
moisture sensitive
 Note: Patents: PCT #CA2018/051078

50mg
250mg

Technical Notes:

1. Asymmetric Imine Hydroboration Catalyzed by Chiral Diazaphospholenes.
2. Enantioselective Imine Reduction Catalyzed by Phosphonium Ions.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Phosphonium-catalyzed
asymmetric imine reductionBorane or silane reductants,
grams of catalyst in 3 steps!17 examples > 90:10 e.r.
> 98% conversion


References:

1. *Angew. Chem. Int. Ed.* **2017**, 56, 16660.
2. *J. Am. Chem. Soc.* **2019**, 141, 36, 14083-14088.

15-1292


NEW

1,3-Bis[(1S)-1-(1-naphthalenyl)ethyl]-2,3-dihydro-1H-1,3,2-diazaphosphol-2-yl trifluoromethanesulfonate, min. 98%
 $C_{27}H_{24}N_2O_3PS$; FW: 544.53; light beige solid
air sensitive, moisture sensitive
 Note: Patents: PCT #CA2018/051078

50mg
250mg

Technical Note:

1. Used as an organocatalyst for enantioselective imine reduction.

Tech. Note (1)
Ref. (1)


References:

1. *J. Am. Chem. Soc.* **2019**, 141, 14083.

15-8585

NEW

(11aS)-3,7-Bis(4-nitrophenyl)-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 95%, (99% ee)
 $C_{33}H_{29}N_2O_8P$; FW: 612.6; off white to light yellow powdr.
 Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

15-8584

NEW

(11aR)-3,7-Bis(4-nitrophenyl)-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee)
 $C_{33}H_{29}N_2O_8P$; FW: 612.6; off white to light yellow powdr.
 Note: Sold in collaboration with Daicel for research purposes only.


25mg
100mg

PHOSPHORUS (Compounds)

15-8460

NEW

[S(R)]-N-[(R)-3,5-Bis(trifluoromethyl)phenyl]
[2-(dicyclohexylphosphanyl)phenyl]-N,2-di-
methyl-2-propanesulfonamide, 95%
C₃₂H₄₂F₆NOPS; FW: 633.7; gray powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for
research purposes only.




50mg

15-8516

NEW


[S(R)]-N-[(S)-[3,5-Bis(trifluoromethyl)phenyl]
[2-(diphenylphosphino)phenyl]meth-
yl]-2-methyl-2-propanesulfonamide, 95%
(2262535-73-7)
C₃₁H₂₈F₆NOPS; FW: 607.6; white to off-white
powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for
research purposes only.



50mg

Technical Note:

- Phosphine ligand for the Ag(I)-catalyzed enantioselective [3+2] cycloaddition of azomethine ylides.



Tech. Note (1)
Ref. (1)


References:

- Org. Biomol. Chem.*, **2019**, *17*, 1395.

15-8581

NEW

(11aS)-3,7-Bis(4-(trifluoromethyl)phenyl)-
10,11,12,13-tetrahydro-10,10,13,13-tetrameth-
yl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg]
[1,3,2]dioxaphosphocin, 95%, (99% ee)
C₃₅H₂₉F₆O₄P; FW: 658.6; off white to light yellow
powdr.
Note: Sold in collaboration with Daicel for
research purposes only.




25mg
100mg

15-7555









NEW

Bis(2,4,4-trimethylpentyl)phosphinic acid,
min. 85%, CYTOP® 501 (83411-71-6)
C₁₆H₃₅O₂P; FW: 290.42; Yellow liq.; f.p. >2226 °F;
d. 0.92
Note: Sold in collaboration with Solvay.



100g
500g

PHOSPHORUS (Compounds)

15-8486 NEW	<p>[S(R)]-N-[(R)-[2-(Diadamantanphosphino)phenyl]phenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95% $C_{38}H_{52}NOPS$; FW: 601.9; white to off-white powdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>		50mg
15-8484 NEW	<p>[S(R)]-N-[(S)-[2-(Diadamantanphosphino)phenyl]phenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95% $C_{38}H_{52}NOPS$; FW: 601.9; white to off-white powdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>		50mg
15-8488 NEW	<p>[S(R)]-N-[(R)-[2-(Diadamantanphosphino)phenyl]phenylmethyl]-2-methyl-2-propanesulfonamide, 95% $C_{37}H_{50}NOPS$; FW: 587.8; white to off-white powdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>		50mg
15-8490 NEW	<p>[S(R)]-N-[(S)-[2-(Diadamantanphosphino)phenyl]phenylmethyl]-2-methyl-2-propanesulfonamide, 95% $C_{37}H_{50}NOPS$; FW: 587.8; white to off-white powdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>		50mg
15-8570 NEW	<p>(11aS)-3,7-Di-9-anthracenyl-10,11,12,13-tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 95%, (99% ee) $C_{48}H_{30}O_5P$; FW: 722.8; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.</p>		25mg 100mg
15-1044 NEW	<p>Di-t-butylphosphine (50% in Toluene) (819-19-2) $[(CH_3)_3C]_2PH$; FW: 146.22; colorless liq.; d. 0.79 <i>air sensitive</i></p>		25g 100g
15-8454 NEW	<p>[S(R)]-N-[(1S)-1-[2-(Dicyclohexylphosphanyl)phenyl]-cyclohexyl]-2-methyl-2-propanesulfonamide, 95% $C_{28}H_{48}NOPS$; FW: 489.7; white to off-white powdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>		50mg
15-8452 NEW	<p>[S(R)]-N-[(1S)-1-[2-(Dicyclohexylphosphanyl)phenyl]-2,2-dimethylpropyl]-2-methyl-2-propanesulfonamide, 95% $C_{27}H_{46}NOPS$; FW: 463.7; white to off-white powdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>		50mg

PHOSPHORUS (Compounds)

15-8450


NEW

[S(R)]-N-[(R)-[2-(Dicyclohexylphosphanyl)phenyl](4-methoxyphenyl)methyl]-N,2-dimethyl-2-propanesulfinamide, 95%

$C_{31}H_{46}NO_2PS$; FW: 511.7; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8448

NEW

[S(R)]-N-[(S)-[2-(Dicyclohexylphosphanyl)phenyl](4-methoxyphenyl)methyl]-N,2-dimethyl-2-propanesulfinamide, 95%

$C_{31}H_{46}NO_2PS$; FW: 511.7; yellow

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.

15-8458


NEW

[S(R)]-N-[(R)-[2-(Dicyclohexylphosphanyl)phenyl](4-(tert-butyl)phenyl)methyl]-N,2-dimethyl-2-propanesulfinamide, 95%

$C_{34}H_{52}NOPS$; FW: 553.8; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8432


NEW

[S(R)]-N-[(S)-1-[2-(Dicyclohexylphosphino)phenyl]-(1,1'-biphenyl)methyl]-2-methyl-2-propanesulfinamide, 95%

$C_{35}H_{46}NOPS$; FW: 559.8; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8444


NEW

[S(R)]-N-[(R)-[2-(Dicyclohexylphosphino)phenyl]-2-naphthalenylmethyl]-N,2-dimethyl-2-propanesulfinamide, 95%

$C_{34}H_{46}NOPS$; FW: 547.8; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8446


NEW

[S(R)]-N-[(R)-[2-(Dicyclohexylphosphino)phenyl]-2-naphthalenylmethyl]-2-dimethyl-2-propanesulfinamide, 95%

$C_{33}H_{44}NOPS$; FW: 533.7; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.




50mg

PHOSPHORUS (Compounds)

15-8440

NEW


[S(R)]-N-[(S)-[2-(Dicyclohexylphosphino)phenyl]-1-naphthalenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95% (2241598-33-2)
 $C_{34}H_{46}NOPS$; FW: 547.8; pink powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



50mg

Technical Note:

- Ligand for the Pd-catalyzed enantioselective reductive Heck reactions to generate 3,3-disubstituted 2,3-dihydrobenzofuran.



Tech. Note (1)
Ref. (1)

References:

- Angew. Chem. Int. Ed.* **2018**, 57, 10373.

15-8442


NEW

[S(R)]-N-[(S)-[2-(Dicyclohexylphosphino)phenyl]-2-naphthalenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{34}H_{46}NOPS$; FW: 547.8; light-yellow to yellow powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg

15-8430

[S(R)]-N-[(R)-[2-(Dicyclohexylphosphino)phenyl]phenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{30}H_{44}NOPS$; FW: 497.7; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8428

NEW

[S(R)]-N-[(S)-[2-(Dicyclohexylphosphino)phenyl]phenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95% (2241598-32-1)
 $C_{30}H_{44}NOPS$; FW: 497.7; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg

Technical Note:

- See 15-8440 (page 57)

15-8424

NEW

[S(R)]-N-[(S)-[2-(Dicyclohexylphosphino)phenyl]phenylmethyl]-2-methyl-2-propanesulfonamide, 95% (1595319-99-5)
 $C_{29}H_{42}NOPS$; FW: 483.7; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg

Technical Note:


- See 15-8426 (page 58)

PHOSPHORUS (Compounds)

15-8426

NEW


[S(R)]-N-[(R)-[2-(Dicyclohexylphosphino)phenyl]phenylmethyl]-2-methyl-2-propanesulfonamide, 95%
 $C_{29}H_{42}NOPS$; FW: 483.7; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.




50mg

Technical Notes:

- Ligand for the Pd-catalyzed enantioselective reductive Heck reactions to generate 3,3-disubstituted 2,3-dihydrobenzofuran.
- Ligand for the Pd-catalyzed enantioselective Heck reaction of aryl triflates and alkynes.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)


References:

- Angew. Chem. Int. Ed. **2018**, 57, 10373.
- J. Am. Chem. Soc. **2019**, 141, 19246.

15-8382

NEW

[S(R)]-N-[(1S)-1-(2',6'-Diisopropyl)-(1,1'-biphenyl)-2-yl]-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95%
 $C_{36}H_{44}NOPS$; FW: 569.8; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.




50mg

15-8410

NEW


[S(R)]-N-[(R)-[4-(1,1-Dimethylethyl)phenyl][5-(diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl]methyl]-2-methyl-2-propanesulfonamide, 95% (2160535-58-8)
 $C_{42}H_{46}NO_2PS$; FW: 659.9; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



50mg

Technical Note:

- Ligand for the Au-catalyzed asymmetric intramolecular cyclization of N-allenamides for the synthesis of chiral tetrahydrocarbolines.



Tech. Note (1)
Ref. (1)

PHOSPHORUS (Compounds)

15-8410 [S(R)]-N-[(R)-[4-(1,1-Dimethylethyl)phenyl][5-(diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl]methyl]-2-methyl-2-propanesulfonamide, 95% (2160535-58-8)
(continued)

References:

1. *Angew. Chem. Int. Ed.* **2017**, 56, 15905.

15-8408

NEW

[S(R)]-N-[(S)-[4-(1,1-Dimethylethyl)phenyl][5-(diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl]methyl]-2-methyl-2-propanesulfonamide, 95%
(2160535-59-9)

C₄₂H₄₆NO₂PS; FW: 659.9; white to off-white powdr.


air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.


50mg

Technical Notes:

1. Ligand for the Au-catalyzed asymmetric intramolecular cyclization of N-allenamides for the synthesis of chiral tetrahydrocarbolines.
2. Ligand for the Pd-catalyzed enantioselective intermolecular denitrogenative cyclization of benzotriazoles with allenes and N-allenamides.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)

References:

1. *Angew. Chem. Int. Ed.* **2017**, 56, 15905.
2. *Angew. Chem. Int. Ed.* **2019**, 58, 11444.

15-4665


NEW

2-[(1,1-Dimethylethyl)thio]-N-[2-(diphenylphosphino)ethyl]-ethanamine (1883429-96-6)

C₂₀H₂₈NPS; FW: 345.48; yellow oil

air sensitive


Note: Patents: US20170088571A1; U.S. Provisional Patent Application No. 62/827,627.



50mg

Technical Note:

1. SNP ligand used in Ru-catalyzed selective ester hydrogenation.



Tech. Note (1)
Ref. (1)

References:

1. *Org. Process Res. Dev.* **2020**, 24, 415.


PHOSPHORUS (Compounds)

15-4670

NEW

2-[(1,1-Dimethylethyl)thio]-N-[2-(diphenylphosphinyl)ethyl]-ethanamine $C_{20}H_{28}NOP$; FW: 361.48; yellow oil

Note: Patents: US20170088571A1; U.S. Provisional Patent Application No. 62/827,627.



250mg

Technical Note:


1. See 15-4665 (page 59)

15-7605

NEW

4,8-Dimethyl-2-phosphabicyclo[3.3.1]nonane, min.**90%, CYTOP® 170 (328952-85-8)** $C_{10}H_{18}P$; FW: 170.2; Colorless to pale yellow liq.; m.p. -22 °C; b.p. 120 °C; f.p. 221 °F; d. 0.95*air sensitive, moisture sensitive*

Note: Sold in collaboration with Solvay.




25g

100g

500g

15-1239

NEW



(9,9-Dimethyl-9H-xanthene-4,5-diyl)bis([o-methoxyphenyl]phenylphosphine), min. 97% (2119686-35-8) $C_{41}H_{36}O_3P_2$; FW: 638.67; white to off-white solid*air sensitive, moisture sensitive, (store cold)*
Note: Sold under license of WO2017/191310 for research purposes only.

25mg

100mg

Technical Notes:

1. P-Chirogenic Xantphos ligand used for the synthesis and application in Rhodium-catalyzed asymmetric hydrogenation.
2. Diastereo- and enantioselective copper catalyzed hydroallylation of disubstituted cyclopropenes.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)

References:


1. *ACS Catal.*, **2017**, 7, 6162–6169
2. *Chem. Sci.*, **2018**, 9, 6503–6508.

15-8245

NEW

((((9,10-Dioxo-9,10-dihydroanthracene-2,6-diyl)bis(oxy))bis(propane-3,1-diyl))bis(phosphonic acid), 98%, 2,6-DPPEAQ $C_{20}H_{22}O_{10}P_2$; FW: 484.33; light yellow to beige powdr.

Note: U.S. Patents: 62/628,599, 62/740,526.

Patent Application: PCT/US19/17479. Product sold under, use subject to, terms and conditions of Limited Use License at www.strem.com/harvard5

1g

5g

Technical Note:

1. A highly stable phosphonate-functionalized anthraquinone used as a redox-active material in a negative potential electrolyte (negolyte) for aqueous redox flow batteries operating at nearly neutral pH.

References:


1. *Adv. Energy Mater.* **2019**, 9(12), 1900039.
2. *Joule*, **2018**, 2(9), 1894–1906.

PHOSPHORUS (Compounds)

15-8396

NEW

[S(R)]-N-[(1S)-2-(Diphenylphosphino)-1-[2-(diadamantanphosphanyl)phenyl]ethyl]-2-methyl-2-propanesulfonamide, 95%
C₄₄H₅₇NOP₂S; FW: 709.9; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




50mg

15-8394

NEW

[S(R)]-N-[(1S)-2-(Diphenylphosphino)-1-[2-(dicyclohexylphosphanyl)phenyl]ethyl]-2-methyl-2-propanesulfonamide, 95%
C₃₆H₄₉NOP₂S; FW: 605.8; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




50mg

15-8398

NEW

[S(R)]-N-[(1S)-2-(Diphenylphosphino)-1-[2-(diphenylphosphino)-4,5-dimethoxyphenyl]ethyl]-2-methyl-2-propanesulfonamide, 95%
C₃₈H₄₁NO₃P₂S; FW: 653.8; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




50mg

15-8392

NEW


[S(R)]-N-[(1S)-2-(Diphenylphosphino)-1-[2-(diphenylphosphino)phenyl]ethyl]-2-methyl-2-propanesulfonamide, 95% (1824731-39-6)
C₃₆H₃₇NOP₂S; FW: 593.7; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




50mg

Technical Notes:

1. Chiral sulfonamide bisphosphine catalysts for highly enantioselective intermolecular cross-Rauhut-Currier reactions.
2. Used in phosphine-catalyzed asymmetric intermolecular cross-vinyllogous Rauhut-Currier Reactions of vinyl Ketones with para-quinone methides.
3. Used in the enantioselective phosphine-catalyzed allylic alkylations of mix-indene with MBH carbonates.
4. Ligand for the Pd-catalyzed fast and living polymerizations of various diazoacetate monomers under mild conditions.
5. Ligand for the direct Pd-catalyzed cross-coupling reaction of easily accessible secondary phosphine oxides and aryl bromides, which provides rapid access to P-chiral phosphine oxides.




Tech. Note (1)
Ref. (1)




Tech. Note (2)
Ref. (2)

PHOSPHORUS (Compounds)

15-8392 [S(R)]-N-[(1S)-2-(Diphenylphosphino)-1-[2-(diphenylphosphino)phenyl]ethyl]-2-methyl-2-propanesulfonamide, 95% (1824731-39-6)




Tech. Note (3)
Ref. (3)



Tech. Note (4)
Ref. (4)

R = Alkyl or Aryl



Tech. Note (5)
Ref. (5)


References:

1. *Angew. Chem. Int. Ed.* **2015**, *54*, 14853.
2. *ACS Catal.* **2017**, *7*, 2805.
3. *Org. Lett.* **2017**, *19*, 6080.
4. *J. Am. Chem. Soc.* **2018**, *140*, 17773.
5. *J. Am. Chem. Soc.* **2019**, *141*, 20556.

15-8358

NEW

[S(R)]-N-[(R)-[6-(Diphenylphosphino)benzo[d][1,3]dioxol-5-yl]-1-naphthalenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%
C₃₅H₃₄NO₃PS; FW: 579.7; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




50mg

15-8354

NEW

[S(R)]-N-[(R)-[6-(Diphenylphosphino)benzo[d][1,3]dioxol-5-yl]-1-naphthalenylmethyl]-2-methyl-2-propanesulfonamide, 95%
C₃₄H₃₂NO₃PS; FW: 565.7; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




50mg

15-8356

NEW

[S(R)]-N-[(S)-[6-(Diphenylphosphino)benzo[d][1,3]dioxol-5-yl]-1-naphthalenylmethyl]-2-methyl-2-propanesulfonamide, 95%
C₃₄H₃₂NO₃PS; FW: 565.7; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




50mg

15-8352

NEW



[S(R)]-N-[(R)-[6-(Diphenylphosphino)benzo[d][1,3]dioxol-5-yl]phenylmethyl]-2-methyl-2-propanesulfonamide, 95%
C₃₀H₃₀NO₃PS; FW: 515.6; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.



50mg


PHOSPHORUS (Compounds)

15-8350 NEW	<p>[S(R)]-N-[(S)-[6-(Diphenylphosphino)benzo[d][1,3]dioxol-5-yl]phenylmethyl]-2-methyl-2-propanesulfonamide, 95% $C_{30}H_{30}NO_3PS$; FW: 515.6; white to off-white powd. <i>air sensitive, light sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>	50mg
15-8360 NEW	<p>[S(R)]-N-[(S)-[2-(Diphenylphosphino)-4,5-dimethoxyphenyl](4-methoxyphenyl)methyl]-N,2-dimethyl-2-propanesulfonamide, 95% $C_{33}H_{38}NO_4PS$; FW: 575.7; light-yellow to yellow powd. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>	50mg
15-8422 NEW	<p>[S(R)]-N-[(1R)-1-[5-(Diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl]-2,2-dimethylpropyl]-2-methyl-2-propanesulfonamide, 95% (2162939-92-4) $C_{36}H_{42}NO_2PS$; FW: 583.8; white to off-white powd. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>	25mg



Technical Note:

1. Ligand for the Pd-catalyzed enantioselective arylation of general sulfonate anions.




Tech. Note (1)
Ref. (1)

References:

1. *J. Am. Chem. Soc.* **2018**, *140*, 3467.

15-8420 NEW	<p>[S(R)]-N-[(1S)-1-[5-(Diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl]-2,2-dimethylpropyl]-2-methyl-2-propanesulfonamide, 95% (2162939-89-9) $C_{36}H_{42}NO_2PS$; FW: 583.8; white to off-white powd. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>	50mg
15-8400	<p>[S(R)]-N-[(R)-[5-(Diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl](4-methoxyphenyl)methyl]-2-methyl-2-propanesulfonamide, 95% (2160535-56-6) $C_{38}H_{40}NO_3PS$; FW: 633.8; white to off-white powd. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>	50mg




Technical Note:

1. Ligand for the Au-catalyzed asymmetric intramolecular cyclization of N-allenamides for the synthesis of chiral tetrahydrocarbolines.

PHOSPHORUS (Compounds)

15-8400
(continued)

[S(R)]-N-[(R)-[5-(Diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl](4-methoxyphenyl)methyl]-2-methyl-2-propanesulfonamide, 95% (2160535-56-6)

Tech. Note (1)
Ref. (1)

References:

- 1.
- Angew. Chem. Int. Ed.*
- 2017**
- , 56, 15905.

15-8402

NEW

[S(R)]-N-[(S)-[5-(Diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl](4-methoxyphenyl)methyl]-2-methyl-2-propanesulfonamide, 95% (2160535-57-7)
 $C_{38}H_{40}NO_3PS$; FW: 633.8; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg

Technical Note:


1. See 15-8400 (page 63)

15-8418

NEW

[S(R)]-N-[(S)-(phenyl)[5-(Diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl](phenyl)methyl]-2-methyl-2-propanesulfonamide, 95%
 $C_{38}H_{38}NO_2PS$; FW: 603.8; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg




15-8512

NEW

N-[(1S)-2-(Diphenylphosphino)-1-(2-(diphenylphosphino)phenyl)ethyl]-3,5-bis(trifluoromethyl)-benzamide, 95%
 $C_{41}H_{31}F_6NOP_2$; FW: 729.6; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg



PHOSPHORUS (Compounds)


15-4645

NEW

N-[2-(Diphenylphosphino)ethyl]-2-[(4-methylphenyl)thio]ethanamine (1883429-99-9)C₂₃H₂₆NPS; FW: 379.50; yellow oil*air sensitive*

Note: Patents: US20170088571A1; U.S. Provisional


Patent Application No. 62/827,627.



100mg

Technical Note:

1. SNP ligand used in Ru-catalyzed selective ester hydrogenation.



Tech. Note (1)
Ref. (1)

References:

1. *Org. Process Res. Dev.* **2020**, 24, 415.


15-4625

NEW

N-[2-(Diphenylphosphino)ethyl]-2-(methylthio)ethanamine (1802182-42-8)C₁₇H₂₂NPS; FW: 303.40; yellow oil*air sensitive*

Note: Patents: US20170088571A1; U.S. Provisional

Patent Application No. 62/827,627.



50mg

Technical Note:


1. See 15-4645 (page 65)

15-8378

NEW

[S(R)]-N-[(1S)-1-[(Diphenylphosphino)methyl]-2,2-dimethylpropyl]-2-methyl-2-propanesulfonamide, 95% (1853342-54-7)C₂₂H₃₂NOPS; FW: 389.5; white to off-white powdr.*air sensitive, moisture sensitive*


Note: Sold in collaboration with Daicel for research purposes only.



100mg

Technical Notes:

1. Used in phosphine-catalyzed asymmetric intermolecular cross-vinylous Rauhut-Currier Reactions of vinyl Ketones with para-quinone methides.
2. Used in the enantioselective phosphine-catalyzed allylic alkylations of mix-indene with MBH carbonates.




Tech. Note (1)
Ref. (1)

PHOSPHORUS (Compounds)

15-8378
(continued)

[S(R)]-N-[(1S)-1-[(Diphenylphosphino)methyl]-2,2-dimethylpropyl]-2-methyl-2-propanesulfonamide, 95% (1853342-54-7)

Tech. Note (2)
Ref. (2)


References:

1. *ACS Catal.* **2017**, 7, 2805.
2. *Org. Lett.* **2017**, 19, 6080.

15-8376

NEW




[S(R)]-N-[(1S)-1-[(Diphenylphosphino)methyl]-2-methylpropyl]-2-methyl-2-Propanesulfonamide, 95% (1803239-46-4)
 C₂₁H₃₀NOPS; FW: 375.5; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



50mg

Technical Notes:

1. Sulfonamide phosphine organocatalyst for the enantioselective intramolecular Rauht-Currier reaction.
2. Used in phosphine-catalyzed asymmetric intermolecular cross-vinylous Rauht-Currier Reactions of vinyl Ketones with para-quinone methides.
3. Used in the enantioselective phosphine-catalyzed allylic alkylations of mix-indene with MBH carbonates.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)

1. *Angew. Chem. Int. Ed.* **2015**, 54, 6874.
2. *ACS Catal.* **2017**, 7, 2805.
3. *Org. Lett.* **2017**, 19, 6080.

PHOSPHORUS (Compounds)

15-8326


NEW

[S(R)]-N-[(S)-[2-(Diphenylphosphino)phenyl]-1,1'-biphenyl)methyl]-2-methyl-2-propanesulfonamide, 95%

C₃₅H₃₄NOPS; FW: 547.7; white to off-white powdr.

air sensitive, moisture sensitive


Note: Sold in collaboration with Daicel for research purposes only.



100mg

Technical Note:

- Ligand for the enantioselective Au-catalyzed [3+3] cycloaddition reaction of 2-(1-alkynyl)-alk-2-en-1-ones with nitrones.



**Tech. Note (1)
Ref. (1)**

References:


- Angew. Chem. Int. Ed. **2014**, 53, 4350.

15-8348

NEW

[S(R)]-N-[(1R)-1-[2-(Diphenylphosphino)phenyl]-2,2-dimethylpropyl]-N,2-dimethyl-2-propanesulfonamide, 95% (2049042-08-0)


C₂₈H₃₆NOPS; FW: 465.6; white to off-white powdr.
moisture sensitive



50mg

Technical Note:



- Ligand used in the Copper(I)-catalyzed asymmetric intermolecular [3+2] cycloaddition of azomethine ylides with α -Trifluoromethyl α,β -unsaturated esters.



PHOSPHORUS (Compounds)

15-8330
(continued)

[S(R)]-N-[(1R)-1-[2-(Diphenylphosphino)phenyl]-2,2-dimethylpropyl]-2-methyl-2-propane-sulfonamide, 95% (1906918-20-4)

Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)


References:

1. *Angew. Chem. Int. Ed.* **2016**, 55, 6324.
2. *ACS Catal.* **2017**, 7, 210.
3. *Org. Biomol. Chem.*, **2019**, 17, 1395.

15-8514

NEW




N-[(1S)-2-(Diphenylphosphino)-1-phenylethyl]-3,5-bis(trifluoromethyl)-benzamide, 95% (1853342-57-0)
 $C_{29}H_{22}F_6NOP$; FW: 545.5; white to off-white powd.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.




50mg

Technical Notes:

1. Organocatalyst used in aza-[4+2] cycloaddition of N-sulfonyl-1-aza-1,3-dienes with vinyl ketones.
2. Catalyst for the enantioselective [3+2] cycloadditions of γ -substituted allenates with β -perfluoroalkyl enones.
3. Organocatalyst used in effective asymmetric vinylogous Mannich reaction of isatin imines with α,α -dicyanoolefins.
4. Chiral phosphine catalyzed α -functionalization of 2-vinylpyridines via enantioselective cross Rauhut-Currier reaction.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)

PHOSPHORUS (Compounds)

15-8514
(continued)N-[(1S)-2-(Diphenylphosphino)-1-phenylethyl]-3,5-bis(trifluoromethyl)-benzamide, 95%
(1853342-57-0)Tech. Note (4)
Ref. (4)

References:

1. *Org. Lett.* **2017**, 19, 1710.
2. *Chem. Sci.*, **2017**, 8, 4660.
3. *Org. Chem. Front.*, **2017**, 4, 101.
4. *Org. Lett.* **2018**, 20, 1304.


15-8366

NEW

[S(R)]-N-[(1S)-2-(Diphenylphosphino)-1-phenylethyl]-2-methyl-2-propanesulfonamide, 95% (1803239-44-2)

C₂₄H₂₈NOPS; FW: 409.5; white to off-white powder.
air sensitive, moisture sensitive




Note: Sold in collaboration with Daicel for research purposes only.



100mg


Technical Notes:

1. Sulfonamide phosphine organocatalyst for the enantioselective intramolecular Rauht-Currier reaction.
2. Organocatalyst used in aza-[4+2] cycloaddition of N-sulfonyl-1-aza-1,3-dienes with vinyl ketones.
3. Used in the enantioselective phosphine-catalyzed allylic alkylations of mix-indene with MBH carbonates.
4. Ligand for the direct Pd-catalyzed cross-coupling reaction of easily accessible secondary phosphine oxides and aryl bromides, which provides rapid access to P-chiral phosphine oxides.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)

PHOSPHORUS (Compounds)

15-8366 [S(R)]-N-[(1S)-2-(Diphenylphosphino)-1-phenylethyl]-2-methyl-2-propanesulfonamide, 95%
(continued) (1803239-44-2)




Tech. Note (4)
Ref. (4)

References:

1. *Angew. Chem. Int. Ed.* **2015**, 54, 6874.
2. *Org. Lett.* **2017**, 19, 1710.
3. *Org. Lett.* **2017**, 19, 6080.
4. *J. Am. Chem. Soc.* **2019**, 141, 20556.

15-8338 [S(R)]-N-[(R)-[2-(Diphenylphosphino)phenyl](4-methoxyphenyl)methyl]-N,2-dimethyl-2-propanesulfonamide, 95%
NEW C₃₁H₃₄NO₂PS; FW: 515.6; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8336 [S(R)]-N-[(S)-[2-(Diphenylphosphino)phenyl](4-methoxyphenyl)methyl]-N,2-dimethyl-2-propanesulfonamide, 95%
NEW C₃₁H₃₄NO₂PS; FW: 515.6; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.


50mg

15-8314 [S(R)]-N-[(S)-[2-(Diphenylphosphino)phenyl](4-methoxyphenyl)methyl]-2-methyl-2-propanesulfonamide, 95% (1616688-62-0)
NEW C₃₀H₃₂NO₂PS; FW: 501.60; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.


100mg

Technical Notes:

1. Ligand for the enantioselective Au-catalyzed [3+3] cycloaddition reaction of 2-(1-alkynyl)-alk-2-en-1-ones with nitrones.
2. Ligand used for Ag(I)-catalyzed enantioselective [3+2] cycloaddition of azomethine ylides.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)

References:

1. *Angew. Chem. Int. Ed.* **2014**, 53, 4350.
2. *Org. Biomol. Chem.*, **2019**, 17, 1395.

PHOSPHORUS (Compounds)


15-8316

NEW

[S(R)]-N-[(R)-[2-(Diphenylphosphino)phenyl] (4-methoxyphenyl)methyl]-2-methyl-2-propanesulfonamide, 95% (1595319-96-2)
 $C_{30}H_{32}NO_2PS$; FW: 501.60; white to off-white powdr.

air sensitive, moisture sensitive


Note: Sold in collaboration with Daicel for research purposes only.




100mg

Technical Notes:


1. Ligand for the enantioselective Au-catalyzed [3+3] cycloaddition reaction of 2-(1-alkynyl)-alk-2-en-1-ones with nitrones.
2. Used in enantioselective Au-catalyzed synthesis of indolyl-substituted cyclopenta[c]furans.
3. Used in the Au-catalyzed enantioselective synthesis of polycyclic indoline skeletons and enantiomerically enriched β -substituted tryptamine-allenes by kinetic resolution.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)

References:

1. *Angew. Chem. Int. Ed.* **2014**, 53, 4350.
2. *Org. Lett.* **2018**, 20, 6403.
3. *Chem. Commun.*, **2019**, 55, 4210.


15-8346

NEW

[S(R)]-N-[(S)-[3,5-Bis(trifluoromethyl)phenyl] [2-(diphenylphosphino)phenyl]methyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{32}H_{30}F_6NOPS$; FW: 621.6; light-yellow to yellow powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.




50mg

PHOSPHORUS (Compounds)

15-8342

NEW

[S(R)]-N-[(R)-[2-(Diphenylphosphino)phenyl]-1-naphthalenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{34}H_{34}NOPS$; FW: 535.7; white to off-white pwdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.




50mg

15-8344

NEW

[S(R)]-N-[(S)-[3,5-Bis(1,1-dimethylethyl)-4-methoxyphenyl][2-(diphenylphosphino)phenyl]methyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{39}H_{50}NO_2PS$; FW: 535.7; white to off-white pwdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8340

NEW

[S(R)]-N-[(S)-[2-(Diphenylphosphino)phenyl]-1-naphthalenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%
 $C_{34}H_{34}NOPS$; FW: 535.7; white to off-white pwdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg

15-8322



NEW

[S(R)]-N-[(S)-[2-(Diphenylphosphino)phenyl]-1-naphthalenylmethyl]-2-methyl-2-propanesulfonamide, 95% (1595319-95-1)
 $C_{33}H_{32}NOPS$; FW: 521.7; white to off-white pwdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

100mg

Technical Notes:


1. Ligand for the enantioselective Au-catalyzed [3+3] cycloaddition reaction of 2-(1-alkynyl)-alk-2-en-1-ones with nitrones.
2. Used in enantioselective Au-catalyzed synthesis of indolyl-substituted cyclopenta[c]furans.
3. Ligand for the enantioselective Au-catalyzed intermolecular exo [4+3]-cycloaddition to generate chiral oxa-bridged benzocycloheptanes.



PHOSPHORUS (Compounds)

15-8322
(continued)

[S(R)]-N-[(S)-[2-(Diphenylphosphino)phenyl]-1-naphthalenylmethyl]-2-methyl-2-propanesulfonamide, 95% (1595319-95-1)

Tech. Note (3)
Ref. (3)

References:

1. *Angew. Chem. Int. Ed.* **2014**, 53, 4350.
2. *Org. Lett.* **2018**, 20, 6403.
3. *Org. Lett.* **2019**, 21, 3018.

15-8334


NEW

[S(R)]-N-[(R)-[2-(Diphenylphosphino)phenyl]phenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%

C₃₀H₃₂NOPS; FW: 485.6; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8332


NEW

[S(R)]-N-[(S)-[2-(Diphenylphosphino)phenyl]phenylmethyl]-N,2-dimethyl-2-propanesulfonamide, 95%

C₃₀H₃₂NOPS; FW: 485.6; light-yellow to yellow powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8312


NEW

[S(R)]-N-[(R)-[2-(Diphenylphosphino)phenyl]phenylmethyl]-2-methyl-2-propanesulfonamide, 95% (1616688-59-5)

C₂₉H₃₀NOPS; FW: 471.60; white to off-white powdr.

air sensitive, moisture sensitive



Note: Sold in collaboration with Daicel for research purposes only.



100mg

Technical Notes:

1. Ligand for the enantioselective Au-catalyzed [3+3] cycloaddition reaction of 2-(1-alkynyl)-alk-2-en-1-ones with nitrones.
2. Ligand used for Ag(I)-catalyzed enantioselective [3+2] cycloaddition of azomethine ylides.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)

References:

1. *Angew. Chem. Int. Ed.* **2014**, 53, 4350.
2. *Org. Biomol. Chem.*, **2019**, 17, 1395.

PHOSPHORUS (Compounds)

- 15-8310** [S(R)]-N-[(S)-[2-(Diphenylphosphino)phenyl]phenylmethyl]-2-methyl-2-propanesulfonamide, 95% (1595319-97-3) 100mg
NEW C₂₈H₃₀NOPS; FW: 471.60; white to off-white powdr
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



Technical Note:

- See 15-8312 (page 73)

- 15-4660** N-[2-(Diphenylphosphinyl)ethyl]-2-[tritylthio]-ethanamine 250mg
NEW C₃₅H₃₄NOPS; FW: 547.70; waxy solid
 Note: Patents: US20170088571A1; U.S. Provisional Patent Application No. 62/827,627.

Technical Note:

- SNP(O) ligand used in Ru-catalyzed selective ester hydrogenation.

Tech. Note (1)
Ref. (1)


References:

- Org. Process Res. Dev. 2020, 24, 415.

- 15-4640** N-[2-(Diphenylphosphinyl)ethyl]-2-[(4-methylphenyl)thio]-ethanamine 250mg
NEW C₂₃H₂₈NOPS; FW: 395.50; white solid
 Note: Patents: US20170088571A1; U.S. Provisional Patent Application No. 62/827,627.

Technical Note:


- See 15-4660 (page 74)



- 15-4620** N-[2-(Diphenylphosphinyl)ethyl]-2-(methylthio)-ethanamine (1839552-30-5) 250mg
NEW C₁₇H₂₂NOPS; FW: 319.4; white solid
 Note: Patents: US20170088571A1; U.S. Provisional Patent Application No 62/827,627.

Technical Note:


- See 15-4660 (page 74)



- 15-4630** N-[2-(Diphenylphosphinyl)ethyl]-2-(phenylthio)-ethanamine 100mg
NEW C₂₂H₂₄NOPS; FW: 381.47; white solid
 Note: Patents: US20170088571A1; U.S. Provisional Patent Application No. 62/827,627.

Technical Note:

- See 15-4660 (page 74)




PHOSPHORUS (Compounds)

15-8364

NEW

[S(R)]-N-[(R)-(3,5-Di-tert-butyl-4-methoxyphenyl)[2-(diphenylphosphino)-4,5-dimethoxyphenyl]-2-methyl-2-propanesulfinamide, 95%
 $C_{40}H_{52}NO_4PS$; FW: 673.9; white to off-white powd.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg




15-8412

NEW

[S(R)]-N-[(R)-[3,5-Di-tert-butyl-4-methoxyphenyl][5-(diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl]methyl]-2-methyl-2-propanesulfinamide, 95%
 $C_{47}H_{56}NO_3PS$; FW: 746.0; white to off-white powd.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

50mg



15-8414

NEW

[S(R)]-N-[(S)-[3,5-Di-tert-butyl-4-methoxyphenyl][5-(diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl]methyl]-2-methyl-2-propanesulfinamide, 95%
 $C_{47}H_{56}NO_3PS$; FW: 746.0; white to off-white powd.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.


50mg

15-8496

NEW

[S(R)]-N-[(S)-(3,5-Di-tert-butyl-4-methoxyphenyl)[2-methyl-2-propanesulfinamide], 95%
 $C_{38}H_{50}NO_2PS$; FW: 573.9; white to off-white powd.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

25mg




15-8468

NEW

[S(R)]-N-[(S)-[3,5-Di-tert-butyl-4-methoxyphenyl][3-benzyloxy-2-(dicyclohexylphosphino)phenyl]methyl]-N,2-dimethyl-2-propanesulfinamide, 95%
 $C_{48}H_{68}NO_3PS$; FW: 746.1; white to off-white powd.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.

25mg



PHOSPHORUS (Compounds)

15-8466


NEW

[S(R)]-N-[(S)-[3,5-Di-tert-butyl-4-methoxyphenyl]][(3-benzyloxy-2-(dicyclohexylphosphino)phenyl)methyl]-2-methyl-2-propanesulfinamide, 95%

C₄₅H₆₆NO₃PS; FW: 732.1; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8504


NEW

[S(R)]-N-[(1S)-1-[2-(Di-tert-butylphosphanyl)phenyl]-2,2-dimethylpropyl]-2-methyl-2-propanesulfinamide, 95%

C₂₃H₄₂NOPS; FW: 411.6; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8502


NEW

[S(R)]-N-[(1R)-1-[2-(Di-tert-butylphosphanyl)phenyl]phenylmethyl]-2-methyl-2-propanesulfinamide, 95%

C₂₅H₃₈NOPS; FW: 431.6; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8500


NEW

[S(R)]-N-[(1S)-1-[2-(Di-tert-butylphosphanyl)phenyl]phenylmethyl]-2-methyl-2-propanesulfinamide, 95%

C₂₅H₃₈NOPS; FW: 431.6; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8492


NEW

[S(R)]-N-[(S)-[2-(Di-tert-butylphosphino)phenyl]adamantylmethyl]-2-methyl-2-propanesulfinamide, 95%

C₂₉H₄₈NOPS; FW: 489.7; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



25mg

15-8494


NEW

[S(R)]-N-[(S)-[2-(Di-tert-butylphosphino)phenyl]([2-naphthalenyl)methyl]-2-methyl-2-propanesulfinamide, 95%

C₂₉H₄₀NOPS; FW: 481.7; white to off-white powdr.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.




25mg

PHOSPHORUS (Compounds)


15-1840 **n-Dodecylphosphonic acid, min. 90% (5137-70-2)** 1g
NEW $\text{CH}_3(\text{CH}_2)_{11}\text{P}(\text{O})(\text{OH})_2$; FW: 250.31; white waxy solid; d. 0.4 5g
 25g

15-8510 **N-[(1S)-(1-Ferrocenyl)-2-(diphenylphosphino)ethyl]-3,5-bis(trifluoromethyl)-benzamide, 95% (2131816-17-4)** 50mg
NEW $\text{C}_{33}\text{H}_{26}\text{F}_6\text{FeNOP}$; FW: 653.4; yellow powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



Technical Note:

1. Ferrocene derived bifunctional phosphine catalyst for asymmetric oxa-[4+2] cycloaddition of α -substituted allenones with enones.




Tech. Note (1)
 Ref. (1)

$\text{R}^1 = \text{CF}_3$, Ar; $\text{R}^2 = \text{CN}$, H, Alkyl; $\text{R}^3 = \text{R}_f$, Ar; $\text{R}^4 = \text{Alkyl}$, $\text{R}^5 = \text{Ar}$, alkyl


References:

1. *Chem. Eur. J.* **2017**, *23*, 13587.


15-8386 **[S(R)]-N-[(1S)-1-Ferrocenylethyl-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95%** 100mg
NEW $\text{C}_{28}\text{H}_{32}\text{FeNOPS}$; FW: 517.4; white to off-white powdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



15-8540 **(11bR)-4-Hydroxy-2,6-bis(adamantan-1-yl)-4-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepinoxide, 95%** 25mg
NEW $\text{C}_{40}\text{H}_{41}\text{O}_4\text{P}$; off white to light yellow powdr.
 Note: Sold in collaboration with Daicel for research purposes only.



15-8544 **(11bR)-4-Hydroxy-2,6-bis(2,6-dimethoxyphenyl)-4-oxide-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin, 98%** 25mg
NEW $\text{C}_{36}\text{H}_{29}\text{O}_8\text{P}$; FW: 620.6; off white to light yellow powdr.
 Note: Sold in collaboration with Daicel for research purposes only.




PHOSPHORUS (Compounds)

15-8534

NEW

(11bR)-4-Hydroxy-2,6-bis(5'-phenyl[1,1':3',1''-terphenyl]-2'-yl)-4-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin oxide, 98%, (99% ee)
 $C_{68}H_{45}O_4P$; FW: 957.1; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

15-8535

NEW

(11bS)-4-Hydroxy-2,6-bis(5'-phenyl[1,1':3',1''-terphenyl]-2'-yl)-4-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin oxide, 98%, (99% ee) (1639447-03-2)


$C_{68}H_{45}O_4P$; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

Technical Note:

- Used in the Lewis base/chiral Brønsted acid catalyzed enantioselective bromocycloetherification.

Tech. Note (1)
Ref. (1)

References:

- Chirality*, **2014**, 26, 344.


15-8530

NEW

(4R)-4-Hydroxy-2,6-bis(3,3'',5,5''-tetrakis(tert-butyl)-[1,1':3',1''-terphenyl]-5'-yl)dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepine-4-oxide, 98%

$C_{88}H_{101}O_4P$; FW: 1253.7; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg


15-8548

NEW

(11bR)-4-Hydroxy-2,6-bis(2,4,6-trimethoxyphenyl)-4-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin, 98%

$C_{38}H_{33}O_{10}P$; FW: 680.6; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

PHOSPHORUS (Compounds)


15-8550

NEW

(11bR)-4-Hydroxy-9,14-dinitro-2,6-bis[2,4,6-tris(1-methylethyl)phenyl]-4-oxide-Dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin, 98% (1535206-22-4)

C₅₀H₅₅N₂O₈P; FW: 848; off white to light yellow powd.


Note: Sold in collaboration with Daicel for research purposes only.




10mg
25mg

Technical Notes:


1. Organocatalyst used for the kinetic resolution of secondary alcohols.
2. Catalyst used in the kinetic resolution of β -amino alcohols via asymmetric substitutions at sp³-hybridized carbon atoms.
3. Brønsted acid catalyst used for asymmetric desymmetrization of α -symmetric acid anhydrides.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)



Tech. Note (3)
Ref. (3)

References:

1. *Angew. Chem. Int. Ed.* **2013**, 52, 10227.
2. *Angew. Chem. Int. Ed.* **2016**, 55, 13137.
3. *Tetrahedron Lett.*, **2016**, 57, 4098.


15-8558

NEW

(6aR)-14-Hydroxy-6,7-diphenyl-14-oxide-dinaphtho[1,2-d:2',1'-f][1,3,2]dioxaphosphepin, 98%, (99% ee) (956610-76-7)

C₃₂H₂₁O₄P; FW: 500.5; off white to light yellow powd.


Note: Sold in collaboration with Daicel for research purposes only.



25mg
100mg

Technical Notes:


1. Organocatalyst for desymmetrization of meso-aziridines.
2. Used in the enantioselective Pd/Brønsted acid-catalyzed carbonylative carbocyclization of enallenes via cross-dehydrogenative coupling with terminal alkynes.



Tech. Note (1)
Ref. (1)

PHOSPHORUS (Compounds)

15-8558 (6aR)-14-Hydroxy-6,7-diphenyl-14-oxide-dinaphtho[1,2-d:2',1'-f][1,3,2]dioxaphosphepin, 98%, (99% ee) (956610-76-7)




Tech. Note (2)
Ref. (2)

References:

1. *J. Am. Chem. Soc.* **2007**, 129, 12084.
2. *Angew. Chem. Int. Ed.* **2017**, 56, 4535.


15-8559 (6aS)-14-Hydroxy-6,7-diphenyl-14-oxide-dinaphtho[1,2-d:2',1'-f][1,3,2]dioxaphosphepin, 98%, (99% ee) (175223-61-7)
C₃₂H₂₁O₄P; FW: 500.5; off white to light yellow powdr.
Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg



Technical Note:

1. Used in the enantioselective racemic resolution of a range of 1,2-amino via liquid-liquid extraction of alcohols.




Tech. Note (1)
Ref. (1)

References:

1. *ChemSusChem* **2018**, 11, 178.

15-8554 (8aR)-18-Hydroxy-8,9-diphenyl-18-oxide-diphenanthro[4,3-d:3',4'-f][1,3,2]dioxaphosphepin, 98%, (99% ee) (871130-18-6)
C₄₀H₂₅O₄P; FW: 600.6; off white to light yellow powdr.
Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg




Technical Notes:


1. Organocatalyst for amidation of imines.
2. Organocatalyst for desymmetrization of meso-aziridines.
3. Chiral Brønsted acid used in catalytic asymmetric aza-Darzens aziridine synthesis.
4. Organocatalyst for the remote α-C–H bond activation of alcohols or amines triggered by the radical trifluoromethylation of alkenes.
5. Used in the enantioselective Pd/Brønsted acid-catalyzed carbonylative carbocyclization of enallenes via cross-dehydrogenative coupling with terminal alkynes.
6. Used in the enantioselective synthesis of tropanes via Brønsted acid catalyzed pseudotransannular desymmetrization.
7. Cu-catalyst used in Cu/achiral pyridine ligand catalyzed enantioselective radical oxytrifluoromethylation of alkenes with alcohols.

PHOSPHORUS (Compounds)


15-8554 (8aR)-18-Hydroxy-8,9-diphenyl-18-oxide-diphenanthro[4,3-d:3',4'-f][1,3,2]dioxaphosphepin, 98%, (99% ee) (871130-18-6)




Tech. Note (1)
Ref. (1)




Tech. Note (2)
Ref. (2)




Tech. Note (3)
Ref. (3)




Tech. Note (4)
Ref. (4)



Tech. Note (5)
Ref. (5)



Tech. Note (6)
Ref. (6)



Tech. Note (7)
Ref. (7)

References:

1. *J. Am. Chem. Soc.* **2005**, 127, 15696.
2. *J. Am. Chem. Soc.* **2007**, 129, 12084.
3. *Org. Lett.* **2011**, 13, 2188.
4. *Org. Chem. Front.*, **2017**, 4, 2139.
5. *Angew. Chem. Int. Ed.* **2017**, 56, 4535.
6. *Angew. Chem. Int. Ed.* **2017**, 56, 8883.
7. *Angew. Chem. Int. Ed.* **2020**, 59, 6780.

PHOSPHORUS (Compounds)


15-8555

NEW

(8aS)-18-Hydroxy-8,9-diphenyl-18-oxide-diphenanthro[4,3-d:3',4'-f][1,3,2]dioxaphosphepin, 98%, (99% ee)
(871130-17-5)

C₄₀H₂₅O₄P; FW: 600.6; light yellow to yellow powdr


Note: Sold in collaboration with Daicel for research purposes only.




25mg
100mg

Technical Notes:


1. Organocatalyst for amidation of imines.
2. Organocatalyst for asymmetric reduction of α -imino esters to generate α -amino esters.
3. Organocatalyst for desymmetrization of meso-aziridines.
4. Organocatalyst for the highly enantioselective addition of indoles to N-acyl imines.
5. Organocatalyst for the enantioselective Friedel-Crafts reaction of pyrrole derivatives with imines.
6. Organocatalyst for the desymmetrization of meso-aziridines with functionalized mercaptans.
7. Organocatalyst for highly enantioselective benzyloxylation of 3-aryloxindoles.
8. Used in for the Ca/Chiral phosphoric acid catalyzed asymmetric chlorination and Michael Reactions of 3-substituted oxindoles.




Tech. Note (1)
Ref. (1-2)




Tech. Note (2)
Ref. (3)




Tech. Note (3)
Ref. (4)



Tech. Note (4)
Ref. (5)





Tech. Note (5)
Ref. (6)



Tech. Note (6)
Ref. (7)

PHOSPHORUS (Compounds)

15-8555 (8aS)-18-Hydroxy-8,9-diphenyl-18-oxide-diphenanthro[4,3-d:3',4'-f][1,3,2]dioxaphosphin, 98%, (99% ee) (871130-17-5)




References:

1. *J. Am. Chem. Soc.* **2005**, 127, 15696.
2. *Chem. Commun.*, **2007**, 4477.
3. *J. Am. Chem. Soc.* **2007**, 129, 5830.
4. *J. Am. Chem. Soc.* **2007**, 129, 12084.
5. *Org. Lett.* **2007**, 9, 2609.
6. *Org. Lett.* **2007**, 9, 4065.
7. *Org. Lett.* **2009**, 11, 5186.
8. *Angew. Chem. Int. Ed.* **2011**, 50, 1135.
9. *J. Am. Chem. Soc.* **2011**, 133, 3339.

15-7560 Isobutylphosphine, min. 95%, CYTOP® 141

NEW

(4023-52-3)
C₄H₉PH₂; FW: 90.11; Colorless liquid; b.p. 78 °C; d. 0.77
air sensitive, moisture sensitive, pyrophoric
Note: Sold in collaboration with Solvay.




25g
100g

15-8498 [S(R)]-N-[(S)-(4-Methoxyphenyl)[2-(di-tert-butylphosphino)phenyl]methyl]-2-methyl-2-propanesulfonamide, 95%

NEW

C₂₆H₄₀NO₂PS; FW: 461.6; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




25mg

15-8416 [S(R)]-N-[(S)-(1-naphthalenyl)[5-(diphenylphosphino)-9,9-dimethyl-9H-xanthen-4-yl]methyl]-2-methyl-2-propanesulfonamide, 95%

NEW

C₄₂H₄₀NO₂PS; FW: 653.8; white to off-white powdr.
air sensitive, moisture sensitive
Note: Sold in collaboration with Daicel for research purposes only.




50mg

PHOSPHORUS (Compounds)

15-0558 (11bS)-8,9,10,11,12,13,14,15-Octahydro-4-hydroxy-2,6-bis[3,5-bis(tri-fluoromethyl)phenyl]phenyl-4-oxide-dinaphtho[2,1-d:1',2'-f][1,3,2]dioxaphosphepin, 98%, (99% ee)
NEW $C_{64}H_{39}F_{24}O_4P$; FW: 1358.9; off white to light yellow pwdr.
 Note: Sold in collaboration with Daicel for research purposes only.



25mg
100mg

15-1279 (1S,1'S)-(-)-(Oxybis(2,1-phenylen))bis((2-methoxyphenyl)(phenyl)phosphine), min. 97% (2119686-55-2)
NEW $C_{38}H_{32}O_3P_2$; FW: 598.61; white to off-white solid
air sensitive, moisture sensitive, (store cold)
 Note: Sold under license of WO2017/191310 for research purposes only.

25mg
100mg

Technical Notes:


1. P-Chirogenic Xantphos ligand used for the synthesis and application in Rhodium-catalyzed asymmetric hydrogenation.
2. Diastereo- and enantioselective copper catalyzed hydroallylation of disubstituted cyclopropenes.




References:

1. *ACS Catal.*, **2017**, 7, 6162–6169.
2. *Chem. Sci.*, **2018**, 9, 6503–6508.

15-7535 9-Phosphabicyclononanes in toluene, mixture of isomers, CYTOP® 282T (13887-02-0)
NEW $C_8H_{15}P$; FW: 142.17; Colorless to pale yellow liq.;
 f.p. 41 °F; d. 0.94
air sensitive, moisture sensitive
 Note: Sold in collaboration with Solvay.


100g
500g

15-8362 [S(R)]-N-[(S)-[2-(Diphenylphosphino)-4,5-dimethoxyphenyl]-1-naphthalenylmethyl]-2-methyl-2-propanesulfonamide, 95%
NEW $C_{35}H_{36}NO_3PS$; FW: 581.7; white to off-white pwdr.
air sensitive, moisture sensitive
 Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-3035 (3-(Tert-butoxy)-2',6'-diisopropyl-6-methoxy-[1,1'-biphenyl]-2-yl)dicyclohexylphosphane
NEW GPhos, 98% (2489243-29-8)
 $C_{35}H_{53}O_2P$; FW: 536.78; white xtl.; m.p. 240
 Note: Patents: US 6,395,916, US 6,307,087


50mg
250mg

PHOSPHORUS (Compounds)

15-3035 (continued) **(3-(Tert-butoxy)-2',6'-diisopropyl-6-methoxy-[1,1'-biphenyl]-2-yl)dicyclohexylphosphane**
GPhos, 98% (2489243-29-8)

Technical Note:

1. Novel mono-phosphine Buchwald ligand that supports a catalyst of improved stability. The GPhos-derived catalyst exhibits improved performance for room temperature C–N cross coupling reactions between a broad scope of primary amines and aryl halides.¹ The ligand is designed to promote C–N cross coupling across a wide range of substrate classes with low catalyst loadings. It shows an improved relative rate of productive catalysis versus catalyst deactivation compared to previous catalysts. The GPhos-supported catalyst exhibits better reactivity than previous catalysts both under ambient conditions and at elevated temperatures. Its use allows for the coupling of a variety of primary amine and aniline nucleophiles, including (1) unhindered, (2) five-membered-ring N-heterocycle-containing, and (3) α -tertiary primary amines, each of which previously required a different catalyst to achieve optimal results. The GPhos-supported catalyst is also well suited for coupling electron-deficient anilines as well as N-heterocycle-containing aryl halides and amines.




References:

1. *J. Am. Chem. Soc.*, **2020**, 142, 35, 15027.

15-6910**NEW**


2-((2R,3R)-3-(tert-Butyl)-4-(3,5-dicyclopentyl-2,6-dimethoxyphenyl)-2,3-dihydrobenzo[d][1,3]oxaphosphol-2-yl)propan-2-ol, min. 97%, (99% ee) (2R,3R)-Baryphos (2416226-97-4)
 C₃₂H₄₅O₃P; FW: 524.68; white solid
 Note: Sold under license from Zejun for research purposes only. Patents: ZL2013105048267, CN104558038.



25mg
100mg
500mg

Technical Note:

1. Ligand/palladium catalyst for asymmetric Suzuki-Miyaura cross-coupling reactions.




References:

1. *J. Am. Chem. Soc.* **2020**, 142, 8036–8043

15-6911**NEW**

2-((2S,3S)-3-(tert-Butyl)-4-(3,5-dicyclopentyl-2,6-dimethoxyphenyl)-2,3-dihydrobenzo[d][1,3]oxaphosphol-2-yl)propan-2-ol, min. 97%, (99% ee) (2S,3S)-Baryphos (2416226-68-9)
 C₃₂H₄₅O₃P; FW: 524.68; white solid
 Note: Sold under license from Zejun for research purposes only. Patents: ZL2013105048267, CN104558038.



25mg
100mg
500mg

Technical Note:

1. See 15-6910 (page 85)

PHOSPHORUS (Compounds)

15-8390 NEW	<p>[S(R)]-N-[(1R)-1-[2-((tert-Butyldiphenylsilyl)oxy)phenyl]]-2-yl-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95% C₄₀H₄₆NO₂PSSi; FW: 663.9; white to off-white powdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>	 <chem>CC(C)(S(=O)(=O)Nc1ccccc1)C[C@H](c2ccccc2)P(c3ccccc3)c4ccccc4</chem>	50mg
15-8388 NEW	<p>[S(R)]-N-[(1S)-1-[2-((tert-Butyldiphenylsilyl)oxy)phenyl]]-2-yl-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95% C₄₀H₄₆NO₂PSSi; FW: 663.9; white to off-white powdr. <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Daicel for research purposes only.</p>		50mg
15-7620 NEW	<p>Tetrabutylphosphonium chloride (68-72 wt% solution in methanol) (2304-30-5) C₁₆H₃₆ClP; FW: 294.88; Pale yellow liq.; f.p. 52 °F; d. 0.9</p>	 <chem>CCCC[P+](CCCC)(CCCC)CCCC.[Cl-]</chem>	100g 500g
15-7625 NEW	<p>Tetrabutylphosphonium chloride (49-51 wt% solution in toluene) (2304-30-5) C₁₆H₃₆ClP; FW: 294.88; Pale yellow liq.; f.p. 43 °F</p>		100g 500g
15-1579 NEW	<p>Tetrabutylphosphonium bromide (75% in water) (3115-68-2) (C₄H₉)₄P+Br-; FW: 339.35; colorless to light brown liq.; f.p. >2200 °F; d. 1.06</p>	 <chem>CCCC[P+](CCCC)(CCCC)CCCC.[Br-]</chem>	100g 500g
15-8596 NEW	<p>(11aR)-10,11,12,13-Tetrahydro-5-hydroxy-3,7-bis[2,4,6-tricyclohexylphenyl]-5-oxide-diindenol[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98% C₆₅H₈₃O₄P; FW: 959.3; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.</p>	 <chem>O=P(O)(c1ccc(C2=CC=CC=C2)cc1)c3c4c5ccc(C6=CC=CC=C6)cc5O[C@H]3C[C@H]4C5</chem>	10mg 25mg
15-8612 NEW	<p>(3aR,8aR)-Tetrahydro-6-hydroxy-2,2-dimethyl-4,4,8,8-tetrakis(3,5-diisopropyl-4-methoxyphenyl)-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98% C₅₉H₈₅O₁₀P; FW: 985.3; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.</p>	 <chem>COc1cc(C=C(C)C)cc(C2=CC=CC=C2)c1C3C4C5C(=C(C)C)C(=C(C)C)C5O[C@H]3C[C@H]4C5</chem>	25mg 100mg
15-8613 NEW	<p>(3aS,8aS)-Tetrahydro-6-hydroxy-2,2-dimethyl-4,4,8,8-tetrakis(3,5-diisopropyl-4-methoxyphenyl)-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98% C₅₉H₈₅O₁₀P; FW: 985.3; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.</p>		25mg 100mg

PHOSPHORUS (Compounds)


15-8617

NEW

(3aR,8aR)-Tetrahydro-6-hydroxy-2,2-dimethyl-4,4,8,8-tetrakis(3,5-di-tert-butylphenyl)-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee)

$C_{63}H_{93}O_6P$; FW: 977.4; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

15-8618

NEW

(3aR,8aR)-Tetrahydro-6-hydroxy-2,2-dimethyl-4,4,8,8-tetrakis(3,5-di-tert-butylphenyl)-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee)

$C_{63}H_{93}O_6P$; FW: 977.4; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg


15-8608

NEW

(3aR,8aR)-Tetrahydro-6-hydroxy-2,2-dimethyl-4,4,8,8-tetrakis(4-(tert-butyl)phenyl)-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee)

$C_{47}H_{61}O_6P$; FW: 753; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

15-8609

NEW

(3aS,8aS)-Tetrahydro-6-hydroxy-2,2-dimethyl-4,4,8,8-tetrakis(4-(tert-butyl)phenyl)-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee)

$C_{47}H_{61}O_6P$; FW: 753; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg


15-8602

NEW

(3aR,8aR)-Tetrahydro-6-hydroxy-2,2-dimethyl-4,4,8,8-tetra-2-naphthalenyl-6-oxide-1,3-dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee) (952649-50-2)


$C_{47}H_{37}O_6P$; FW: 728.8; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.

25mg
100mg

Technical Note:

1. Ligand for the Zn-catalyzed asymmetric Simmons-Smith cyclopropanation of both functionalized and unfunctionalized olefins.




Tech. Note (1)
Ref. (1)

References:




1. *Adv. Synth. Catal.* **2006**, *348*, 2363.

PHOSPHORUS (Compounds)

15-8603 NEW	(3aS,8aS)-Tetrahydro-6-hydroxy-2,2-dimethyl-4,4,8,8-tetra-2-naphthalenyl-6-oxide-1,3-dioxolo[4,5-e][1,3,2]dioxaphosphepin, 98%, (99% ee) C ₄₇ H ₃₇ O ₆ P; FW: 728.8; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
15-8566 NEW	(11aR)-10,11,12,13-Tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-3,7-di-2-naphthalenyl-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee) C ₄₁ H ₃₅ O ₄ P; FW: 622.7; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
15-8567 NEW	(11aS)-10,11,12,13-Tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-3,7-di-2-naphthalenyl-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee) C ₄₁ H ₃₅ O ₄ P; FW: 622.7; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
15-8572 NEW	(11aR)-10,11,12,13-Tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-3,7-di-9-phenanthrenyl-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee) C ₄₉ H ₃₉ O ₄ P; FW: 722.8; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
15-8573 NEW	(11aS)-10,11,12,13-Tetrahydro-10,10,13,13-tetramethyl-5-hydroxy-3,7-di-9-phenanthrenyl-5-oxide-diindeno[7,1-de:1',7'-fg][1,3,2]dioxaphosphocin, 98%, (99% ee) C ₄₉ H ₃₉ O ₄ P; FW: 722.8; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
15-8629 NEW	(3aR,8aR)-4,4,8,8-Tetrakis(3,5-di-tert-butylphenyl)-2,2-bis(4-chlorophenyl)-6-hydroxy-tetrahydro-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 95%, (99% ee) C ₇₃ H ₉₅ Cl ₂ O ₆ P; FW: 1170.4; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg
15-8630 NEW	(3aS,8aS)-4,4,8,8-Tetrakis(3,5-di-tert-butylphenyl)-2,2-bis(4-chlorophenyl)-6-hydroxy-tetrahydro-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 95%, (99% ee) C ₇₃ H ₉₅ Cl ₂ O ₆ P; FW: 1170.4; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.	25mg 100mg



PHOSPHORUS (Compounds)

15-8625	(3aR,8aR)-4,4,8,8-Tetrakis(3,5-di-tert-butylphenyl)-2,2-bis(4-fluorophenyl)-6-hydroxy-tetrahydro-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee) $C_{73}H_{95}F_2O_6P$; FW: 1137.5; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.		25mg 100mg
15-8626	(3aS,8aS)-4,4,8,8-Tetrakis(3,5-di-tert-butylphenyl)-2,2-bis(4-fluorophenyl)-6-hydroxy-tetrahydro-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98% $C_{73}H_{95}F_2O_6P$; FW: 1137.5; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.		25mg 100mg
15-8634	(3aR,8aR)-4,4,8,8-Tetrakis(3,5-di-tert-butylphenyl)-6-hydroxy-2,2-bis(4-(trifluoromethyl)phenyl)tetrahydro-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98% $C_{75}H_{95}F_6O_6P$; FW: 1237.5; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.		25mg 100mg
15-8635	(3aS,8aS)-4,4,8,8-Tetrakis(3,5-di-tert-butylphenyl)-6-hydroxy-2,2-bis(4-(trifluoromethyl)phenyl)tetrahydro-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98% $C_{75}H_{95}F_6O_6P$; FW: 1237.5; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.		25mg 100mg
15-8638	(3aR,8aR)-4,4,8,8-Tetrakis(3,5-di-tert-butylphenyl)-6-hydroxy-2,2-di-p-tolyltetrahydro-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee) $C_{75}H_{101}O_6P$; FW: 1129.6; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.		25mg 100mg
15-8639	(3aS,8aS)-4,4,8,8-Tetrakis(3,5-di-tert-butylphenyl)-6-hydroxy-2,2-di-p-tolyltetrahydro-6-oxide-[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee) $C_{75}H_{101}O_6P$; FW: 1129.6; off white to light yellow powdr. Note: Sold in collaboration with Daicel for research purposes only.		25mg 100mg

PHOSPHORUS (Compounds)


15-8621

NEW

(3a'R,8a'R)-4',4',8',8'-Tetrakis(3,5-di-tert-butylphenyl)-6'-hydroxy-tetrahydrospiro[cyclohexane-1,2'-6'-ox[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee)

C₆₆H₉₇O₆P; FW: 1017.4; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.



25mg
100mg

15-8622

NEW

(3a'S,8a'S)-4',4',8',8'-Tetrakis(3,5-di-tert-butylphenyl)-6'-hydroxy-tetrahydrospiro[cyclohexane-1,2'-6'-ox[1,3]dioxolo[4,5-e][1,3,2]dioxaphosphepine, 98%, (99% ee)

C₆₆H₉₇O₆P; FW: 1017.4; off white to light yellow powdr.

Note: Sold in collaboration with Daicel for research purposes only.


25mg
100mg

15-7630

NEW

Tetraoctylphosphonium bromide (45-55 wt% in toluene) (23906-97-0)

C₃₂H₆₈BrP; FW: 563.76; Colorless to pale yellow liq.; f.p. 43 °F
hygroscopic



100g
500g

Technical Notes:

1. Thermally stable, soluble quaternary phosphonium derivative used in metal separation¹ and lubricant additive applications².
2. Catalyst for the oxidation of alkyl aromatic compounds³ and preparation of mononitrated aromatic compounds⁴.

References:

1. *Green Chemistry*, **2012**, 14, 1657-1665.
2. US20150232777A1.
3. US9932287B2.
4. WO2012156540A3.


15-7580

NEW

Tributyl(ethyl)phosphonium bromide, min. 95% (45-55% in water), CYPHOS® IL 676W (7392-50-9)

C₁₄H₃₂BrP; FW: 311.29; Colorless to pale yellow liq.


Note: Sold in collaboration with Solvay.




100g
500g

Technical Notes:

1. CYPHOS® IL 676W has been used in phase transfer catalysis [1].
2. CYPHOS® IL 676W is an effective in aluminum electroplating [2].
3. CYPHOS® IL 676W has also been exploited in organic transformations as a cocatalyst [3].
4. CYPHOS® IL 676W has been observed to increase the SN2 reactions rate of sulfonate esters with bromide ion [5].
5. CYPHOS® IL 676W has been demonstrated effective at capturing SO2 from flue gas [5].



Tech. Note (3)
Ref. (3)









Tech. Note (4)
Ref. (4)

PHOSPHORUS (Compounds)

15-7580 Tributyl(ethyl)phosphonium bromide, min. 95% (45-55% in water), CYPHOS® IL 676W
(continued) (7392-50-9)

References:

1. JP Patent 02304061, Dec. 17, **1990**.
2. JP Patent 04021794, Jan. 24, **1992**.
3. *Tetrahedron Lett.* **2004**, 45, 7359.
4. *J. Phys. Org. Chem.* **2006**, 19, 281.
5. *AIChE*, **2015**, 6, 2028.

15-7570 NEW	Tributylmethylphosphonium p-toluenesulfonate, min. 95% (45-55% in water), CYPHOS® IL 120W (55767-12-9) $C_{13}H_{30}P(C_7H_7O_3S)$; FW: 388.55; Colorless to pale yellow liq. Note: Sold in collaboration with Solvay.		100g 500g
15-7575 NEW	Tributyl(octyl)phosphonium chloride, min. 93%, CYPHOS® IL 253 (56315-19-6) $C_{26}H_{44}ClP$; FW: 351.00; Colorless to pale yellow viscous liq.; b.p. >2200 °C; f.p. 306 °F; d. 0.92 Note: Sold in collaboration with Solvay.		100g 500g
15-7615 NEW	Tributylphosphine sulfide/tributylphosphine (50:50 mixture), CYTOP® 506 (3084-50-2) $C_{12}H_{27}PS/C_{12}H_{27}P$; Colorless to amber liq.; f.p. >2212 °F <i>air sensitive, moisture sensitive</i> Note: Sold in collaboration with Solvay.		25g 100g 500g
15-7585 NEW	Tributyltetradecyl phosphonium chloride, (ca. 50% in water), CYPHOS® 4345W (81741-28-8) $C_{26}H_{56}ClP$; FW: 435.16; Yellow green viscous liq.; d. 0.95 Note: Sold in collaboration with Solvay.		100g 500g
15-7565 NEW	Triethyl(octyl)phosphonium chloride (45-55% solution), CYPHOS® IL 541W (482647-71-2) $C_{14}H_{32}ClP$; FW: 266.83; Colorless to pale yellow liq. Note: Sold in collaboration with Solvay.		100g 500g
15-7590 NEW	Trihexyl(octyl)phosphonium chloride, min. 95%, CYPHOS® IL 356 (850134-85-9) $C_{26}H_{56}ClP$; FW: 435.16; Colorless to pale yellow liq.; f.p. >230 °F Note: Sold in collaboration with Solvay.		100g 500g

PHOSPHORUS (Compounds)

15-8374


NEW

[S(R)]-N-[(1S)-1-(2',4',6'-Triisopropyl)-(1,1'-bi-phenyl)-2-yl-2-(diphenylphosphino)ethyl]-2-methyl-2-propanesulfonamide, 95%

$C_{39}H_{50}NOPS$; FW: 611.9; white to off-white powder.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



50mg

15-8456


NEW

[S(R)]-N-[(R)-[2,4,6-(Triisopropylphenyl)methyl]-4-methoxyphenyl][2-(dicyclohexylphosphino)phenyl][methyl]-2-methyl-2-propanesulfonamide, 95%

$C_{38}H_{60}NOPS$; FW: 609.9; white to off-white powder.

air sensitive, moisture sensitive

Note: Sold in collaboration with Daicel for research purposes only.



50mg


15-6665

NEW

Tri-n-octylphosphine, min. 90% (4731-53-7)

$(C_8H_{17})_3P$; colorless to pale yellow liq.; f.p. >2200 °F; d. 0.83

air sensitive



100g

500g

Technical Notes:

1. Ligand used in the manufacture of nanocrystals, specifically quantum dots.¹
2. A General Phosphorus Source for the Low-Temperature Conversion of Metals into Metal Phosphides.²
3. Raw material for chalcogenides TOP=E (E = O, S, Se), used for II-VI semiconductor nanoparticles such as CdSe and CdTe, and IV-VI nanoparticles, such as PbSe.

References:

1. *Physika E*, **2017**, 90, 175-182.
2. *Chem. Mater.* **2007**, 19, 17, 4234-4242.
3. *Mater. Chem.* **2010**, 20, 5797-5809.


15-6659

NEW

Trioctylphosphine selenide min. 80% (20612-73-1)

$(C_8H_{17})_3PSe$; FW: 449.59; colorless to pale yellow liq.; f.p. >230 °F; d. 0.9

STENCH



25g

100g

500g

15-7550


NEW

Trioctylphosphine sulfide, CYTOP® 505

(2551-53-3)

$C_{24}H_{51}PS$; FW: 402.71; colorless to pale yellow liq.; f.p. 230 °F; d. 0.89

Note: Sold in collaboration with Solvay.



25g

100g

500g

15-7610

NEW

amp


Tripropylphosphine, min. 98%, CYTOP® 330

(2234-97-1)

$C_9H_{21}P$; FW: 160.24; Colorless liq.; d. 0.8

air sensitive, moisture sensitive, pyrophoric



Note: Sold in collaboration with Solvay.




50g

250g



PHOTOCHEMICAL EQUIPMENT (Other)

98-7850	EvoluChem™ Light proofing upgrade Note: This item is compatible with the EvoluChem™ PhotoRedOx Box (98-7500 & 98-7550). Sold in collaboration with HepatoChem.	1 ea
NEW		
98-7975	EvoluChem™ PhotoRedOx Box Light Source - Electric Power PAR20-18W, LG, Wavelength 380nm, 25 degrees, 100-240V AC Note: This item is compatible with the EvoluChem™ PhotoRedOx Box (98-7500 & 98-7550). A light source is required for use of the PhotoRedOx Box. For related products, see www.strem.com/evoluchem . Sold in collaboration with HepatoChem.	1 ea
NEW		
98-7950	EvoluChem™ PhotoRedOx Box Light Source - Wavelength 365nm, Electric Power 18W Note: This item is compatible with the EvoluChem PhotoRedOx Box (98-7500 & 98-7550). A light source is required for use of the PhotoRedOx Box. For related products see www.strem.com/evoluchem . Sold in collaboration with HepatoChem.	1ea
NEW		
98-7780	EvoluChem™ PhotoRedOx Box Photochemistry Holder - 2 x 20 ml vials Note: This item is compatible with the EvoluChem PhotoRedOx Box (98-7500 & 98-7550). A photochemistry holder is required for use of the PhotoRedOx Box. For additional holder options, please see www.strem.com/evoluchem . Sold in collaboration with HepatoChem.	1ea
NEW		
98-7790	EvoluChem™ PhotoRedOx Box Photochemistry Holder - 8 x 4 ml vials Note: This item is compatible with the EvoluChem™ PhotoRedOx Box (98-7500 & 98-7550). A photochemistry holder is required for use of the PhotoRedOx Box. For additional holder options, please see www.strem.com/evoluchem . Sold in collaboration with HepatoChem.	1 ea
NEW		
98-7550	EvoluChem™ PhotoRedOx Box (Temperature controlled) Note: Sold in collaboration with HepatoChem.	1 ea
NEW		

RHENIUM (Compounds)

75-4200	Trimethylsilyl perrhenate, 98% (99.999%-Re) PURATREM (16687-12-0) $\text{ReO}_5\text{OSi}(\text{CH}_3)_3$; FW: 323.39; white solid; m.p. 79.5-80.5 <i>air sensitive, light sensitive, moisture sensitive, hygroscopic</i>	250mg 1g 5g
NEW		

RUTHENIUM (Compounds)

44-3210	Dichloro[N-[2-(diphenylphosphino-κP)ethyl]-2-(methylthio-κS)ethanamine-κN](triphenylphosphine) ruthenium $\text{C}_{35}\text{H}_{37}\text{Cl}_2\text{NP}_2\text{RuS}$; FW: 737.67; yellow solid	50mg
NEW		
Technical Note: 1. Catalyst for selective ester hydrogenation.		
$\text{R}-\text{C}(=\text{O})-\text{O}-\text{R}^1 + \text{H}_2 \xrightarrow[\text{KO}-t\text{C}_4\text{H}_9, 40^\circ\text{C}]{[\text{Ru}] \text{ SNP(O) Ligand (50 bar)}} \text{R}-\text{CH}_2-\text{OH} + \text{R}^1\text{OH}$		
References: 1. <i>Org. Process Res. Dev.</i> 2020 , 24, 415.		
44-3245	Dichloro[rel-[N(S)]-N-[2-(diphenylphosphinyl-κO)ethyl]-2-[(R)-benzylthio-κS]ethanamine-κN](triphenylphosphine) ruthenium $\text{C}_{41}\text{H}_{41}\text{Cl}_2\text{NOP}_2\text{RuS}$; FW: 829.76; magenta powdr. <i>air sensitive</i> Note: Patents: US20170088571A1; U.S. Provisional Patent Application No. 62/827,627.	50mg
NEW		
Technical Note: 1. See 44-3210 (page 93)		

RUTHENIUM (Compounds)

44-3230


NEW

Dichloro[rel-[N(S)]-N-[2-(diphenylphosphinyl-κO)ethyl]-2-[(R)-methylthio-κS]ethanamine-κN](triphenylphosphine) ruthenium (1839552-42-9)

$C_{35}H_{37}Cl_2NOP_2RuS$; FW: 753.66; brick red solid
moisture sensitive

Note: Patents: US20170088571A1; U.S.

Provisional Patent Application No. 62/827,627.



100mg

Technical Note:

1. See 44-3210 (page 93)

44-3240


NEW

Dichloro[rel-[N(S)]-N-[2-(diphenylphosphinyl-κO)ethyl]-2-[(R)-phenylthio-κS]ethanamine-κN](triphenylphosphine) ruthenium

$C_{40}H_{39}Cl_2NOP_2RuS$; FW: 815.74; dark red xtl.
air sensitive

Note: Patents: US20170088571A1; U.S.

Provisional Patent Application No. 62/827,627.



50mg

Technical Note:

1. See 44-3210 (page 93)

44-3235


NEW

Dichloro[rel-[N(S)]-N-[2-(diphenylphosphinyl-κO)ethyl]-2-[(R)-tert-butylthio-κS]ethanamine-κN](triphenylphosphine) ruthenium

$C_{39}H_{43}Cl_2NOP_2RuS$; FW: 795.75; black brown solid
air sensitive

Note: Patents: US20170088571A1; U.S.

Provisional Patent Application No. 62/827,627.



100mg

Technical Note:

1. See 44-3210 (page 93)


STRONTIUM (Compounds)

38-1200

NEW

Bis(tri-isopropylcyclopentadienyl)strontium 98% (147658-82-0)

$[(C_3H_7)_3C_5H_2]_2Sr$; FW: 470.29; yellow liq.
air sensitive, moisture sensitive

500mg
2g

Technical Note:

1. Homoleptic strontium Cp precursor for ALD of Sr containing chalcogenide thin films.

STRONTIUM (Compounds)

38-1200 Bis(tri-isopropylcyclopentadienyl)strontium 98% (147658-82-0)
(continued)

Target Deposit	Deposition Technique	Delivery Temperature	Pressure	Co-reactants	Deposition Temperature	Ref.
SrO	ALD	100 °C	7.5 Torr	H ₂ O	250 °C	[1]
SrO	ALD	140-145 °C	1 Torr	H ₂ O	150-390 °C	[7]
SrTiO ₃	ALD	100 °C	7.5 Torr	H ₂ O, Ti(OPr) ₄	250-325 °C	[1]
SrTiO ₃	ALD	150 °C	0.2 Torr	H ₂ O, Ti(NMe ₂) ₄	120-330 °C	[10]
SrTiO ₃	ALD	90 °C	0.7 Torr	O ₃ , H ₂ O, Ti(Me ₃ Cp)(OMe) ₃	370 °C	[6]
SrTiO ₃	ALD	90 °C	0.7 Torr	O ₃ , H ₂ O, Ti(OPr) ₂ (tmhd) ₂	370 °C	[6]
SrTiO ₃	ALD	130 °C		O ₃ , Ti(NEtMe) ₄	300 °C	[4]
SrTiO ₃	PEALD	120 °C		O ₂ plasma, Ti(OPr) ₂ (tmhd) ₂	250 °C	[5]
La _x Sr _{1-x} TiO ₃	ALD	130 °C	1 Torr	H ₂ O, Ti(OPr) ₄ , La(fmd) ₃	225 °C	[9]
SrSnO ₃	ALD	130 °C	1 Torr	Sn(NMe ₂) ₄ , H ₂ O	180 °C	[12]
SrZrO ₃	ALD	130 °C	1 Torr	H ₂ O, Zr(Nme ₂) ₄	225 °C	[11]
SrRuO ₃	CVD/ALD			O ₂ , RuO ₄	230 °C	[8]
SrS	ALD	100 °C	7.5 Torr	H ₂ S	120-460 °C	[2]
M:SrS (M=Cu, Pb, Ce, Mn, Eu)	ALD	100 °C	7.5 Torr	H ₂ S, Cu(Cp)(PEt) ₃ , Pb(Bu ₃ Cp) ₃ , Ce(Cp) ₃ , Mn(thmd) ₃ , Eu(thmd) ₃	120-380 °C	[3]





Thermal Behavior:

- TGA available in [2]
- Decomposition in process ~250 °C [1, 7]
- Vapor pressure ~0.5 Torr at 120 °C [5]


References:

1. *Chem. Vap. Deposition* **2001**, 7, 75.
2. *Chem. Mater.* **2002**, 14, 1937.
3. *J. Electrochem. Soc.* **2004**, 151, H221.
4. *J. Appl. Phys.* **2011**, 109, 094101.
5. *J. Electrochem. Soc.* **2011**, 158, G34.
6. *Chem. Mater.* **2013**, 25, 953.
7. *J. Phys. Chem. C* **2013**, 117, 11578.
8. *ECS Trans.* **2013**, 58, 171.
9. *J. Appl. Phys.* **2014**, 115, 224108.
10. *Thin Solid Films* **2014**, 550, 53.
11. *J. Appl. Phys.* **2018**, 124, 044102.
12. *J. Vac. Sci. Technol. A*, **2019**, 37, 050902.

BORON (Compounds)

05-0110 HAZ	CALLERY™ Dimethylsulfide borane, min. 93% (13292-87-0) (CH ₃) ₂ S·BH ₃ ; FW: 75.97; colorless to slightly yellow liq.; m.p. -40°; f.p. 65°F; d. 0.801 <i>moisture sensitive, (store cold)</i> Note: Product and trademark of Ascensus Specialties LLC		25g 100g
05-1000 HAZ	CALLERY™ (R)-Methyl oxazaborolidine, 1M in toluene (112022-83-0) C ₁₈ H ₂₀ BNO; FW: 277.17; colorless to amber liq.; f.p. 40°F (toluene); d. 0.925 <i>moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		10ml 50ml
05-1001 HAZ	CALLERY™ (S)-Methyl oxazaborolidine, 1M in toluene (112022-81-8) C ₁₈ H ₂₀ BNO; FW: 277.17; colorless to amber liq.; f.p. 40°F (toluene); d. 0.925 <i>moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		10ml 50ml
93-0540 HAZ 	CALLERY™ Triethylborane, min 99.8% (97-94-9) B(C ₂ H ₅) ₃ ; FW: 98.00; clear to slightly amber liq.; m.p. -92.9°; b.p. 95°; f.p. 32.8°F; d. 0.6961 (23°) <i>pyrophoric</i> Note: Product and trademark of Ascensus Specialties LLC		100g



LITHIUM (Compounds)

03-1275 NEW HAZ	CALLERY™ Lithium hexamethyldisilazane, 24% solution in tetrahydrofuran (4039-32-1) C ₆ H ₁₈ LiNSi ₂ ; FW: 167.33; yellow to amber; f.p. -6.16; d. 0.867 <i>moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		100g 500g 2.5kg
------------------------------	--	---	-----------------------



POTASSIUM (Elemental Forms)

19-1900 amp HAZ	CALLERY™ Potassium Metal, min. 99% (breakseal ampoule) (7440-09-7) K; FW: 39.09; under argon; m.p. 63.7°; b.p. 760°; d. 0.86 <i>air sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		1g 5g
93-1900 amp HAZ	CALLERY™ Potassium Metal, min. 99% (prescored ampoule) (7440-09-7) K; FW: 39.09; under argon; m.p. 63.7°; b.p. 760°; d. 0.86 <i>air sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		1g 5g 25g
19-1910 amp HAZ	CALLERY™ Sodium-Potassium Alloy 22:78 (min. 99%) (11135-81-2) NaK; liq.; 78% K, 22% Na (under argon); m.p. -12.6°; b.p. 785°; d. 0.855 <i>air sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		1g 5g 25g 50g





POTASSIUM (Compounds)

19-1045 NEW HAZ	CALLERY™ Potassium hexamethyldisilazane, 15% in toluene (40949-94-8) [(CH ₃) ₃ Si] ₂ NK; FW: 199.49; colorless to amber liquid; f.p. -6.88; d. 0.879 <i>air sensitive, moisture sensitive, STENCH</i> Note: Product and trademark of Ascensus Specialties LLC		100g 500g 2.5kg
19-1000 NEW HAZ	CALLERY™ Potassium tert-amylate, 25% solution in toluene (41233-93-6) C ₅ H ₁₁ KO; FW: 126.24; clear to amber; f.p. 39.92 <i>air sensitive, heat sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		250g 1kg 5kg














POTASSIUM (Compounds)

19-1020 NEW HAZ	CALLERY™ Potassium tert-butoxide, 20% solution in tetrahydrofuran (865-47-4) C_4H_9KO ; FW: 112.21; colorless to amber; f.p. 0.14; d. 0.906 <i>air sensitive, heat sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		100g 500g 2.5kg
93-1910 HAZ	CALLERY™ Potassium tert-butoxide, min. 99% (865-47-4) $KOC(CH_3)_3$; FW: 112.21; white to off-white powdr.; m.p. 256° dec. <i>moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		100g 500g 2.5kg














SODIUM (Compounds)

11-1275 NEW HAZ	CALLERY™ Sodium hexamethyldisilazane, 40% solution in tetrahydrofuran (1070-89-9) $C_6H_{18}NNaSi_2$; FW: 183.38; yellow to brown; f.p. -23.98 <i>air sensitive, heat sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		100g 500g 2.5kg
11-1630 NEW HAZ	CALLERY™ Sodium tert-amylate, 35% in tetrahydrofuran (14593-46-5) $(C_5H_{11})(CH_3)_2CONa$; FW: 110.13; yellow to amber liquid; f.p. -6.16 <i>air sensitive, heat sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		500g 2.5kg
93-1022 HAZ	CALLERY™ Sodium tert-butoxide, min. 98% (865-48-5) C_4H_9ONa ; FW: 96.11; White to light yellow powdr. <i>moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		100g 500g
11-0595 HAZ	CALSELECT™ Sodium triacetoxymethylborohydride, min. 97% (56553-60-7) $Na(OOCCH_3)_3BH$; FW: 211.94; white powdr.; m.p. 106 °C <i>air sensitive, moisture sensitive, (store cold)</i> Note: Product and trademark of Ascensus Specialties LLC		25g 100g 500g 2.5kg
11-1117 NEW HAZ	VENPURE™ AF CAPLETS (sodium borohydride), min. 99% (16940-66-2) $NaBH_4$; FW: 37.83; White caplets; m.p. > 360 °C <i>air sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		50g 250g 1kg
11-1120 NEW HAZ	VENPURE™ AF GRANULE (sodium borohydride), min. 99% (16940-66-2) $NaBH_4$; FW: 37.83; White granules; m.p. > 360 °C <i>air sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		50g 250g 1kg
11-1119 NEW HAZ	VENPURE™ SF GRANULE (sodium borohydride), min. 98.5% (16940-66-2) $NaBH_4$; FW: 37.83; White granules; m.p. > 360 °C <i>air sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		50g 250g 1kg
93-1118 HAZ	VENPURE™ SF POWDER (sodium borohydride), min. 98.5% (16940-66-2) $NaBH_4$; FW: 37.83; white powdr.; m.p. > 360 °C; d. 1.074 <i>moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		50g 250g 1kg
11-1121 NEW HAZ	VENPURE™ SOLUTION (sodium borohydride 12%; sodium hydroxide 40%) in water (16940-66-2) $NaBH_4$; FW: 37.83; Clear colorless; b.p. 130-135 °C; d. 1.4 <i>air sensitive, moisture sensitive</i> Note: Product and trademark of Ascensus Specialties LLC		500g 2.5kg













Coming Soon

05-1785	CALLERY™ 9-Borabicyclo-[3.3.1]-nonane in tetrahydrofuran 0.5M (rep.monomer) (280-64-8)	
05-0175	CALLERY™ Borane tetrahydrofuran, 1M in tetrahydrofuran (14044-65-6)	
05-0156	CALLERY™ (+)-Diisopinocampheylchloroborane, 60-65% solution in hexanes ((+)-DPC) (112246-73-8)	
05-0155	CALLERY™ (-)-Diisopinocampheylchloroborane, 60-65% solution in heptanes ((-)-DPC) (85116-73-6)	
05-0130	CALLERY™ 5-Ethyl-2-methylpyridine borane, min. 93% (1014979-56-6)	
05-1655	CALLERY™ Isopropyl Pinacol Borate, min. 98% (61676-62-8)	
03-1430	CALLERY™ Lithium tert-Butoxide, 20% solution in tetrahydrofuran (1907-33-1)	
03-1410	CALSELECT™ Lithium triethylborohydride, 1M in tetrahydrofuran (22560-16-3)	
03-1400	CALSELECT™ Lithium tri-sec-butylborohydride, 1M in tetrahydrofuran (38721-52-7)	
03-4870	CALSELECT™ Lithium Tri-tert-butoxy Aluminum Hydride, 1M in tetrahydrofuran (17476-04-9)	
05-1750	CALLERY™ Methoxydiethylborane, 50% in tetrahydrofuran (7397-46-8)	
05-1700	CALLERY™ Methoxydiethylborane, min. 97% (7397-46-8)	
05-1635	CALLERY™ Methyl Boronic Acid, min. 96.5% (13061-96-6)	

Coming Soon


05-1780	CALLERY™ N,N-Diethylaniline borane, min. 97% (13289-97-9)	
05-0145	CALLERY™ Pinacolborane with 1% triethylamine, min. 97.5% (25015-63-8)	
19-1035	CALLERY™ Potassium ethoxide, 24% in ethanol (DN w/ 5% toluene) (917-58-8)	
19-1090	CALLERY™ Potassium hexamethyldisilazane, 20% solution in tetrahydrofuran (40949-94-8)	
19-1040	CALLERY™ Potassium isopropoxide, 19% in isopropanol (6831-82-9)	
19-1095	CALLERY™ Potassium Methoxide, 25% solution in Methanol (865-33-8)	
19-1070	CALLERY™ Potassium tert-amylate, 5% solution in cyclohexane (41233-93-6)	
19-1060	CALLERY™ Potassium tert-amylate, 15% solution in cyclohexane (41233-93-6)	
19-1085	CALLERY™ Potassium tert-amylate, 15% solution in hexanes (41233-93-6)	
19-1080	CALLERY™ Potassium tert-butoxide, 12% solution in tetrahydrofuran (865-47-4)	
19-1030	CALLERY™ Potassium 3,7-Dimethyl-3-octylate, 50% in n-heptane (263148-42-1)	
11-1615	CALLERY™ Sodium mentholate, 50% in hexanes (19321-38-1)	
11-1725	CALLERY™ Sodium n-butoxide, 20% solution in n-Butanol (2372-45-4)	

Coming Soon

11-1625	CALLERY™ Sodium tert-amylate (14593-46-5)	
11-1740	CALLERY™ Sodium tert-amylate, 25% in Cyclohexane (14593-46-5)	
11-1635	CALLERY™ Sodium tert-amylate, 25% in Toluene (14593-46-5)	
11-1735	CALLERY™ Sodium tert-butoxide, 20% solution in tetrahydrofuran (865-48-5)	
05-0565	CALLERY™ Tri-n-butylborane-3-methoxypropylamine, min. 98% (345269-15-0)	
05-0140	CALLERY™ Triethylamine borane, min. 95% (1722-26-5)	
05-0550	CALLERY™ Triethylborane, 14% in tetrahydrofuran (97-94-9)	
05-1665	CALLERY™ Triethylborane (85%); Triethylaluminum (15%) (97-93-8; 97-94-9)	
05-0555	CALLERY™ Triethylborane-1,3-diaminopropane, min. 97% (148861-07-8)	
05-0560	CALLERY™ Triethylborane-diethylenetriamine, min. 97% (1187733-83-0)	
05-1675	CALLERY™ Trimethyl Borate Azeotrope, 72% in Methanol (121-43-7; 67-56-1)	
05-1645	CALLERY™ Trimethylboroxin, 50% solution in tetrahydrofuran (823-96-1)	

Coming Soon

11-1600 CALSELECT™ Sodium tri-sec-butylborohydride, 1M in tetrahydrofuran (67276-04-4)



03-1450 VENPURE™ Lithium borohydride 10%, in tetrahydrofuran (16949-15-8)

11-1122 VENPURE™ 930 Solution (sodium borohydride 9%; sodium hydroxide 30%) (16940-66-2)



cGMP & Custom Synthesis

Providing reliable solutions for over 30 years

Ascensus delivers consistent quality and expertise in cGMP and custom synthesis production. All chemists at the Newburyport site are regularly trained on GMP standards and participate in GMP production to provide customers with reliable products and services.

cGMP Production

Custom synthesis services and process development are conducted under the strictest confidence. A broad range of chemical transformations are possible at our Newburyport facility. For pharmaceutical applications, manufacturing is conducted under current Good Manufacturing Practices (cGMP) in kilo-lab suites. This capability has been available at the Newburyport site since 1988. Complete documentation is provided, including validation and stability studies. Both Active Pharmaceutical Ingredients (APIs) and late stage raw materials for APIs are manufactured. Active Drug Master Files (DMFs) are also maintained in North America and Europe. Quality agreements are in place with major US and International drug companies.

Site Capabilities

GMP Suites & Facilities

- 3 ISO Class 7 rooms
- 2 controlled, non-classified GMP rooms
- Dedicated glassware
- Dedicated production suite
- Verified Line Clearance prior to project change-over

Analytical Support

- Method Development
- Methods Transfer
- Methods Verification/ Validation
- Data Integrity
- Release and Stability Testing

Equipment & Scale

- Kilo-lab production
- Up to 22L dedicated glassware
- Capable of 50-75L
- Expansion capacity (150L)
- 10s to few 100s of kilos of API
- Cryogenic (down to -78°C)
- Max temperature 265°C

Cleaning

- Dedicated glassware
- Cleaning SOPs
 - Rinse release
 - Depyrogenated
 - Swab-released
- Cleaning Validation


















Current APIs

- | | | |
|------------------------|-------------------------------|------------------------------|
| • Cisplatin | • Isoproterenol HCl* | • Chromic Chloride |
| • Carboplatin | • Selenious Acid** | • Pamidronic Acid* |
| • Dexmedetomidine HCl* | • Cupric Sulfate Pentahydrate | • Sodium Sulfide Nonahydrate |

* Marketed exclusively through Chemwerth

** Available through Chemwerth if a DMF is required

Available Booklets



Visit strem.com for new product announcements.

THE STREM CHEMIKER

**Ascensus Specialties
Strem Chemicals, Inc.**

7 Mulliken Way
Newburyport, MA 01950-4098 U.S.A.
Tel.: (978) 499-1600 Fax: (978) 465-3104
(Toll-free numbers below US & Canada only)
Tel.: (800) 647-8736 Fax: (800) 517-8736

CHEMISTRY MADE FOR WHAT MATTERS

APIs
BORANES
BOROHYDRIDES
CATALYSTS & LIGANDS
CVD/ALD PRECURSORS
OTHER SPECIALTY CHEMICALS
SPECIALTY ALCOHOLATES
STREM CATALOG

Follow us and view our blogs on strem.com

