


Kit Catalog Number: 96-0350

Active Catalog number **05-0058** Fluorescein bis(benzyl boronic ester) **FBBBE**
Ingredients: Catalog number **05-0054** Coumarin benzyl boronic ester **CBBE**

Sold under license from The Regents of the University of California.
 US Patent application 61/762,706.



<p>Introduction:</p>	<p>FBBBE and CBBE are robust fluorescent probes, effective for imaging hydrogen peroxide in biological settings. In the supplied form, the fluorescence signal is effectively quenched by the addition of the benzyl ether boronic ester protecting groups to the fluorophores. In the presence of biologically relevant levels of hydrogen peroxide (10-200 :M), the protecting groups cleave and a strong fluorescent signal is observed. (FBBBE: $I_{exc} = 480$ nm, $I_{em} = 512$ nm; CBBE: $I_{exc} = 370$ nm, $I_{em} = 450$ nm) 1. <i>ChemBioChem</i>, 2013, 14, 593.</p>
<p>Contents:</p>	<p>FBBBE: 3x5mg CBBE: 3x5mg Dimethylsulfoxide (ACS spectrophotometric grade): 6 x 1.0mL</p>
<p>SDS:</p>	<p>The Safety Data Sheets for the three products contained in this kit can be downloaded from the Strem Chemicals Web Site at www.strem.com. Locate the SDS using the following catalog numbers: FBBBE: 05-0058 CBBE: 05-0054 Dimethylsulfoxide (ACS spectrophotometric grade): 97-4940</p>
<p>Storage Conditions:</p>	<p>Both compounds are stable at ambient conditions for several days. For long term storage, the kit should be stored at -20°C and protected from light.</p>
<p>Preparation of FBBBE:</p> <p>Preparation of CBBE:</p>	<p>Step 1: Add 654 μL of DMSO to 5.0 mg FBBBE (MW: 764.3 g/mol) (resulting concentration = 10 mM). FBBBE is readily soluble in DMSO. The solution can be portioned into 50 μL aliquots as needed. These solutions must be stored in the freezer at < -20°C. Step 2: Add 50 μL of 10 mM stock to 50 μL DMSO (resulting concentration = 5 mM). This stock solution can be stored in the freezer at < -20°C. Step 3: This solution can be further diluted to the desired concentration in buffer. For example add 10 μL of the 5 mM stock solution to 990 μL 1X PBS, resulting in a 50 μM solution (this solution should be used the same day in which it was prepared). This will result in a solution containing only 1% DMSO (v/v). This solution can then be used in the biological setting to image H₂O₂.</p> <p>Step 1: Add 661 L of DMSO to 5.0 mg CBBE (MW: 378.2 g/mol) (resulting concentration = 20 mM). Vortex until all solid has been thoroughly dissolved. The solution can be portioned into 50 μL aliquots as needed. These solutions must be stored in the freezer at < -20°C. Step 2: Add 50 μL of 20 mM stock to 150 μL DMSO (resulting concentration = 5 mM). This stock solution can be stored in the freezer at < -20°C. Step 3: The solution can be further diluted to the desired concentration in buffer. For example add 10 μL of the 5 mM stock solution to 990 μL 1X PBS, resulting in a 50 μM solution (this solution should be used the same day in which it was prepared). This will result in a solution containing only 1% DMSO (v/v). This solution can then be used to image H₂O₂ in the biological setting.</p>

Visit www.strem.com for new product information and a searchable catalog.

Strem Chemicals, Inc.
 7 Mulliken Way
 Newburyport, MA 01950
 U.S.A
 Tel: 978.499.1600
 Fax: 978.465.3104
 Email: info@strem.com

Strem Chemicals, Inc.
 15, rue de l'Atome
 Zone Industrielle
 67800 BISCHHEIM France
 Tel: (33) 03 88 62 52 60
 Fax: (33) 03 88 62 26 81
 Email: info.europe@strem.com

Strem Chemicals, Inc.
 Postfach 1215
 77672 KEHL
 Germany
 Tel: 0 78 51/ 7 58 79
 Email: info.europe@strem.com

Strem Chemicals UK Ltd.
 An Independent Distributor of Strem Chemicals Products
 Newton Hall, Town Street
 Newton, Cambridge
 England CB22 7ZE
 Tel: 0845 643 7263
 Fax: 0845 643 7362
 Email: enquiries@strem.co.uk