


27-3015



28-1720

Renewable energy sources such as photovoltaic power, wind power and biomass conversions to carbon-free energy carriers are growing sources for the decarbonization of the global energy systems. In this respect, hydrogen produced from renewable resources is a potentially promising source of clean energy.

Erwin Reisner and associates from the University of Cambridge have developed cobalt and nickel molecular catalysts that can be used for sunlight driven hydrogen production in water when immobilised on a semiconductor [1].

When semiconducting materials (e.g. TiO₂), co-functionalized with a catalyst and a dye (e.g. Ru or organic dye), are then exposed to light, the dye gets excited and oxidized. Generated electrons are transferred through the conduction band of the semiconductor and are injected directly into the H₂ evolution catalyst. The oxidized dye can in principle be regenerated by a water oxidation catalyst, thus closing the fuel-making cycle (Fig. 1).




Fig. 1. Schematic representation of dye-sensitized photocatalysis (DSP) for water splitting.
cat_{ox} = water oxidation catalyst, cat_{red} = H₂ evolution catalyst [1]

©[Chem. Soc. Rev., 2015, DOI: 10.1039/c5cs00733j] - Published by The Royal Society of Chemistry

27-3015: Triethylammonium {chlorobis(dimethylglyoximate)(4-hydrogenphosphonatepyridinyl)cobaltate(III)}, also known as a cobaloxime (CoP), is a catalyst that leads to the evolution of H₂ either electro- or photocatalytically in pH-neutral water, even in the presence of atmospheric O₂. The phosphonic acid group in (Et₃NH)[CoP] enables the complex to dissolve in water and allows for its immobilization on metal oxide surfaces for heterogeneous applications [2-4, 6, 9-10].

28-1720: Bis{P,P'-1,5-diphenyl-3,7-bis[(4-hydrogenphosphonate)phenyl]-1,5,3,7-diazadiphosphocine}nickel(II) bromide (hydrogen bromide adduct), also known as Ni bis(diphosphine)(NiP), is an alternative H₂-evolution catalyst. The NiP inner core of the complex is responsible for catalytic activity. The phosphonic acid group in the outer-sphere provides the catalyst with good solubility in aqueous solutions. In addition, the product can be anchored onto semiconductor particles [5-8].

As mentioned, the CoP catalyst shows a tolerance towards O₂, but its activity suffers in the presence of CO. In contrast, the NiP catalyst is strongly inhibited by O₂, but can produce H₂ under high CO concentrations [6].

Visit www.strem.com for new product information and a searchable catalog.

Strem Chemicals, Inc.
7 Mulliken Way
Newburyport, MA 01950
U.S.A
Tel: 978.499.1600
Fax: 978.465.3104
Email: info@strem.com

Strem Chemicals, Inc.
15, rue de l'Atome
Zone Industrielle
67800 BISCHHEIM France
Tel: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.
Postfach 1215
77672 KEHL
Germany
Tel: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.
An Independent Distributor of Strem Chemicals Products
Newton Hall, Town Street
Newton, Cambridge
England CB22 7ZE
Tel: 0845 643 7263
Fax: 0845 643 7362
Email: enquiries@strem.co.uk

27-3015 Triethylammonium {chlorobis(dimethylglyoximate)(4-hydrogenphosphonatepyridinyl) cobaltate(III)} 5mg
CAS# 1280199-86-1 25mg
 $C_{19}H_{35}ClCoN_6O_7P$; FW: 584.88; light-brown solid

Technical Notes:

1. Synthetic cobaloxime catalyst for reduction of aqueous protons to H_2 in the presence of atmospheric oxygen.
2. Catalysts contains phosphonate anchor groups for immobilisation on metal oxide semiconductor surface enabling light-driven H_2 evolution.

28-1720 Bis{P,P'-1,5-diphenyl-3,7-bis[(4-hydrogenphosphonate)phenyl]-1,5,3,7-diazadiphosphocine} nickel (II) bromide (hydrogen bromide adduct) 5mg
CAS# 1514896-39-9 25mg
 $C_{60}H_{69}Br_3N_4NiO_{12}P_8$; FW: 1584.41; red-brown xtls.

Technical Notes:

1. Synthetic Ni catalyst for reduction of aqueous protons to H_2 .
2. Catalyst tolerates CO.
3. Catalyst contains phosphonate anchor groups for immobilisation on metal oxide semiconductor surface enabling light-driven H_2 evolution.

References:

1. *Chem. Soc. Rev.*, **2015**, 45, 9.
2. *Angew. Chem. Int. Ed.*, **2012**, 51, 9381.
3. *Angew. Chem. Int. Ed.*, **2012**, 51, 12749.
4. *Chem. Commun.*, **2014**, 50, 12768.
5. *J. Am. Chem. Soc.*, **2014**, 136, 356.
6. *Chem. Commun.*, **2014**, 50, 15995.
7. *Angew. Chem. Int. Ed.*, **2014**, 53, 11538.
8. *J. Amer. Chem. Soc.*, **2015**, 137, 6018.
9. *Chem. Commun.*, **2011**, 47, 1695-1697
10. *Chem. Sci.*, **2015**, 6, 2727-2736

Visit www.strem.com for new product information and a searchable catalog.

Strem Chemicals, Inc.
7 Mulliken Way
Newburyport, MA 01950
U.S.A
Tel: 978.499.1600
Fax: 978.465.3104
Email: info@strem.com

Strem Chemicals, Inc.
15, rue de l'Atome
Zone Industrielle
67800 BISCHHEIM France
Tel: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.
Postfach 1215
77672 KEHL
Germany
Tel: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.
An Independent Distributor of Strem Chemicals Products
Newton Hall, Town Street
Newton, Cambridge
England CB22 7ZE
Tel: 0845 643 7263
Fax: 0845 643 7362
Email: enquiries@strem.co.uk