


# STREM Nitric Oxide Sensor (Extracellular) Kit ("NO-ON") (FL2A) (Cell-impermeable NO fluorescent probe)

metals • inorganics • organometallics • catalysts • ligands • custom synthesis • cGMP facilities • nanomaterials

## Kit Catalog number: 96-0397

**Active Ingredient:** 2-{4,5-Bis[(6-(2-ethoxy-2-oxoethoxy)-2-methylquinolin-8-ylamino)methyl]-6-hydroxy-3-oxo-3H-xanthen-9-yl}benzoic acid FL2A

**Catalog number - 07-0287**



**Introduction:** The copper complex of **FL2A** is a novel, cell-impermeable fluorescent NO probe, that allows direct imaging of nitric oxide produced through fluorescence turn on. A solution of the copper (II) complex of **FL2A** can be readily prepared using this kit. For additional information, consult the following references:

1. *Inorg. Chem.*, **2010**, *49*, 7464.
2. *PNAS*, **2010**, *107*, 8525.

**Contents:** Ligand **FL2A**: 5 x 0.5mg  
Dimethylsulfoxide (ACS spectrophotometric grade): 5 x 1.0ml  
Copper (II) chloride (1.0 mM solution in water): 5 x 1.0ml

**MSDS:** The Material Safety Data Sheets for the three products contained in this kit can be downloaded from the Strem Chemicals Web Site at [www.strem.com](http://www.strem.com). Locate the MSDS using the following catalog numbers:  
**FL2A**: 07-0291  
Dimethylsulfoxide (ACS spectrophotometric grade): 97-4940  
Copper (II) chloride as a 1.0 mM solution in water: 97-3060

**Storage conditions:** The kit should be stored at -20°C and protected from light.

**Preparation of the active copper complex of FL2A:** **Step I**  
Allow the kit to warm to room temperature. Add 589 microliters of DMSO to a 0.5mg vial of **FL2A** (resulting concentration - 1.0 mM). The **FL2A** is readily soluble in the DMSO. The solution can be partitioned into aliquots of 40-300 µL as required. These solutions must be stored in the freezer at <-20°C. The DMSO solution of **FL2A** is stable for three months at -80°C. It is advisable to check the extinction coefficient of the solution before preparing the copper complex. (log  $\epsilon$ (499 nm) = 4.66)

**Step II**  
A Cu<sub>2</sub>**FL2A** solution should be freshly prepared by adding 1:2 **FL2A** solution (1.0 mM) to the copper (II) solution (1.0 mM) at room temperature.  
Note: The prepared DMSO/water stock solution of Cu**FL2A** solution [log  $\epsilon$ (495 nm) = 4.19] can be kept at room temperature, but should be protected from light. The solution can be diluted with media to provide the concentration required for cell sensing experiments. When the extinction coefficient of the red solution of Cu**FL2A** diminishes by 20% of the original value, the solution should be discarded. Do not use the solution after 1 hour, and do not freeze the solution.

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.

**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950-4098  
U.S.A.  
Tel.: (978) 499-1600  
Fax: (978) 465-3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
15, rue de l'Atome  
Zone Industrielle  
67800 BISCHHEIM France  
Tel.: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Telefon: 0 78 51/ 7 58 79  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK Ltd.**  
An Independent Distributor of  
Strem Chemicals Products  
Newton Hall, Town Street  
Newton, Cambridge CB22 7ZE UK  
Tel.: 0845 643 7263  
Fax: 0845 643 7362  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)




# **STREM** Nitric Oxide Sensor (Intracellular) Kit ("NO-ON") (FL) (Cell-permeable NO fluorescent probe)

metals • inorganics • organometallics • catalysts • ligands • custom synthesis • cGMP facilities • nanomaterials

## Kit Catalog number: 96-0293

**Active Ingredient:** 2-{2-Chloro-6-hydroxy-5-[2-methylquinolin-8-ylamino)methyl]-3-oxo-3H-xanthen-9-yl}benzoic acid (FL)

Catalog number - 07-0293



**Introduction:** The copper complex of FL is a novel, cell-permeable fluorescent NO probe that allows direct imaging of nitric oxide produced in living cells through fluorescence turn on. The probe is best suitable for immediate intracellular detection of NO. Imaging of the cells with CuFL be performed at the latest one hour after incubation. A solution of the copper (II) complex of FL can be readily prepared using this kit. For additional information, consult the following references:

1. *Nature Chemical Biology*, **2006**, 2, 375.
2. *Nature Protocols*, **2007**, 2, 408.
3. *J. Am. Chem. Soc.*, **2006**, 128, 14364.

**Contents:** Ligand FL: 5 x 0.5mg  
Dimethylsulfoxide (ACS spectrophotometric grade): 5 x 1.0ml  
Copper (II) chloride (1.0 mM solution in water): 5 x 1.0ml

**MSDS:** The Material Safety Data Sheets for the three products contained in this kit can be downloaded from the Strem Chemicals Web Site at [www.strem.com](http://www.strem.com). Locate the MSDS using the following catalog numbers:  
FL: 07-0293  
Dimethylsulfoxide (ACS spectrophotometric grade): 97-4940  
Copper (II) chloride as a 1.0 mM solution in water: 97-3060

**Storage conditions:** The kit should be stored at -20°C and protected from light.

**Preparation of the active copper complex of FL:**

**Step I**

Allow the kit to warm to room temperature. Add 931 microliters of DMSO to a 0.5mg vial of FL (resulting concentration - 1.0 mM). The FL is readily soluble in the DMSO. The solution can be partitioned into aliquots of 40-300 µL as required. These solutions must be stored in the freezer at <-20°C. The DMSO solution of FL is stable for three months at -80°C. It is advisable to check the extinction coefficient of the solution before preparing the copper complex.  
(log  $\epsilon$ (504 nm) = 4.62)

**Step II**

A CuFL solution should be freshly prepared by adding 1:1 FL solution (1.0 mM) to the copper (II) solution (1.0 mM) at room temperature.  
Note: The prepared DMSO/water stock solution of CuFL solution [log  $\epsilon$ (499 nm) = 4.60] can be kept at room temperature, but should be protected from light. The solution can be diluted with media to provide the concentration required for cell sensing experiments. When the extinction coefficient of the red solution of CuFL diminishes by 20% of the original value, the solution should be discarded. Do not use the solution after 1 hour, and do not freeze the solution.

metals • inorganics • organometallics • catalysts • ligands • custom synthesis • cGMP facilities • nanomaterials

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.

**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950-4098  
U.S.A.  
Tel.: (978) 499-1600  
Fax: (978) 465-3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
15, rue de l'Atome  
Zone Industrielle  
67800 BISCHHEIM France  
Tel.: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Telefon: 0 78 51/ 7 58  
79  
Email:  
[info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK Ltd.**  
An Independent Distributor of  
Strem Chemicals Products  
Newton Hall, Town Street  
Newton, Cambridge CB22 7ZE UK  
Tel.: 0845 643 7263  
Fax: 0845 643 7362  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)




# STREM Nitric Oxide Sensor (Intracellular) Kit ("NO-ON") (FL2E) (Cell-trappable NO fluorescent probe)

metals • inorganics • organometallics • catalysts • ligands • custom synthesis • cGMP facilities • nanomaterials

Kit Catalog number: 96-0396

Active Ingredient: 2-{4,5-Bis[(6-(2-ethoxy-2-oxoethoxy)-2-methylquinolin-8-ylamino)methyl]-6-hydroxy-3-oxo-3H-xanthen-9-yl}benzoic acid FL2E

Catalog number - 07-0291



**Introduction:** The copper complex of **FL2E** is a novel, cell-trappable fluorescent NO probe, that allows direct imaging of nitric oxide produced in living cells through fluorescence turn on. The FL2E ligand is trapped in the cell. A solution of the copper (II) complex of **FL2E** can be readily prepared using this kit.  
For additional information, consult the following references:  
1. *Inorg. Chem.*, **2010**, *49*, 7464.

**Contents:** Ligand **FL2E**: 5 x 0.5mg  
Dimethylsulfoxide (ACS spectrophotometric grade): 5 x 1.0ml  
Copper (II) chloride (1.0 mM solution in water): 5 x 1.0ml

**MSDS:** The Material Safety Data Sheets for the three products contained in this kit can be downloaded from the Strem Chemicals Web Site at [www.strem.com](http://www.strem.com). Locate the MSDS using the following catalog numbers:  
**FL2E**: 07-0291  
Dimethylsulfoxide (ACS spectrophotometric grade): 97-4940  
Copper (II) chloride as a 1.0 mM solution in water: 97-3060

**Storage conditions:** The kit should be stored at -20°C and protected from light.

**Preparation of the active copper complex of FL2E:** **Step I**  
Allow the kit to warm to room temperature. Add 571 microliters of DMSO to a 0.5mg vial of **FL2E** (resulting concentration - 1.0 mM). The **FL2E** is readily soluble in the DMSO. The solution can be partitioned into aliquots of 40-300 µL as required. These solutions must be stored in the freezer at <-20°C. The DMSO solution of **FL2E** is stable for three months at -80°C. It is advisable to check the extinction coefficient of the solution before preparing the copper complex. (log ε(500 nm) = 4.25)

**Step II**

A Cu<sub>2</sub>**FL2E** solution should be freshly prepared by adding 1:2 **FL2E** solution (1.0 mM) to the copper (II) solution (1.0 mM) at room temperature.

Note: The prepared DMSO/water stock solution of Cu**FL2E** solution [log ε(496 nm) = 4.05] can be kept at room temperature, but should be protected from light. The solution can be diluted with media to provide the concentration required for cell sensing experiments. When the extinction coefficient of the red solution of Cu**FL2E** diminishes by 20% of the original value, the solution should be discarded. Do not use the solution after 1 hour, and do not freeze the solution.

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.

**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950-4098  
U.S.A.  
Tel.: (978) 499-1600  
Fax: (978) 465-3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
15, rue de l'Atome  
Zone Industrielle  
67800 BISCHEIM France  
Tel.: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Telefon: 0 78 51/ 7 58 79  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK Ltd.**  
An Independent Distributor of  
Strem Chemicals Products  
Newton Hall, Town Street  
Newton, Cambridge CB22 7ZE UK  
Tel.: 0845 643 7263  
Fax: 0845 643 7362  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)

