



STREM Phosphines for Colloidal Synthesis of Semiconducting Quantum Dots

METALS • INORGANICS • ORGANOMETALLICS • CATALYSTS • LIGANDS • NANOMATERIALS • CUSTOM SYNTHESIS • cGMP FACILITIES



Tri-n-octylphosphine (TOP) [15-6655] and Trioctylphosphine oxide (TOPO) [15-6660/15-6661] are used as surfactant/capping agents for the colloidal synthesis of CdE (E=S, Se, Te) type semiconductor nanocrystals (Quantum Dots).^[1] TOP and TOPO prevent nanoparticles from aggregation and uncontrolled growth as well as stabilize them in solution. In addition, TOPO serves as a reaction medium. Its high boiling temperature (411.2°C) facilitates reactions leading to homogeneous nucleation and controlled growth of the particles. TOPO/TOP coated Quantum Dots (QDs) are compatible with organic solvents such as chloroform or toluene.

TOP with a lone electron pair on the phosphorus atom serves as a common reagent in the chemical synthesis of II-VI, III-V and IV-VI QDs.^[1] These pair electrons exhibit strong solvent ability toward selenium (93-3416; 34-0090), sulfur (93-1616; 93-1617; 93-1618) and tellurium (93-5222; 52-5200; 93-5220) powders forming thus TOP-E anion precursors for selenation, sulfonation and telluration of metal cations.

In TOP/TOPO media, CdSe QDs are synthesized starting with Me_2Cd (48-5040) and $(TMS)_2Se$.^[1] For selenation, use a TOPSe solution prepared by mixing elemental Se (93-3416 or 34-0090) and TOP.^[2] Alternatively, CdSe QDs can be prepared in TOPO by injection of $Me_2Cd:Se$ powder solution in tributylphosphine^[3] or using CdO (93-4835, 93-4817, 48-0800) as a precursor.^[4] TOPO is used for the preparation of core-shell structured QDs, e.g. CdSe/ZnS^[5] or CdSe/CdS/ZnS^[6] and Perovskite type QDs.^[7]

For other long-chained alkylphosphonic acids used in the colloidal synthesis of QDs, or other types of nanoparticles, please review the following Strem product numbers: 96-1525 (Kit); 15-0958; 15-1835; 15-2400; 15-2410; 15-3510; 15-3520; 15-5145.

References:

1. *J. Am. Chem. Soc.*, **1993**, 115, 8706.
2. *J. Mater. Chem. B*, **2013**, 1, 1381.
3. *Nature*, **2002**, 420, 800.
4. *J. Am. Chem. Soc.*, **2001**, 123, 183.
5. *Nano Lett.*, **2001**, 1, 207.
6. *J. Phys. Chem. B*, **2004**, 108, 18826.
7. *ACS Nano*, **2017**, 11, 10373.

Visit www.strem.com for new product information and a searchable catalog.

Strem Chemicals, Inc.
7 Mulliken Way
Newburyport, MA 01950
U.S.A
Tel: 978.499.1600
Fax: 978.465.3104
Email: info@strem.com

Strem Chemicals, Inc.
15, rue de l'Atome
Zone Industrielle
67800 BISCHHEIM France
Tel: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.
Postfach 1215
77672 KEHL
Germany
Tel: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.
An Independent Distributor of Strem Chemicals Products
Newton Hall, Town Street
Newton, Cambridge
England CB22 7ZE
Tel: 0845 643 7263
Fax: 0845 643 7362
Email: enquiries@strem.co.uk