



What is Laser Ablation?

Laser ablation is a physical process used to generate nanoparticles. It involves the use of short pulses of laser energy focused on a target in a solvent. The target absorbs the energy from the laser pulse and is vaporized. The vaporized material then condenses as **nanoparticles**.

Why Laser Ablation?

Laser ablation generates **nanoparticle** dispersions that are free of any contaminants, such as unreacted starting materials.

The use of **short laser pulses** makes it possible to carry out the **ablation process** in a wide range of carriers, including volatile organic solvents and reactive monomers. The laser technology also allows the conversion of almost any solid or powder material into a colloid. And that means that nearly endless material combinations are possible by varying the target and the solvent for screening studies and individual applications.




In addition, the released **nanoparticles** have a positive charge on the surface, and thus electrostatic stabilization by the surrounding medium occurs without the need, in most cases, for any type of chemical stabilizer.

The positive surface charge also provides a means for interaction with molecules that bear electron-donating substituents. In-situ functionalization is also possible with monomers, polymers and biomolecules. Deposition can also be achieved using electrophoresis.

In fact, **laser-generated nanoparticles** can be readily and permanently immobilized on micro- and nanoparticulate supports under ambient conditions without heat treatment. And because the **nanoparticles** are not covered with hydrophobic ligands, even water-insoluble supports like oxides, sulfates and phosphates can be used. Their high purity also means that supported catalysts prepared from laser-generated **nanoparticles** don't suffer from poisoning.

Advantages of Laser Ablation Generated Nanoparticles

- ❖ Very high purity – free of particle-associated residual chemicals
- ❖ High stability due to electrostatic stabilization (12 mo., 2mo. in H₂O)
- ❖ High dispersibility
- ❖ Direct use of organic solvents
- ❖ Adjustable ion release rates
- ❖ Increase adsorption efficiency and reduced cost
- ❖ High bio-compatibility in different applications
- ❖ Higher reactivity and lower toxicity for medical/pharma applications



Courtesy University of Duisburg-Essen.



Possible Applications for Ablation Generated Nanoparticles

- ❖ Biomedical, immunoassay, and toxicology studies
- ❖ Targeted drug delivery
- ❖ Supported heterogeneous catalysts
- ❖ Nanocoatings
- ❖ Nanocomposites

Visit www.strem.com for new product information and searchable catalog.

Strem Chemicals, Inc.

7 Mulliken Way
Newburyport, MA 01950-4098
U.S.A.
Tel.: (978) 499-1600
Fax: (978) 465-3104
Email: info@strem.com
[particular_nano0213.doc](#)

Strem Chemicals, Inc.

15, rue de l'Atome
Zone Industrielle
67800 BISCHHEIM France
Tel.: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.

Postfach 1215
77672 KEHL
Germany
Telefon: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.

41 Hills Road
Cambridge
England CB2 1NT
Tel.: 0845 643 7263
Fax: 0845 643 7362
Email: enquiries@strem.co.uk



Selection Matrix for Nanoparticles prepared using Laser Ablation

Strem offers a wide range of residue-free **nanoparticles** prepared via **Laser Ablation** by Particular GmbH.

Selection Matrix- pick element, size, concentration, solvent from available list

| Element | Size | | Concentration | | Solvent | | | |
|---------|-------|---------|---------------|---------|---------|---------|-----------------|--------------|
| | <20nm | 50-70nm | 100mg/L | 500mg/L | water | acetone | Ethylene Glycol | iso propanol |
| Ag | x | x | x | | x | x | | 50-70nm only |
| | x | x | | x | x | | | |
| Au | x | x | x | | x | x | x | x |
| | x | x | | x | x | | | |
| Pt | x | x | x | | x | x | 50-70nm only | x |
| | x | x | | x | x | | | |
| Pd | x | x | x | | x | x | 50-70nm only | |
| | x | x | | x | x | | | |
| Rh | x | x | x | | x | x | 50-70nm only | x |
| | x | x | | x | x | | | |
| Ru | x | x | x | | x | x | 50-70nm only | x |
| | x | x | | x | x | | | |
| Cu | x | x | | | | x | 50-70nm only | |
| Ti | x | x | | | | x | x | 50-70nm only |

x = available Stability: 12 months (2 month in water)



Nanoparticles from metals (silver, titanium, platinum, gold) in water

Visit www.strem.com for new product information and searchable catalog.

Strem Chemicals, Inc.
 7 Mulliken Way
 Newburyport, MA 01950-4098
 U.S.A.
 Tel.: (978) 499-1600
 Fax: (978) 465-3104
 Email: info@strem.com
 particular_nano0213.doc

Strem Chemicals, Inc.
 15, rue de l'Atome
 Zone Industrielle
 67800 BIRSCHHEIM France
 Tel.: (33) 03 88 62 52 60
 Fax: (33) 03 88 62 26 81
 Email: info.europe@strem.com

Strem Chemicals, Inc.
 Postfach 1215
 77672 KEHL
 Germany
 Telefon: 0 78 51/ 7 58 79
 Email: info.europe@strem.com

Strem Chemicals UK Ltd.
 41 Hills Road
 Cambridge
 England CB2 1NT
 Tel.: 0845 643 7263
 Fax: 0845 643 7362
 Email: enquiries@strem.co.uk



Reactant & Surfactant-free pure Nanoparticles via Laser Ablation



Particular
customized material

Cat. No. Description

| | |
|----------------|--|
| 79-0410 | Gold nanoparticles, pure, (<20nm) in water at 100mg/L (surfactant and reactant-free, OD>1, stabilized with < 0.01 mmol/l of citrate) |
| 79-0412 | Gold nanoparticles, pure, (<20nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 79-0414 | Gold nanoparticles, pure, (<20nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 79-0416 | Gold nanoparticles, pure, (<20nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 79-0418 | Gold nanoparticles, pure, (<20nm) in water at 500mg/L (surfactant and reactant-free, OD>5, stabilized with < 0.01 mmol/l of citrate) |
| 79-0420 | Gold nanoparticles, pure, (50-70nm) in water at 100mg/L (surfactant and reactant-free, OD>1, stabilized with < 0.01 mmol/l of citrate) |
| 79-0422 | Gold nanoparticles, pure, (50-70nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 79-0424 | Gold nanoparticles, pure, (50-70nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 79-0426 | Gold nanoparticles, pure, (50-70nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 79-0428 | Gold nanoparticles, pure, (50-70nm) in water at 500mg/L (surfactant and reactant-free, OD>5, stabilized with < 0.01 mmol/l of citrate) |
| 47-0710 | Silver nanoparticles, pure, (<20nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 47-0712 | Silver nanoparticles, pure, (<20nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 47-0718 | Silver nanoparticles, pure, (<20nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 47-0720 | Silver nanoparticles, pure, (50-70nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 47-0722 | Silver nanoparticles, pure, (50-70nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 47-0726 | Silver nanoparticles, pure, (50-70nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 47-0728 | Silver nanoparticles, pure, (50-70nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 78-1402 | Platinum nanoparticles, pure, (<20nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 78-1404 | Platinum nanoparticles, pure, (<20nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 78-1408 | Platinum nanoparticles, pure, (<20nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 78-1410 | Platinum nanoparticles, pure, (<20nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 78-1412 | Platinum nanoparticles, pure, (50-70nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 78-1414 | Platinum nanoparticles, pure, (50-70nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 78-1416 | Platinum nanoparticles, pure, (50-70nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 78-1418 | Platinum nanoparticles, pure, (50-70nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 78-1422 | Platinum nanoparticles, pure, (50-70nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |

Visit www.strem.com for new product information and searchable catalog.

Strem Chemicals, Inc.

7 Mulliken Way
Newburyport, MA 01950-4098
U.S.A.
Tel.: (978) 499-1600
Fax: (978) 465-3104
Email: info@strem.com
[particular_nano0213.doc](#)

Strem Chemicals, Inc.

15, rue de l'Atome
Zone Industrielle
67800 BISCHHEIM France
Tel.: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.

Postfach 1215
77672 KEHL
Germany
Telefon: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.

41 Hills Road
Cambridge
England CB2 1NT
Tel.: 0845 643 7263
Fax: 0845 643 7362
Email: enquiries@strem.co.uk

Strem Nanoparticles prepared via Laser Ablation

| Cat. No. | Description |
|----------|---|
| 46-4010 | Palladium nanoparticles, pure, (<20nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 46-4012 | Palladium nanoparticles, pure, (<20nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 46-4018 | Palladium nanoparticles, pure, (<20nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 46-4020 | Palladium nanoparticles, pure, (50-70nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 46-4022 | Palladium nanoparticles, pure, (50-70nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 46-4024 | Palladium nanoparticles, pure, (50-70nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 46-4028 | Palladium nanoparticles, pure, (50-70nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 45-1322 | Rhodium nanoparticles, pure, (<20nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 45-1324 | Rhodium nanoparticles, pure, (<20nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 45-1328 | Rhodium nanoparticles, pure, (<20nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 45-1330 | Rhodium nanoparticles, pure, (<20nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 45-1332 | Rhodium nanoparticles, pure, (50-70nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 45-1334 | Rhodium nanoparticles, pure, (50-70nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 45-1336 | Rhodium nanoparticles, pure, (50-70nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 45-1338 | Rhodium nanoparticles, pure, (50-70nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 45-1340 | Rhodium nanoparticles, pure, (50-70nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 44-2810 | Ruthenium nanoparticles, pure, (<20nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 44-2812 | Ruthenium nanoparticles, pure, (<20nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 44-2816 | Ruthenium nanoparticles, pure, (<20nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 44-2818 | Ruthenium nanoparticles, pure, (<20nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 44-2820 | Ruthenium nanoparticles, pure, (50-70nm) in water at 100mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 44-2822 | Ruthenium nanoparticles, pure, (50-70nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 44-2824 | Ruthenium nanoparticles, pure, (50-70nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 44-2826 | Ruthenium nanoparticles, pure, (50-70nm) in isopropanol at 100mg/L (surfactant and reactant-free) |
| 44-2828 | Ruthenium nanoparticles, pure, (50-70nm) in water at 500mg/L (surfactant and reactant-free, stabilized with < 0.01 mmol/l of citrate) |
| 29-0092 | Copper nanoparticles, pure, (<20nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 29-0096 | Copper nanoparticles, pure, (50-70nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 29-0098 | Copper nanoparticles, pure, (50-70nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 22-0192 | Titanium nanoparticles, pure, (<20nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 22-0194 | Titanium nanoparticles, pure, (<20nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 22-0196 | Titanium nanoparticles, pure, (50-70nm) in acetone at 100mg/L (surfactant and reactant-free) |
| 22-0198 | Titanium nanoparticles, pure, (50-70nm) in ethylene glycol at 100mg/L (surfactant and reactant-free) |
| 22-0203 | Titanium nanoparticles, pure, (50-70nm) in isopropanol at 100 mg/L (surfactant and reactant-free) |

References:

1. Wagener, P.; Barcikowski, S.; Bärsch, N.: Fabrication of nanoparticles and nanomaterials using laser ablation in liquids. In: [Photonik international 1/2011](#).
2. Wagener, P.; Schwenke, A.; Barcikowski, S.: How Citrate Ligands Affect Nanoparticle Adsorption to Microparticle Supports. In: [Langmuir 28 \(2012\), pp. 6132-6140](#).
3. Sajti, C. L.; Sattari, R.; Chichkov, B. N.; Barcikowski, S.: Gram Scale Synthesis of Pure Ceramic Nanoparticles by Laser Ablation in Liquid. In: [Journal of Physical Chemistry C 114 \(2010\), 6, pp. 2421-2427](#).
4. Petersen, S.; Barcikowski, S.: Conjugation Efficiency of Laser-Based Bioconjugation of Gold Nanoparticles with Nucleic Acids. In: [The Journal of Physical Chemistry C 113 \(2009\), 46, pp. 19830-19835](#).

Visit www.strem.com for new product information and searchable catalog.

Strem Chemicals, Inc.

7 Mulliken Way
Newburyport, MA 01950-4098
U.S.A.
Tel.: (978) 499-1600
Fax: (978) 465-3104
Email: info@strem.com
[particular_nano0213.doc](#)

Strem Chemicals, Inc.

15, rue de l'Atome
Zone Industrielle
67800 BISCHHEIM France
Tel.: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.

Postfach 1215
77672 KEHL
Germany
Telefon: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.

41 Hills Road
Cambridge
England CB2 1NT
Tel.: 0845 643 7263
Fax: 0845 643 7362
Email: enquiries@strem.co.uk