


STREM Catkit Single Use Vials for low catalyst loading experiments

metals · inorganics · organometallics · catalysts · ligands · custom synthesis · cGMP facilities · nanomaterials

Limited substrate availability, expense, convenience or safety may drive the need or preference for small scale reactions. However, weighing less than 1mg of catalysts and ligands for such small scale reactions can be difficult, especially if precise ratios of metal precursor and ligand are required eg. Suzuki Miyaura Coupling.




The **Table 1** below shows the quantity of metal precursor and ligand needed for a 2mol% loading of catalyst using 25mg of substrate with a molecular weight of 350.

Metal Precursor	Weight	Ligand	Weight
Pd ₂ (dba) ₃	0.65mg	SPhos	0.59mg
Pd(OAc) ₂	0.32mg	P(t-Bu) ₃	0.29mg
CuI	0.27mg	dtbpf	0.68mg

These quantities cannot be consistently weighed out. The practical issue of accurately weighing out small quantities of components can be resolved by the use of Admixtures (See Sharpless AD-Mix, *J. Org. Chem.*, **1992**, 57, 2771-3)

The metal precursor and ligand can be "diluted" with a suitable base by grinding the mixture of components until homogeneous. This mixture can then be easily weighed out into single use vials so that each contains a known but low catalyst loading.



Kit of Catkits

96-3790 **NEW** → **Kit of Catkits – Single-Use Vials for low catalyst loading experiments**
The kit contains 4 different 5 x 1 vials of the items listed below.
Components available for individual sale.

Available Catkits

46-2040 NEW →	trans-Dichlorobis(tricyclohexylphosphino)palladium(II)/potassium phosphate admixture [Catkit single-use vials – 6.62 wt% Pd complex] [29934-17-6] PdCl ₂ [(C ₆ H ₁₁) ₃ P] ₂ /K ₃ PO ₄ ; off-white powder. Note: Each vial contains 453mg of admixture, which will deliver 4 mole% of palladium catalyst and 2 equivalents of base, to a reaction using 1 mmole of substrate. Weight-percent of components: 6.62wt% palladium complex 93.38 wt% potassium phosphate	5x1 vial 25x1 vial
46-2038 NEW →	trans-Dichlorobis(triphenylphosphine)palladium(II)/potassium phosphate admixture [Catkit single-use vials – 6.32wt % complex] [13965-03-2] PdCl ₂ [P(C ₆ H ₅) ₃] ₂ /K ₃ PO ₄ ; off-white powder. Note: Each vial contains 453mg of admixture, which will deliver 4 mole% of palladium catalyst and 2 equivalents of base, to a reaction using 1 mmole of substrate. Weight-percent of components: 6.32wt% palladium complex 93.68 wt% potassium phosphate	5x1 vial 25x1 vial
46-2030 NEW →	Palladium(II) acetate/1,1'-bis(di-t-butylphosphino)ferrocene/potassium phosphate admixture [Catkit single-use vials – 2.02 wt% Pd(OAc) ₂] [95408-45-0] off-white powder. Note: Each vial contains 453mg of admixture, which will deliver 4 mole% of palladium catalyst and 2 equivalents of base, to a reaction using 1 mmole of substrate. Weight-percent of components: 2.02 wt% palladium acetate 4.28 wt% phosphine ligand 93.7 wt% potassium phosphate	5x1 vial 25x1 vial
46-2033 NEW →	Palladium(II) acetate/2-dicyclohexylphosphino-2,6-dimethoxy-1,1'-biphenyl (SPhos)/potassium phosphate admixture [Catkit single-use vials – 1.96 wt% Pd(OAc) ₂] [1028206-58-7] off-white powder. Note: Each vial contains 453mg of admixture, which will deliver 4 mole% of palladium catalyst and 2 equivalents of base, to a reaction using 1 mmole of substrate. Weight-percent of components: 1.96wt% palladium acetate 7.17 wt% phosphine ligand 90.8 wt% potassium phosphate	5x1 vial 25x1 vial

For detailed description of the Catalyst Kit approach ("Catkit") see J.C. Ruble "Making Catalysis Simple", *New Horizons in Catalysis*, 4 December 2009.

Also Available

29-0490 NEW →	Copper(I) iodide/cesium carbonate admixture [5.50 wt% CuI] [7681-65-4] CuI/Cs ₂ CO ₃ ; off-white powder. Weight-percent of components: 5.50 wt% copper(I) iodide 94.50 wt% cesium carbonate	5g 25g
--------------------------------	--	-----------

Visit www.strem.com for new product information and searchable catalog.

Strem Chemicals, Inc.
7 Mulliken Way
Newburyport, MA 01950-4098
U.S.A.
Tel.: (978) 499-1600
Fax: (978) 465-3104
Email: info@strem.com

Strem Chemicals, Inc.
15, rue de l'Atome
Zone Industrielle
67800 BIRSCHHEIM France
Tel.: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.
Postfach 1215
77672 KEHL
Germany
Telefon: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK, Ltd.
41 Hills Road
Cambridge
England CB2 1NT
Tel.: 0845 643 7263
Fax: 01223 368021
Email: enquiries@strem.co.uk