






Catalog # 05-0175 CALLERY™ Borane tetrahydrofuran, 1M in tetrahydrofuran, BTHF




Borane-tetrahydrofuran complex (BTHF) is a valuable reagent for the reduction of functional groups and for hydroboration reactions with carbon-carbon double and triple bonds. Functional groups that are readily reduced by BTHF include aldehyde, ketone, carboxylic acid, amide, oxime, imine, and nitrile. The carboxylic acid group is reduced at a faster rate than most groups including non-conjugated alkene. Conjugated α,β -unsaturated carboxylic acids give saturated alcohols as the major products [1, Review].


Technical Notes:

- Used in Au-catalyzed intermolecular hydroaminations with less-reactive internal alkynes and unprotected aliphatic amines giving excellent yields with low catalyst loading.
- Reductant used in In-catalyzed reduction of aromatic and aliphatic nitriles to primary amines.
- Dopant/Reducing agent used for generation of B-doped graphene nanoplatelets (borane-reduced graphene oxide, B-rG-O) using a solution process.
- Used for the reduction of functionalized tertiary phosphine oxides with BH_3 .
- Component of *N*-heterocyclic carbene borane organocatalyst used for asymmetric reduction of tert-butanesulfinyl ketimines.
- Catalyst used for the hydroboration of alkynes and alkenes.
- Used for the modification of alkoxy-functionalized silicon QD surfaces via ligand exchange of the alkoxy-surface groups with alkyl or alkenyl-surface groups.
- Ligand component for Cu-catalyzed transfer hydrogenation of *N*-heteroaromatics with an oxazaborolidine complex.
- Used in Ir-catalyzed photocatalytic C–F bond borylation of polyfluoroarenes with $\text{NHC}\cdot\text{BH}_3$.
- Reductant used in functionalization of [2.2]paracyclophanes via a reductive sulfanylation reaction.


Tech Note (1)
Ref. (2)Tech Note (2)
Ref. (3)Tech Note (3)
Ref. (4)Tech Note (4)
Ref. (5)




Tech Note (5)
Ref. (6)




Tech Note (6)
Ref. (7)




Tech Note (7)
Ref. (8)



Tech Note (8)
Ref. (9)



Tech Note (9)
Ref. (10)



Tech Note (10)
Ref. (11)

References:

1. [Chem. Rev. 2006, 106, 2617.](#)
2. [J. Am. Chem. Soc. 2009, 131, 12100.](#)
3. [J. Org. Chem. 2012, 77, 221.](#)
4. [ACS Nano, 2013, 7, 19.](#)
5. [J. Org. Chem. 2015, 80, 1672.](#)
6. [J. Org. Chem. 2015, 80, 11441.](#)
7. [J. Am. Chem. Soc. 2016, 138, 7114.](#)
8. [Synthesis 2018, 50, 803.](#)
9. [ACS Omega 2019, 4, 8487.](#)
10. [Org. Lett. 2020, 22, 1742.](#)
11. [J. Org. Chem. 2021, 86, 507.](#)