







BINAP Catalyst & Ligand Kits for asymmetric hydrogenation & other catalytic applications

BINAP is a chiral ligand that has been used in a variety of asymmetric reactions since 1980. Especially BINAP-Ru complex affords high enantioselectives in asymmetric hydrogenations of a wide range of olefins and functionalized ketones.

ToIBINAP has high solubility in organic solvents. This

property solves the problems caused by low solubility intrinsic to BINAP and its complexes. **XyIBINAP** is usually used in the Noyori's Ru-Diphosphine-Diamine complexes for hydrogenation of simple ketones. In many cases this ligand gives higher ee values than those with BINAP or ToIBINAP. **H₈-BINAP**-Ru complexes hydrogenate unsaturated carboxylic acids in higher enantioselectives than BINAP-Ru complexes do.



These kits contain 250mg of each of the items listed below it. All items are available for individual sale in 250mg, 1g & 5g unit sizes

96-6950 Takasago BINAP Ligand Kit

Item#	Acronym	CAS#	Item#	Acronym	CAS#
15-0150	(R)-(+)-BINAP	76189-55-4	15-0151	(S)-(-)-BINAP	76189-56-5
15-0152	(R)-(+)-ToIBINAP	99646-28-3	15-0153	(S)-(-)-ToIBINAP	100165-88-6
15-0476	(R)-(+)-XyIBINAP	137219-86-4	15-0477	(S)-(-)-XyIBINAP	135139-00-3

96-6953 Takasago BINAP Ru Acetate Catalyst Kit

Metal catalysts in this kit are from the following four catalyst types:

[RuCl(p-cymene)(P^P)]Cl (I)	[NH ₂ Me ₂][{RuCl(P^P)} ₂ (μ-Cl) ₃] (II)	[Ru(OAc) ₂ (P^P)] (III)	[RuCl ₂ (P^P)(N^N)] (IV)		
Item#	Acronym (Ligand P^P)	CAS#	Item#	Acronym (Ligand P^P)	CAS#
[Ru(OAc) ₂ (P^P)](III)					
44-0152	(R)-BINAP	325146-81-4	44-0153	(S)-BINAP	261948-85-0
44-0162	(R)-ToIBINAP	116128-29-1	44-0163	(S)-ToIBINAP	106681-15-6
44-0164	(R)-XyIBINAP	374067-50-2	44-0165	(S)-XyIBINAP	374067-49-9

96-6951 Takasago BINAP Ru Cymene Catalyst Kit

Item#	Acronym (Ligand P^P)	CAS#	Item#	Acronym (Ligand P^P)	CAS#
[RuCl[(p-cymene)(P^P)]Cl(I)					
44-0084	(R)-BINAP	145926-28-9	44-0086	(S)-BINAP	130004-33-0
44-0088	(R)-ToIBINAP	131614-43-2	44-0089	(S)-ToIBINAP	228120-95-4
44-0092	(R)-XyIBINAP	944451-24-5	44-0093	(S)-XyIBINAP	944451-25-6

Sold in collaboration with Takasago

Visit www.strem.com for new product information and a searchable catalog.

Strem Chemicals, Inc.
 7 Mulliken Way
 Newburyport, MA 01950
 U.S.A.
 Tel: 978.499.1600
 Fax: 978.465.3104
 Email: info@strem.com

Strem Chemicals, Inc.
 15, rue de l'Atome
 Zone Industrielle
 67800 BISCHHEIM France
 Tel: (33) 03 88 62 52 60
 Fax: (33) 03 88 62 26 81
 Email: info.europe@strem.com

Strem Chemicals, Inc.
 Postfach 1215
 77672 KEHL
 Germany
 Tel: 0 78 51/ 7 58 79
 Email: info.europe@strem.com

Strem Chemicals UK Ltd.
 An Independent Distributor of Strem Chemicals Products
 Newton Hall, Town Street
 Newton, Cambridge
 England CB22 7ZE
 Tel: +44 (0)1223 873 028
 Fax: +44 (0)1223 870207
 Email: enquiries@strem.co.uk

96-6954 Takasago BINAP Ru Diamine Catalyst Kit

Item#	Acronym (Ligand P [∧] P)	CAS#	Item#	Acronym (Ligand P [∧] P)	CAS#
<i>[RuCl₂(P[∧]P)(N[∧]N)](IV)</i>					
44-0212	RuCl ₂ [(R)-XylBINAP][(R)-daipen]	220114-32-9	44-0213	RuCl ₂ [(S)-XylBINAP][(S)-daipen]	220114-01-2
44-0217	(R)-RUCY [®] XylBINAP	1384974-38-2	44-0218	(S)-RUCY [®] XylBINAP	1312713-89-5
44-0226	RuCl ₂ [(S)-XylBINAP][(S,S)-dpen]	220114-38-5	44-0224	RuCl ₂ [(S)-XylBINAP][(S,S)-dpen]	220114-03-4

96-6952 Takasago BINAP Ru Dimer Catalyst Kit

Item#	Acronym (Ligand P [∧] P)	CAS#	Item#	Acronym (Ligand P [∧] P)	CAS#
<i>[NH₂Me₂][{RuCl(P[∧]P)}₂(μ-Cl)₃]Cl (II)</i>					
44-0510	(R)-BINAP	199684-47-4	44-0511	(S)-BINAP	199541-17-8
44-0512	(R)-TolBINAP	749935-02-2	44-0513	(S)-TolBINAP	309735-86-2
44-0514	(R)-XylBINAP	944451-08-5	44-0515	(S)-XylBINAP	944451-10-9

H₈-BINAP - Also Available (not included in the above kits) - Unit Sizes: 50mg, 250mg

15-2972	(R)-(+)-H ₈ -BINAP	139139-86-9	15-2973	(S)-(-)-H ₈ -BINAP	139139-93-8
44-0094	[RuCl(p-cymene)((R)-H ₈ -BINAP)]Cl	944451-26-7	44-0095	[RuCl(p-cymene)((S)-H ₈ -BINAP)]Cl	944451-27-8
44-0516	[NH ₂ Me ₂][{RuCl((R)-H ₈ -BINAP)} ₂ (μ-Cl) ₃]	204933-84-6	44-0517	[NH ₂ Me ₂][{RuCl((S)-H ₈ -BINAP)} ₂ (μ-Cl) ₃]	944451-12-1
44-0166	Ru(OAc) ₂ [(R)-H ₈ -BINAP]	374067-51-3	44-0167	Ru(OAc) ₂ [(S)-H ₈ -BINAP]	142962-95-6

Visit www.strem.com for new product information and a searchable catalog.

Strem Chemicals, Inc.
7 Mulliken Way
Newburyport, MA 01950
U.S.A
Tel: 978.499.1600
Fax: 978.465.3104
Email: info@strem.com

Strem Chemicals, Inc.
15, rue de l'Atome
Zone Industrielle
67800 BISCHEIM France
Tel: (33) 03 88 62 52 60
Fax: (33) 03 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.
Postfach 1215
77672 KEHL
Germany
Tel: 0 78 51/ 7 58 79
Email: info.europe@strem.com

Strem Chemicals UK Ltd.
An Independent Distributor of Strem Chemicals Products
Newton Hall, Town Street
Newton, Cambridge
England CB22 7ZE
Tel: +44 (0)1223 873 028
Fax: +44 (0)1223 870207
Email: enquiries@strem.co.uk