




# Vertical Electropolished Stainless Steel Bubblers



All valves come with Kel-F stem tip rated to 200°F. High temperature valves also available. Please inquire.

## Vertical Electropolished Stainless Steel Bubblers with Fill Port



Catalog #	Description	Cylinder Vol. (ml)	A mm	B mm	C mm	D mm	E mm	Temp. Valve
96-4151	150 ml, vertical, electropolished with fill-port	150	51	108	140	83	13	Standard
96-4149	150 ml, vertical, electropolished with fill-port, a w/special valve	150	51	108	140	83	13	Standard
96-4290	300 ml, vertical, electropolished with fill-port	300	51	187	140	83	13	Standard
96-4598	600ml, vertical, electropolished, with fill-port	600	76	163	140	83	13	Standard
96-4998	1000ml, vertical, electropolished with fill-port	1000	76	254	140	83	13	Standard
96-5002	1200ml, vertical, electropolished with fill-port	1200	101	165	140	83	13	Standard
96-5003	1500ml, vertical, electropolished with fill-port	1500	101	203	140	83	13	Standard
96-5001	2000ml, vertical, electropolished with fill-port	2000	101	254	140	83	13	Standard
96-5010	2500ml, vertical, electropolished with fill-port	2500	127	203	140	83	13	Standard
96-5011	3000ml, vertical, electropolished with fill-port	3000	127	254	140	83	13	Standard
<b>Catalog # - Cylinders with High Temperature Valve – Stellite, 600°F (or 315° C)</b>								
98-0271	150 ml, vertical, electropolished with fill-port, valve has a metal seat (high temp. – rating 315°C or 600°F)	150	51	108	140	83	13	High
§98-0276	150cc, vertical, electropolished with fill port, high temperature metal seat (high temp. – rating 315°C or 600°F), dip tube cut to 1" from the head	150	51	108	140	83	13	High
98-0270	300ml, vertical, electropolished with fill-port, valve has a metal seat (high temp. – rating 315°C or 600°F)	300	51	187	140	83	13	High
96-4295	300 ml, high-temperature valves (high temp. – rating 315°C or 600°F), fill-port, no dip tube, handles rotated 90° front	300	51	187	140	83	13	High
96-4599	600ml, vertical, electropolished with fill-port, metal seat (high temp. – rating 315°C or 600°F) valve	600	76	163	140	83	13	High
96-4996	1000cc, vertical, electropolished with fill-port, metal seat (high temp. – rating 315°C or 600°F) valve	1000	76	254	140	83	13	High

§ Dip Tube cut to 1" into cylinder

See reverse side for Horizontal Stainless Steel Bubblers

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.


**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950-4098  
U.S.A.  
Tel.: (978) 499-1600  
Fax: (978) 465-3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
15, rue de l'Atome  
Zone Industrielle  
67800 BISCHHEIM France  
Tel.: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Telefon: 0 78 51/ 7 58 79  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK, Ltd.**  
41 Hills Road  
Cambridge  
England CB2 1NT  
Tel.: 0845 643 7263  
Fax: 01223 368021  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)

# STREM Horizontal Electropolished Stainless Steel Bubblers



All valves come with Kel-F stem tip rated to 200°F. High temperature valves also available. Please inquire.

## Horizontal Electropolished Stainless Steel Bubblers with Fill Port

Catalog #	Description	Cylinder Vol. (ml)	A mm	B mm	C mm	D mm	Temp. Valve
96-5151	150ml, horizontal in line, electropolished with fill-port	150	51	108	24	222	Standard
98-0273	150ml, electropolished with fill port, horizontal with valve turned 90° away from fill port	150	51	108	24	222	Standard
96-5298	300ml, horizontal in line, electropolished with fill-port	300	51	187	24	222	Standard
96-5599	600ml, horizontal in line, electropolished with fill-port	600	76	163	24	222	Standard
96-5998	1000ml, horizontal in line, electropolished with fill-port	1000	76	254	24	222	Standard
96-5004	2000ml, horizontal in line, electropolished with fill-port	2000	102	254	24	222	Standard
Catalog #	Cylinders with High Temperature Valve – Stellite, 600°F (or 315° C)						
98-0272	150ml, horizontal in line, electropolished with fill-port, valve has a metal seat (high temp. – rating 315°C or 600°F)	150	51	108	24	222	High
96-5006	300ml, horizontal in line, electropolished with fill-port valve has a metal seat (high temp. – rating 315°C or 600°F)	300	51	187	24	222	High
96-5007	600ml, horizontal in line, electropolished with fill-port valve has a metal seat (high temp. – rating 315°C or 600°F)	600	76	163	24	222	High
96-5008	1000ml, horizontal in line, electropolished with fill-port valve has a metal seat (high temp. – rating 315°C or 600°F)	1000	76	254	24	222	High



See reverse side for Vertical Stainless Steel Bubblers

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.

**Strem Chemicals, Inc.**  
 7 Mulliken Way  
 Newburyport, MA 01950-4098  
 U.S.A.  
 Tel.: (978) 499-1600  
 Fax: (978) 465-3104  
 Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
 15, rue de l'Atome  
 Zone Industrielle  
 67800 BISCHHEIM France  
 Tel.: (33) 03 88 62 52 60  
 Fax: (33) 03 88 62 26 81  
 Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
 Postfach 1215  
 77672 KEHL  
 Germany  
 Telefon: 0 78 51/ 7 58 79  
 Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK, Ltd.**  
 41 Hills Road  
 Cambridge  
 England CB2 1NT  
 Tel.: 0845 643 7263  
 Fax: 01223 368021  
 Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)