

**Stainless Steel Bubblers, vertical, electropolished with fill-port, PCTFE valve stem tip (121°C), DOT 4B, UN stamped**

Catalog #	Vol. (mL)	A mm	B mm	C mm	D mm	E mm	Temp. Valve	Special Configuration
95-4151	150	51	108	143	83	13	Standard	
95-4290	300	51	187	143	83	13	Standard	
95-4598	600	76	164	143	83	13	Standard	
95-4998	1000	76	254	143	83	13	Standard	
95-5002	1200	102	184	143	83	13	Standard	
95-5003	1500	102	223	143	83	13	Standard	
95-5001	2000	102	292	143	83	13	Standard	
95-5011	3000	152	213	140	83	13	Standard	
95-3000	150	51	108	140	83	13	Standard	replaceable-seat valves with rotated handles
95-4153	150	51	108	143	83	13	Standard	with rotated handles

**Stainless Steel Bubblers, vertical, electropolished with fill-port, high temperature valves (315°C), DOT 4B, UN stamped**

Catalog #	Vol. (mL)	A mm	B mm	C mm	D mm	E mm	Temp. Valve	Special Configuration
95-0271	150	51	108	140	83	13	High	
95-0270	300	51	187	140	83	13	High	
95-4599	600	76	164	140	83	13	High	
95-4002	1200	102	184	143	83	13	High	
95-4155	150	51	108	140	83	13	High	with rotated handles
95-0276	150	51	108	140	83	13	High	with short dip tube, cut to 1 inch from top of bubbler
95-4157	200	51	150	140	83	13	High	with rotated handles
95-0280	300	51	187	140	83	13	High	with rotated handles
95-4295	300	51	187	140	83	13	High	with rotated handles, <b>no dip tube</b>
95-4610	600	76	164	127	83	13	High	with rotated handles, <b>no dip tube</b> , stems same height

See reverse side for Horizontal Stainless Steel Bubblers


Visit [www.strem.com](http://www.strem.com) for new product information and a searchable catalog.

**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950  
U.S.A  
Tel: 978.499.1600  
Fax: 978.465.3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
15, rue de l'Atome  
Zone Industrielle  
67800 BISCHHEIM France  
Tel: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Tel: 0 78 51/ 7 58 79  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK Ltd.**  
An Independent Distributor of Strem Chemicals Products  
Newton Hall, Town Street  
Newton, Cambridge  
England CB22 7ZE  
Tel: +44 (0)1223 873 028  
Fax: +44 (0)1223 870207  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)



Stainless Steel Bubblers, horizontal, electropolished with fill-port, PCTFE valve stem tip (121°C), DOT 4B, UN stamped							
Catalog #	Vol. (mL)	A mm	B mm	C mm	D mm	Temp. Valve	Special Configuration
95-5151	150	51	108	24	222	Standard	
95-5298	300	51	187	24	222	Standard	
95-5599	600	76	164	24	222	Standard	
95-5998	1000	76	254	24	222	Standard	
95-5004	2000	152	152	24	222	Standard	
95-0273	150	51	108	24	222	Standard	with reversed handles

Stainless Steel Bubblers, horizontal, electropolished with fill-port, high temperature valves (315°C), DOT 4B, UN stamped							
Catalog #	Vol. (mL)	A mm	B mm	C mm	D mm	Temp. Valve	Special Configuration
95-0272	150	51	108	24	222	High	
95-5006	300	51	187	24	222	High	
95-5007	600	76	164	24	222	High	
95-5008	1000	76	254	24	222	High	
95-5009	1000	76	254	24	222	High	no dip tube

See reverse side for Vertical Stainless Steel Bubblers

Visit [www.strem.com](http://www.strem.com) for new product information and a searchable catalog.

**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950  
U.S.A.  
Tel: 978.499.1600  
Fax: 978.465.3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
15, rue de l'Atome  
Zone Industrielle  
67800 BISCHHEIM France  
Tel: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Tel: 0 78 51/ 7 58 79  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK Ltd.**  
An Independent Distributor of Strem Chemicals Products  
Newton Hall, Town Street  
Newton, Cambridge  
England CB22 7ZE  
Tel: +44 (0)1223 873 028  
Fax: +44 (0)1223 870207  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)