


Catalog # 19-1045 CALLERY™ Potassium hexamethyldisilazane, 15% in toluene







Technical Notes:

Lithium hexamethyldisilazane (KHMDS) is a strong non-nucleophilic, hindered amine base, with higher base strength than alkali metal alkoxides ($pK_a=26$). KHMDS is a useful reagent for wide variety of chemical reactions and transformations. Applications include alkylation, arylation, acylation, ring formation, isomerization, rearrangements, aldol condensations, Wittig and Horner-Emmons reactions and polymerization. In addition, KHMDS is able to catalyze transition metal-free reactions and to act as a ligand while reacting with a wide range of metal halides forming $M(\text{HMDS})_x$ catalysts


Transition metal-free application

1. Catalyst for the addition of primary and secondary phosphines to carbodiimides to generate substituted phosphaguanidines
2. **Total synthesis.** Catalyst used in the synthesis of 2,6-cis-disubstituted tetrahydropyran systems via intramolecular amide enolate alkylation
3. **Cross-couplings.** Catalyst for the cross-dehydrogenative couplings of hydrosilanes with amines
4. **1,4-Addition reactions.** Catalyst for the 18-crown-6 assisted direct-type 1,4-addition reactions of alkylazaarenes
5. Catalyst for the selective benzylic C–H bond addition of alkylpyridines to styrenes
6. Catalyst for the synthesis of aryl ethers by reacting alcohols/phenols (ROH) with aryl ammonium salts (ArNMe_3^+), which are readily prepared from anilines (ArNR'_2 , $\text{R}'=\text{H}$ or Me)
7. **Silylation of (poly)azines.** Catalyst for the mild and direct site-selective sp^2 C–H silylation of (poly)azines
8. Used in NHC-catalyzed enantioselective decarboxylative annulations to access dihydrobenzoxazinones and quinolones
9. **Olefination.** Catalyst for the stereoselective synthesis of ribofuranoid *exo*-glycals by olefination using ribofuranosyl sulfones
10. **Mannich Reactions.** Catalyst for the asymmetric Mannich reactions with simple amides


Tech Note (1)
Ref. (1)Tech Note (2)
Ref. (2)Tech Note (3)
Ref. (3)




Tech Note (4)
Ref. (4)




Tech Note (5)
Ref. (5)




Tech Note (6)
Ref. (6)




Tech Note (7)
Ref. (7)



Tech Note (8)
Ref. (8)



Tech Note (9)
Ref. (9)



Tech Note (10)
Ref. (10)


References:

1. [Chem. Commun., 2006, 3812](#)
2. [J. Org. Chem. 2015, 80, 3315](#)
3. [ChemCatChem 2016, 8, 1373](#)
4. [Angew. Chem. Int. Ed. 2017, 56, 4520](#)
5. [Angew. Chem. Int. Ed. 2018, 57, 1650](#)
6. [Angew. Chem. Int. Ed. 2018, 57, 3641](#)
7. [J. Am. Chem. Soc. 2019, 141, 127](#)
8. [Angew. Chem. Int. Ed. 2019, 58, 5941](#)
9. [J. Org. Chem. 2021, 86, 657](#)
10. [J. Am. Chem. Soc. 2021, 143, 5598](#)


Application with transition metals

1. **Hydrosilylation of carbonyls.** Ligand for the Fe-catalyzed hydrosilylation of carbonyl compounds
2. **Arylation.** Base additive used in Pd-Catalyzed C(sp³)-H arylation of diarylmethanes at room Temperature
3. **Selective hydrogenation.** Base additive used in Fe-catalyzed selective hydrogenation of activated amides to amines and alcohols
4. **Hydrosilylation of alkenes.** Ligand for the Co-catalyzed hydrosilylation of alkenes with tertiary silanes


- Alkylation of alcohols.** Base additive Co-catalyzed alkylation of secondary alcohols with primary alcohols via borrowing hydrogen/hydrogen autotransfer
- Hydroboration.** Ligand for the REM catalyzed hydroboration of imines and nitriles (REM = Er, Y, Dy or Gd)
- Asymmetric allylic alkylation.** Base additive in the direct intermolecular Pd-catalyzed asymmetric allylic alkylation of α -aryl cyclic vinylogous esters, enabling a straightforward enantioselective synthesis of 6-allyl-6-aryl-3-ethoxycyclohex-2-en-1-ones, common motifs embedded in numerous structurally diverse natural products




Tech Note (1)
Ref. (1)




Tech Note (2)
Ref. (2)




Tech Note (3)
Ref. (3)




Tech Note (4)
Ref. (4)



Tech Note (5)
Ref. (5)



Tech Note (6)
Ref. (6)



Tech Note (7)
Ref. (7)

References:

- [Angew. Chem. Int. Ed. 2010, 49, 10186](#)
- [J. Am. Chem. Soc. 2012, 134, 13765](#)
- [Chem. Commun., 2016, 52, 5285](#)
- [J. Am. Chem. Soc. 2017, 139, 1798](#)
- [Chem. Eur. J. 2017, 23, 12110](#)
- [Inorg. Chem. 2018, 57, 15069](#)
- [Org. Lett. 2021, 23, 920](#)

Strem Chemicals, Inc.

strem.com