



Catalog # 46-3010 Tris(dibenzylideneacetone)dipalladium(0) chloroform adduct




See also [46-0210](#) and [46-3000](#)

Technical Notes:


1. Used for Pd-catalyzed asymmetric arylation, vinylation, and Allenylation of tert-cyclobutanols via enantioselective C-C Bond cleavage.
2. Used for synthesis of chiral chromans through the Pd-catalyzed asymmetric allylic alkylation (AAA).
3. Catalyst for double N-arylation of primary amines to synthesize multisubstituted carbazoles from 2,2'-biphenylene ditriflates.
4. Palladium catalyst for regio- and enantioselective allylic alkylation of ketones through allyl enol carbonates.
5. Used for Pd-catalyzed enantioselective C-3 allylation of 3-substituted-1*H*-indoles using trialkylboranes.
6. Used for enantioselective construction of spirocyclic oxindolic cyclopentanes by Pd-catalyzed trimethylenemethane-[3+2]-cycloaddition.
7. Used for Pd-catalyzed insertion of  $\alpha$ -diazoesters into vinyl halides to generate  $\alpha,\beta$ -unsaturated  $\gamma$ -amino Esters.
8. Pd catalyst for decarboxylative asymmetric allylic alkylation of enol carbonates.
9. Palladium catalyst for asymmetric addition of oxindoles and allenes.
10. Catalyst for diastereo- and enantioselective formal [3+2]-cycloaddition between substituted vinylcyclopropanes and electron-deficient olefins.
11. Used for Pd-catalyzed asymmetric decarboxylative cycloaddition of vinyl ethylene carbonates with Michael acceptors.
12. Catalyst for enantioselective [6+4] cycloaddition of vinyl oxetanes with azadienes to access ten-membered heterocycles.




**Tech. Note (1)**  
**Ref. (1)**




**Tech. Note (2)**  
**Ref. (2)**




Tech. Note (3)  
Ref. (3)




Tech. Note (4)  
Ref. (4)




Tech. Note (5)  
Ref. (5)




Tech. Note (6)  
Ref. (6)




Tech. Note (7)  
Ref. (7)




Tech. Note (8)  
Ref. (8)




Tech. Note (9)  
Ref. (9)



Tech. Note (10)  
Ref. (10)



**Tech. Note (11)**  
**Ref. (11)**



**Tech. Note (12)**  
**Ref. (12)**

References:

1. *J. Am. Chem. Soc.*, **2003**, 125, 8862.
2. *J. Am. Chem., Soc.* **2004**, 126, 11965.
3. *J. Org. Chem.*, **2005**, 70, 413.
4. *J. Am. Chem. Soc.*, **2005**, 127, 2846.
5. *J. Am. Chem. Soc.*, **2006**, 128, 6314.
6. *J. Am. Chem. Soc.*, **2007**, 129, 12396.
7. *Angew. Chem. Int. Ed.*, **2009**, 48, 3677.
8. *J. Am. Chem. Soc.*, **2009**, 131, 18343.
9. *J. Am. Chem. Soc.*, **2011**, 133, 20611.
10. *J. Am. Chem. Soc.*, **2012**, 134, 17823.
11. *Angew. Chem. Int. Ed.*, **2014**, 53, 11257.
12. *Angew. Chem. Int. Ed.*, **2018**, 57, 1596.