




metals • inorganics • organometallics • catalysts • ligands • custom synthesis • cGMP facilities • nanomaterials




16-0215



16-0650



16-1700, 16-1706,



16-0210

Sold in collaboration with PhosphonicS Ltd. for research purposes only.

Sold in 10g and 50g units.

Metal scavengers are increasingly used as an effective way to solve metal removal problems associated with Active Pharmaceutical Ingredients (API's) and synthetic intermediates. This is an often complex technical challenge, influenced by a number of factors including the significant structural variations with API's, the polar functional groups which APIs tend to contain, the environment (solvent and pH) and potential incompatibilities within the API. PhosphonicS has designed a portfolio of silica-based materials containing a diverse range of functional groups to address effective metal scavenging.

### Properties of PhosphonicS' Broad Portfolio of Metal Scavengers

- High Affinity for a **wide range of metals**, in **different oxidation states**
- **Fast kinetics** - highly active at ambient temperatures
- **High selectivity** for the metal, meaning **minimal loss of the API**/compound
- **Broad solvent** and **pH compatibility** with both organic and aqueous formulations, reducing the requirements of time-consuming and costly solvent switches
- **No swelling** of the materials are required
- **Excellent stability** - thermal, physical, chemical and mechanical
- **Very high purity** - **eliminating issues of extractable impurities**
- **Enhanced performance**, due to **multiple functional groups** and **higher effective loadings**
- **Metal recycling** options
- **Availability on Process scale**

### Performance Benefits from using PhosphonicS' Metal Scavengers

- **Faster:** purification, batch processing, problem-solving, market introduction
- **Reduced:** Development time, valuable API losses, total manufacturing costs, environmental burden
- **Enhanced:** process, productivity, compliance

**96-6750 PhosphonicS Metals Scavenging Kit** - Contains a 10g unit of each of the following nine products

**16-0215 2-Aminoethyl sulfide ethyl Silica (PhosphonicS SEA)** white to pale yellow solid  
 Silica size: 60-200 microns Loading: 0.8 mmol/g Pore size: 60Å  
 Scavenges: Fe(III), Ni(II), Cu(II), Ru(III), Rh(I), Rh(III), Ag(I), Cd(II), Sn(II), Sn(IV), Pt(II), Hg(II), Pb(II)

**16-0650 2-Mercaptoethyl ethyl sulfide Silica (PhosphonicS SEM26)** white solid  
 Silica size: 60-200 microns Loading: 0.8 mmol/g Pore size: 60Å  
 Scavenges: Fe(III), Ni(0), Cu(II), Ru(III), Rh(I), Rh(III), Pd(0), Pd(II), Ag(I), Pt(II), Hg(II)

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.

**Strem Chemicals, Inc.**  
 7 Mulliken Way  
 Newburyport, MA 01950-4098  
 U.S.A.

Tel.: (978) 499-1600  
 Fax: (978) 465-3104  
 Email: [info@strem.com](mailto:info@strem.com)

phosphonicS0121

**Strem Chemicals, Inc.** 15,  
 rue de l'Atome  
 Zone Industrielle  
 67800 BISCHHEIM France

Tel.: (33) 03 88 62 52 60  
 Fax: (33) 03 88 62 26 81  
 Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
 Postfach 1215  
 77672 KEHL  
 Germany  
 Telefon: 0 78 51/ 7 58 79

Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK Ltd.**  
 An Independent Distributor of Strem  
 Chemicals Products  
 Newton Hall, Town Street  
 Newton, Cambridge  
 England CB22 7ZE  
 Tel: 0845 643 7263  
 Fax: 0845 643 7362  
 Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)

<b>16-1700</b>	<b>3-Mercaptopropyl ethyl sulfide Silica (60Å, high-cross linking) (PhosphonicS SPM36)</b> white solid Silica size: 60-200 microns      Loading: 0.8 mmol/g      Pore size: 60Å Scavenges: Fe(III), Ni(II), Cu(II), Ru(III), Rh(I), Rh(III), Pd(0), Pd(II), Ag(I), Cd(II), Sn(II), Sn(IV), Pt(II), Au, Hg(II), Pb(II)
<b>16-1706</b>	<b>3-Mercaptopropyl ethyl sulfide Silica (60Å, low-cross linking) (PhosphonicS SPM32)</b> white to pale yellow solid Silica size: 60-200 microns      Loading: 0.8 mmol/g      Pore size: 60Å Scavenges: Fe(III), Ni(0), Ni(II), Cu(II), Ru(III), Rh(I), Rh(III), Pd(0), Pd(II), Ag(I), Cd(II), Sn(II), Sn(IV), Pt(II), Au, Hg(II), Pb(II)
<b>16-0210</b>	<b>Triamine ethyl sulfide amide Silica (PhosphonicS STA3)</b> white solid Silica size: 60-200 microns      Loading: 0.9 mmol/g      Pore size: 60Å Scavenges: Al(III), VO(II), Cr(III), Mn(II), Fe(III), Co(II), Ni(II), Cu(II), Zn(II), Ga(III), Sr(II), Ru(III), Rh(I), Rh(III), Pd(0), Pd(II), Ag(I), Cd(II), Sn(II), Sn(IV), Pt(II), Hg(II), Pb(II), UO <sub>2</sub> (II)

Scavenger	Pd(II)	Pd(0)	Pt	Rh(I)	Rh(III)	Ru(II)	Ru(III)	Ag(I)
16-0210	✓	✓	✓	✓	✓	✓	✓	✓
16-0650	✓	✓	✓	✓	✓		✓	✓
16-1706	✓	✓	✓	✓	✓		✓	✓
16-1704								
16-1700	✓	✓	✓	✓	✓		✓	✓
16-1702								
16-1540	✓	✓	✓	✓		✓	✓	✓
16-0215	✓	✓	✓	✓	✓		✓	✓
16-0200	✓	✓	✓	✓	✓		✓	

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.

**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950-4098  
U.S.A.

Tel.: (978) 499-1600  
Fax: (978) 465-3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**,  
rue de l'Atome  
Zone Industrielle  
67800 BISCHEIM France

Tel.: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Telefon: 0 78 51/ 7 58 79

Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK Ltd.**  
An Independent Distributor of Strem  
Chemicals Products  
Newton Hall, Town Street  
Newton, Cambridge  
England CB22 7ZE  
Tel: 0845 643 7263  
Fax: 0845 643 7362  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)