

Immobilized Enzyme Kits

96-4050: CalB immo KIT

each kit contains 3x10g of each of the following items:

Strem Catalog #	Description	Immobilized on	Immobilization	Enzyme activity (PLU/g dry)
07-3130	CalB immo Plus™	DVB/methacrylate	Adsorption	>9,000
07-3148	CalB immo 8806™	Octadecyl methacrylate	Adsorption	>10,000

Appearance: White to slightly yellow spherical beads, free from foreign matter

Optimal Storage Condition: Keep in dry conditions (2-8°C).

Shelf Life: 1 year

Contents are available for individual sale.

Principal Applications

- Screening of immobilized lipases for process development
- Esterifications (regio- and stereo-selective)
- Transesterfication

- Amidation
- · Fats and oils modification

Advantages

· Fast screening in process development

• Wide selection of enzyme carriers for different applications

Availability & Shipping

Our products are available globally in various sizes including bulk shipping.

For a complete listing of the Strem Catalog products and catalysis applications, please visit strem.com or contact our team today info@ascensusspecialties.com.

Collaboration

Sold in collaboration with Purolite for research purposes only.





• CORPORATE OFFICE

Ascensus Specialties 2821 Northup Way, Suite 275 Bellevue, WA 98004

• P: +1-424-448-1679
E: ask@ascensusspecialties.com
www.ascensusspecialties.com

© 2022 ASCENSUS SPECIALTIES LLC. ALL RIGHTS RESERVED

The ASCENSUS word, design marks and all branded products are trademarks of Ascensus Specialties Inc. in the United States and elsewhere.